

Amplify and Tag Targets

For Purified DNA

- 1 Dilute 1 ng purified DNA input material to 0.2 ng/μl with nuclease-free water.
- 2 Create a master mix in the Master Mix tube.
 - ▶ PCR1 (4.7 μl)
 - ▶ FEM (0.3 μl)
 - ▶ DPMA or DPMB (5.0 μl)
- 3 Pipette to mix and then centrifuge briefly.
- 4 Add 10 μl master mix to FSP plate.
- 5 Dilute 2 μl 2800M with 98 μl nuclease-free water in a 1.5 ml microcentrifuge tube. Gently flick and then centrifuge briefly.
- 6 Add 5 μl diluted 2800M to the appropriate well.
- 7 Add 5 μl nuclease-free water to the appropriate well.
- 8 Add 5 μl diluted purified DNA to each well. Pipette to mix.
- 9 Centrifuge at 1000 × g for 30 seconds.
- 10 Transport to the post-PCR area.
- 11 Place on the thermal cycler and run the PCR1 program.

SAFE STOPPING POINT

If you are stopping, seal the plate and store at 2°C to 8°C for up to 2 days. Alternatively, leave on the thermal cycler overnight.

For Crude Lysate

- 1 Create a master mix in the Master Mix tube.
 - ▶ PCR1 (4.7 μl)
 - ▶ FEM (0.3 μl)
 - ▶ DPMA or DPMB (5.0 μl)
 - ▶ Nuclease-free water (3.0 μl)
- 2 Pipette to mix and then centrifuge briefly.

- 3 Add 13 μl master mix to FSP plate.
- 4 Dilute 2 μl 2800M with 38 μl nuclease-free water in a 1.5 ml microcentrifuge tube.
- 5 Vortex and then centrifuge briefly.
- 6 Add 2 μl diluted 2800M to the appropriate wells.
- 7 Add 2 μl nuclease-free water to the appropriate wells.
- 8 Add 2 μl diluted crude lysate sample.
- 9 Centrifuge at 1000 × g for 30 seconds.
- 10 Transport to the post-PCR area.
- 11 Place on the thermal cycler and run the PCR1 program.
- 10 Add 5 μl nuclease-free water to the wells containing reagents from step 7. Pipette to mix.
- 11 Create FTA sample master mix in the FTA Master Mix tube.
 - ▶ PCR1 (4.7 μl)
 - ▶ FEM (0.3 μl)
 - ▶ DPMA or DPMB (5.0 μl)
 - ▶ Nuclease-free water (5.0 μl)
- 12 Pipette to mix and then centrifuge briefly.
- 13 Add 15 μl FTA master mix to FTA punch in the FSP plate.
- 14 Centrifuge at 1000 × g for 30 seconds.
- 15 Transport to the post-PCR area.
- 16 Place on the thermal cycler and run the PCR1 program.

SAFE STOPPING POINT

If you are stopping, seal the plate and store at 2°C to 8°C for up to 2 days. Alternatively, leave on the thermal cycler overnight.

For FTA Card Input

- 1 Place a 1.2 mm FTA card punch into the FSP plate.
- 2 Add 100 μl 1X TBE buffer.
- 3 Place on a PCR tube storage rack.
- 4 Shake at 1800 rpm for 2 minutes.
- 5 Centrifuge at 1000 × g for 30 seconds.
- 6 Remove and discard all supernatant.
- 7 Add the following reagents to the FSP plate intended for positive and negative controls:
 - ▶ PCR1 (4.7 μl)
 - ▶ FEM (0.3 μl)
 - ▶ DPMA or DMPB (5.0 μl)
- 8 Dilute 2 μl 2800M with 98 μl nuclease-free water in a new 1.5 ml microcentrifuge tube. Gently flick and then centrifuge briefly.
- 9 Add 5 μl diluted 2800M to the wells containing reagents from step 7. Pipette to mix.

SAFE STOPPING POINT

If you are stopping, seal the plate and store at 2°C to 8°C for up to 2 days. Alternatively, leave on the thermal cycler overnight.

Enrich Targets

- 1 Centrifuge the FSP at 1000 × g for 30 seconds.
- 2 Arrange Index 1 (i7) adapters in columns 1–12.
- 3 Arrange Index 2 (i5) adapters in rows A–H.
- 4 Place the plate on the ForenSeq Index Plate Fixture.
- 5 Using a multichannel pipette, add 4 µl Index 1 (i7) adapters to each column.
- 6 Using a multichannel pipette, add 4 µl Index 2 (i5) adapters to each row.
- 7 Vortex PCR2 and then centrifuge briefly.
- 8 Add 27 µl PCR2.
- 9 Centrifuge at 1000 × g for 30 seconds.
- 10 Place on the thermal cycler and run the PCR2 program.

SAFE STOPPING POINT

If you are stopping, seal the plate and store at 2°C to 8°C for up to 7 days. Alternatively, leave on the thermal cycler overnight.

Purify Libraries

- 1 Prepare SPB according to the number of libraries you are preparing.

| Libraries | Procedure |
|-----------|---|
| < 16 | Add 50 µl SPB × the number of libraries to a 1.7 ml microcentrifuge tube. |
| 16–96 | Add [50 µl SPB × (the number of libraries/8)] + 5 µl SPB to each well of a column of a new midi plate or reagent reservoir. |
| > 96 | Add (50 µl SPB × the number of libraries) + 200 µl SPB to a multichannel reagent reservoir. |

- 2 Add 45 µl SPB to the PBP plate.
- 3 Centrifuge the FSP plate at 1000 × g for 30 seconds.
- 4 Transfer 45 µl to the PBP plate.
- 5 Shake at 1800 rpm for 2 minutes.
- 6 Incubate at room temperature for 5 minutes.
- 7 Place on the magnetic stand until liquid is clear.
- 8 Remove and discard all supernatant.
- 9 Wash two times with 200 µl 80% EtOH.
- 10 Centrifuge at 1000 × g for 30 seconds.
- 11 Place on the magnetic stand.
- 12 Use a 20 µl pipette to remove residual EtOH.
- 13 Remove from the magnetic stand.
- 14 Add 52.5 µl RSB.
- 15 Shake at 1800 rpm for 2 minutes.
- 16 Incubate at room temperature for 2 minutes.
- 17 Place on the magnetic stand until liquid is clear.
- 18 Transfer 50 µl to the PLP plate.
- 19 Centrifuge at 1000 × g for 30 seconds.

SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to one year.

Normalize Libraries

- 1 Add 46.8 µl LNA1 and 8.5 µl LNB1 to the LNA1/LNB1 Master Mix tube.
- 2 Vortex and then invert several times to mix.
- 3 Pour into a reagent reservoir.
- 4 Transfer 45 µl to the NWP plate.
- 5 Place the PLP plate on the magnetic stand until liquid is clear.
- 6 Transfer 20 µl from the PLP plate to the NWP plate.
- 7 Shake at 1800 rpm for 30 minutes.
- 8 Combine 1.8 µl HP3 and 33.3 µl nuclease-free water in a 1.5 ml microcentrifuge tube.
- 9 Add 30 µl LNS2 to the NLP plate.
- 10 Place the NWP plate on the magnetic stand and wait until the liquid is clear (~2 minutes).
- 11 Remove and discard all supernatant.
- 12 Remove from the magnetic stand.
- 13 Wash two times with 45 µl LNWI.
- 14 Remove from the magnetic stand.
- 15 Centrifuge at 1000 × g for 30 seconds.
- 16 Place on the magnetic stand until liquid is clear.
- 17 Use a 20 µl pipette to remove supernatant.
- 18 Remove from the magnetic stand.
- 19 Add 32 µl freshly prepared 0.1 N HP3.
- 20 Shake at 1800 rpm for 5 minutes.
- 21 Place on the magnetic stand until liquid is clear.
- 22 Transfer 30 µl to the NLP plate. Pipette to mix.
- 23 Centrifuge at 1000 × g for 30 seconds.

SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 30 days.

Pool Libraries

- 1 Transfer 5 μ l to a new eight-tube strip.
- 2 Transfer the contents to the PNL tube.
- 3 Vortex and then centrifuge briefly.

SAFE STOPPING POINT

If you are stopping, cap the tube and store at -25°C to -15°C for up to 30 days.

Denature and Dilute Libraries

- 1 Create an HSC denaturation reaction in the HSC mixture tube.
 - ▶ HSC (2 μ l)
 - ▶ HP3 (2 μ l)
 - ▶ Nuclease-free water (36 μ l)
- 2 Vortex and then centrifuge briefly.
- 3 Incubate at room temperature for 5 minutes.
- 4 Add 591 μ l HT1 to the DNL tube.
- 5 Transfer 7 μ l from the PNL tube to the DNL tube. Pipette to mix.
- 6 Transfer 2 μ l HSC mixture to the DNL tube. Pipette to mix.
- 7 Vortex and then centrifuge briefly.
- 8 Place on the 96°C microheating system for 2 minutes.
- 9 Invert several times to mix.
- 10 Immediately place in the ice-water bath or on the benchtop cooler for 5 minutes.
- 11 Immediately load the entire contents onto the reagent cartridge.

Acronyms

| Acronym | Definition |
|---------|--|
| 2800M | Control DNA 2800M |
| DNL | Diluted Normalized Libraries |
| DPMA | DNA Primer Mix A |
| DPMB | DNA Primer Mix B |
| FEM | Enzyme Mix |
| FSP | ForenSeq Sample Plate |
| HP3 | 2N NaOH |
| HSC | Human Sequencing Control |
| HT1 | Hybridization Buffer |
| LNA1 | Library Normalization Additives 1 |
| LNB1 | Library Normalization Beads 1 |
| LNS2 | Library Normalization Storage Buffer 2 |
| LNW1 | Library Normalization Wash 1 |
| NLP | Normalized Library Plate |
| NWP | Normalization Working Plate |
| PBP | Purification Bead Plate |
| PCR1 | PCR1 Reaction Mix |
| PCR2 | PCR2 Reaction Mix |

| Acronym | Definition |
|---------|-----------------------------|
| PLP | Purified Library Plate |
| PNL | Pooled Normalized Libraries |
| RSB | Resuspension Buffer |
| SPB | Sample Purification Beads |