

# ForenSeq™ DNA Signature Prep

## Reference Guide

For Research, Forensic, or Paternity Use Only

Revision History	3
Introduction	5
DNA Input Recommendations	6
Additional Resources	7
Protocol Introduction	8
Tips and Techniques	9
Library Prep Workflow	10
Amplify and Tag Targets	11
Enrich Targets	14
Purify Libraries	16
Normalize Libraries	18
Pool Libraries	21
Denature and Dilute Libraries	22
Supporting Information	24
Technical Assistance	35



This document and its contents are proprietary to Illumina, Inc. and its affiliates ("Illumina"), and are intended solely for the contractual use of its customer in connection with the use of the product(s) described herein and for no other purpose. This document and its contents shall not be used or distributed for any other purpose and/or otherwise communicated, disclosed, or reproduced in any way whatsoever without the prior written consent of Illumina. Illumina does not convey any license under its patent, trademark, copyright, or common-law rights nor similar rights of any third parties by this document.

The instructions in this document must be strictly and explicitly followed by qualified and properly trained personnel in order to ensure the proper and safe use of the product(s) described herein. All of the contents of this document must be fully read and understood prior to using such product(s).

FAILURE TO COMPLETELY READ AND EXPLICITLY FOLLOW ALL OF THE INSTRUCTIONS CONTAINED HEREIN MAY RESULT IN DAMAGE TO THE PRODUCT(S), INJURY TO PERSONS, INCLUDING TO USERS OR OTHERS, AND DAMAGE TO OTHER PROPERTY.

ILLUMINA DOES NOT ASSUME ANY LIABILITY ARISING OUT OF THE IMPROPER USE OF THE PRODUCT(S) DESCRIBED HEREIN (INCLUDING PARTS THEREOF OR SOFTWARE).

© 2015 Illumina, Inc. All rights reserved.

**Illumina, 24sure, BaseSpace, BeadArray, BlueFish, BlueFuse, BlueGnome, cBot, CSPPro, CytoChip, DesignStudio, Epicentre, ForenSeq, Genetic Energy, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSelect, MiSeq, MiSeqDx, MiSeq FGx, NeoPrep, NextBio, Nextera, NextSeq, Powered by Illumina, SureMDA, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode, verifi, VeriSeq**, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

For Research, Forensic, or Paternity Use Only– not for any clinical or therapeutic use in humans or animals. This product includes 2800M Control DNA manufactured by Promega Corporation for distribution by Illumina, Inc.



## Revision History

Document	Date	Description of Change
Material # 20000923 Document # 15049528 v01	September 2015	<p>Updated introduction to indicate that DPMA contains primer pairs for 58 STRs and 94 identity-informative SNPs</p> <p>Updated this document to current format for library prep documentation. Revised instructions to be more succinct</p> <p>Changed reference from PCR Product to library following amplification</p> <p>Removed reference to obsolete Experienced User Card and added reference to new protocol guide and checklist</p> <p>Removed kit box and tube part numbers</p> <p>Removed pipettes from <i>Consumables</i> as they are standard lab items</p> <p>Removed thermal cycler from pre-PCR <i>Equipment</i></p> <p>Corrected 2800M Control Alleles in the <i>Loci</i> tables for the following.</p> <ul style="list-style-type: none"> <li>• rs279844, DYS612, Y-GATA-H4, DXS10103</li> <li>• rs1805007, rs1294331, rs1413212, rs993934</li> <li>• rs1355366, rs2399332, rs1979255, rs2046361</li> <li>• rs251934, rs338882, rs1336071, rs214955</li> <li>• rs727811, rs763869, rs1463729, rs3780962</li> <li>• rs735155, rs2111980, rs2920816, rs1335873</li> <li>• rs1886510, rs354439, rs722290, rs1821380</li> <li>• rs1382387, rs729172, rs1024116, rs1736442</li> <li>• rs719366, rs1800407, rs2814778, rs3737576</li> <li>• rs1876482, rs3827760, rs1229984, rs3811801</li> <li>• rs870347, rs1871534, rs2196051, rs3814134</li> <li>• rs1079597, rs1572018, rs1800414, rs2593595</li> <li>• rs4411548, rs2042762, rs3916235, rs310644</li> </ul> <p>Removed loci rs7520386 from Identity Informative SNPs <i>Loci</i> table</p> <p>Corrected amplicon lengths for DXS8378 in the X Haplotype Markers <i>Loci</i> table</p> <p>Changed <i>Loci</i> table headers from Target Start to Amplicon Start Position and Target End to Amplicon End Position and defined positions</p> <p>Incorporated alternate procedures to prepare FTA Card into <i>Amplify and Tag Targets</i> procedures</p> <p>Added 2800M as a positive template control to prepare FTA Card</p>
Part # 15049528 Rev. D	February 2015	<ul style="list-style-type: none"> <li>• Updated introduction to indicate that DNA Primer Mix A supports 7 X haplotype markers</li> <li>• Removed 1000 µl pipettes and tips from <i>Consumables and Equipment</i></li> <li>• In the <i>Loci</i> tables: <ul style="list-style-type: none"> <li>• Moved SNPs rs16891982 and rs12913832 from phenotypic-informative SNP to ancestry-informative SNP and indicated that they are used for both predictions.</li> <li>• Corrected the vWA minimum and maximum amplicon length</li> </ul> </li> </ul>


Document	Date	Description of Change
Part # 15049528 Rev. C	January 2015	<ul style="list-style-type: none"> <li>• Updated <i>Loci</i> tables: <ul style="list-style-type: none"> <li>• Revised autosomal STR, Y haplotype marker, and X haplotype marker amplicon lengths</li> <li>• Removed X haplotype marker DXS10148</li> <li>• Changed right column heading of identity, phenotypic, and ancestry-informative SNPs to 2800M Control Alleles</li> </ul> </li> <li>• Added number of reactions supported to <i>Kit Contents</i></li> <li>• Changed catalog numbers for kit, guide, and experienced user card</li> <li>• Changed MiSeq FGx Reagent Kit name and catalog number</li> </ul>
Part # 15049528 Rev. B	September 2014	<ul style="list-style-type: none"> <li>• Modified the reagent volumes in the <i>Amplify and Tag Targets</i> and <i>Prepare FTA Card</i> master mix tables to actual reagent volumes without overage</li> <li>• Corrected locus D5S818 2800M Control alleles to 12,12</li> <li>• Updated <i>Additional Resources</i> to remove updated support page url and remove web navigation instructions and written urls</li> <li>• Separate seal and shake as separate substeps</li> <li>• Updated SDS link to <a href="http://support.illumina.com/sds.html">support.illumina.com/sds.html</a></li> </ul>
Part # 15049528 Rev. A	August 2014	Initial release.

## Introduction

This protocol explains how to prepare DNA libraries using the reagents provided in the Illumina® ForenSeq™ DNA Signature Prep Kit to genotype database or casework reference samples in a single sequencing run.

A primer mix containing a pair of tagged oligos for each target sequence is mixed with the DNA sample. PCR cycles link the tags to copies of each target to form DNA templates consisting of the regions of interest flanked by universal primer sequences. The tags are used to attach indexed adapters, which are then amplified using PCR, purified, pooled into a single tube, and then sequenced. The index sequences allow the sequencing system to separate and isolate the data generated from each sample.

Figure 1 ForenSeq DNA Signature Prep Overview



Targeted primer mixes enable analysis of autosomal, Y- and X-chromosome Short Tandem Repeat (STR) targets, identity-informative SNPs, with the option to include ancestry-informative and phenotypic-informative SNPs depending on which primer mix is used. The ForenSeq DNA Signature Prep enables analysis of these markers on gDNA ranging from high-quality single source to difficult samples. This process is done within a single reaction with integrated indexing to support sequencing of up to 96 database or 32 casework samples per run. ForenSeq DNA Signature Prep applies the long paired-end read capability and high data quality of your Illumina sequencing system.

The ForenSeq DNA Signature Prep Kit offers:

- ▶ Multiplexing—Amplify STR and SNP amplicons in a single reaction, and sequence up to 96 samples in a single sequencing run.
- ▶ Two different primer mixes:
  - ▶ DNA Primer Mix A—Contains primer pairs for 58 STRs (including 27 autosomal STRs and 7 X and 24 Y haplotype markers) and 94 identity-informative SNPs.
  - ▶ DNA Primer Mix B—Contains all markers in DNA Primer Mix A, plus primer pairs for 56 ancestry-informative SNPs and 22 phenotypic-informative SNPs (2 ancestry-informative SNPs are also used for phenotype prediction).
- ▶ Library generation—Allows for simultaneous preparation of 96 samples to generate libraries of PCR products within a single plate. Each library is a collection of amplified DNA fragments from a single sample.

## DNA Input Recommendations

It is important to quantify the input DNA and assess the DNA quality before beginning the ForenSeq DNA Signature Prep protocol. Follow these DNA input recommendations:

- ▶ 1 ng of human genomic DNA (gDNA) input is recommended.
- ▶ Use a fluorometric based method for quantification, such as qPCR.
- ▶ The ForenSeq DNA Signature Prep Kit is compatible with lysates from buccal swabs and FTA card stains as DNA input.
  - ▶ If using crude lysates, 2  $\mu$ l input material is required per sample. See *Equipment* on page 26 for recommended lysis buffers.
  - ▶ If using FTA paper, a 1.2 mm FTA card punch per sample is required.

## Additional Resources

The following documentation is available for download from the Illumina website.

Resource	Description
<i>ForenSeq DNA Signature Prep Protocol Guide</i> (document # 1000000001629)	Provides only protocol instructions. The protocol guide is intended for experienced users.
<i>ForenSeq DNA Signature Prep Checklist</i> (document # 1000000001630)	Provides a checklist of the protocol steps. The checklist is intended for experienced users.
<i>MiSeq FGx Instrument Reference Guide</i> (document # 15050524)	Provides an overview of MiSeq FGx™ instrument components and software, instructions for performing sequencing runs, and procedures for proper instrument maintenance and troubleshooting.
<i>ForenSeq Universal Analysis Software User Guide</i> (document # 15053876)	Provides an overview of ForenSeq Universal Analysis Software (UAS), which performs secondary analysis such as demultiplexing, alignment, variant calling, and report generation.

Visit the ForenSeq DNA Signature Prep Kit support page on the Illumina website for access to requirements and compatibility, additional documentation, frequently asked questions, and best practices.

## Protocol Introduction

- ▶ Review best practices before proceeding. See *Additional Resources* on page 7 for information about ForenSeq DNA Signature Prep Kit best practices on the Illumina website.
- ▶ Processing fewer than 8 samples at the same time, including positive and negative controls, can cause problems with pipetting accuracy, due to the small volumes used when preparing the master mix.
- ▶ Create a sample sheet to record the positions of each sample and index adapter. For more information, see the *ForenSeq Universal Analysis Software User Guide* (document # 15053876).



### NOTE

For indexing guidelines when pooling 16 or fewer libraries, see *Additional Resources* on page 7 for information about ForenSeq DNA Signature Prep Kit best practices on the Illumina website.

- ▶ Follow the protocol in the order shown using the specified volumes and incubation parameters.
- ▶ Confirm kit contents and make sure that you have the required equipment and consumables. For more information, see *Supporting Information* on page 24.



## Tips and Techniques

Unless a safe stopping point is specified in the protocol, proceed immediately to the next step.

### Avoiding Cross-Contamination

- ▶ When adding or transferring samples, change tips between *each sample*.
- ▶ When adding adapters or primers, change tips between *each row* and *each column*.
- ▶ Remove unused index adapter tubes from the working area.

### Sealing the Plate

- ▶ Always seal the 96-well plate before the following steps in the protocol:
  - ▶ Shaking steps
  - ▶ Vortexing steps
  - ▶ Centrifuge steps
  - ▶ Thermal cycling steps
- ▶ Apply the adhesive seal to cover the plate and seal with a rubber roller.
- ▶ Microseal 'B' adhesive seals are effective at -40°C to 110°C, and suitable for skirted or semiskirted PCR plates. Use Microseal 'B' for shaking, centrifuging, and long-term storage.
- ▶ Microseal 'A' adhesive film is effective for thermal cycling and easy to cut when using fewer than 96 wells.

### Plate Transfers

- ▶ When transferring volumes between plates, transfer the specified volume from each well of a plate to the corresponding well of the other plate.
- ▶ If beads are aspirated into the pipette tips, dispense back to the plate on the magnetic stand and wait until the liquid is clear (~2 minutes).

### Centrifugation

- ▶ Centrifuge at any step in the procedure to consolidate liquid or beads in the bottom of the well, and to prevent sample loss.

# Library Prep Workflow

Figure 2 ForenSeq DNA Signature Prep Workflow



## Amplify and Tag Targets

This process amplifies and tags the gDNA using a ForenSeq oligonucleotide primer mix with regions specific to DNA sequences upstream and downstream of STRs and SNPs.

This protocol requires Control DNA 2800M and a negative PCR amplification control (nuclease-free water) in each experiment. If these controls are not included, troubleshooting support is limited.



### NOTE

Processing fewer than 8 samples at the same time, including positive and negative controls, can affect pipetting accuracy due to the small volumes used when preparing the master mix.

### Consumables

- ▶ 2800M (Control DNA 2800M)
- ▶ One of the following:
  - ▶ DPMA (DNA Primer Mix A)
  - ▶ DPMB (DNA Primer Mix B)
- ▶ FEM (Enzyme Mix)
- ▶ PCR1 (PCR1 Reaction Mix)
- ▶ 1.5 ml microcentrifuge tubes (2)
- ▶ 96-well 0.3 ml PCR plate, skirted or semiskirted
- ▶ Human gDNA:
  - ▶ Purified DNA (1 ng per sample)
  - ▶ Crude lysate (2  $\mu$ l per sample)
  - ▶ FTA card (1.2 mm punch per sample)
- ▶ [For FTA card] 1X TBE buffer (100  $\mu$ l per FTA card punch)
- ▶ Microseal 'A' film
- ▶ Microseal 'B' adhesive seal
- ▶ Nuclease-free water
- ▶ [Optional] RNase/DNase-free 8-tube strip and caps



### NOTE

Use Microseal 'A' when sealing the plate before placing on the thermal cycler. Use Microseal 'B' for other steps that require a sealed plate.

### About Reagents

- ▶ For information on the loci detected with DPMA and DPMB, see *Loci* on page 27.

## Preparation

- 1 Prepare the following consumables.

Item	Storage	Instructions
DPMA or DPMB	-25°C to -15°C	Thaw at room temperature.
FEM	-25°C to -15°C	Thaw at room temperature. Return to storage after use.
PCR1	-25°C to -15°C	Thaw at room temperature.
1X TBE buffer	-25°C to -15°C	Thaw at room temperature.
2800M	2°C to 8°C	Let stand for 30 minutes to bring to room temperature.

- 2 Create a sample sheet to record the positions of each sample and index adapter.
- 3 Save the following PCR1 program on the thermal cycler in the post-amplification area.



**CAUTION**

Failure to use the thermal ramping mode for your thermal cycler can have an adverse effect on results. See ramping modes for selected *Thermal Cyclers* on page 26.

- ▶ Choose the preheat lid option and set to 100°C
- ▶ 98°C for 3 minutes
- ▶ 8 cycles of:
  - ▶ 96°C for 45 seconds
  - ▶ 80°C for 30 seconds
  - ▶ 54°C for 2 minutes, with specified ramping mode
  - ▶ 68°C for 2 minutes, with specified ramping mode
- ▶ 10 cycles of:
  - ▶ 96°C for 30 seconds
  - ▶ 68°C for 3 minutes, with specified ramping mode
- ▶ 68°C for 10 minutes
- ▶ Hold at 10°C



**NOTE**

The PCR1 program takes approximately 3.5 hours and can be run overnight.

- 4 Label tube and plates with a marker as follows.
  - ▶ [For Purified DNA or Crude lysate] Master Mix - 1.5 ml microcentrifuge tube
  - ▶ [For FTA card] FTA Master Mix - 1.5 ml microcentrifuge tube
  - ▶ FSP - PCR plate

## Procedure

- 1 Quantify gDNA using a fluorometric-based method.
- 2 **For purified DNA or crude lysate**—Perform the following steps:
  - a Dilute input material as follows.
    - ▶ [Purified DNA] Dilute 1 ng to 0.2 ng/μl with nuclease-free water.
    - ▶ [Crude lysate] Dilute 2 μl with 3 μl nuclease-free water.
  - b Create a master mix for 8 or more reactions in the Master Mix tube. Multiply each reagent volume by the number of reactions being prepared. Make 10% extra reagent for overage.
    - ▶ PCR1 (4.7 μl)
    - ▶ FEM (0.3 μl)
    - ▶ DPMA or DPMB (5.0 μl)
  - c Pipette to mix and then centrifuge briefly.
  - d If processing more than 8 samples, evenly distribute the master mix into each well of an 8-tube strip, and then use a multichannel pipette to dispense.
  - e Add 10 μl master mix to each well of the FSP plate.
  - f Dilute 2 μl 2800M with 98 μl nuclease-free water in a new 1.5 ml microcentrifuge tube. Gently flick the tube and then centrifuge briefly.
  - g Add 5 μl diluted 2800M as a positive template control to the appropriate well according to the sample sheet.

- h Add 5  $\mu$ l nuclease-free water as a negative PCR amplification control to the appropriate well according to the sample sheet.
  - i Add samples to the wells according to the sample sheet. Pipette to mix.
    - ▶ [Purified DNA] Add 5  $\mu$ l diluted DNA (0.2 ng/ $\mu$ l) sample to each well.
    - ▶ [Crude lysate] Add 5  $\mu$ l diluted crude lysate sample to each well.
- 3 **For FTA card input material**—Perform the following steps:
- a Place a 1.2 mm FTA card punch into each well of the FSP plate according to the sample sheet.
  - b Add 100  $\mu$ l 1X TBE buffer.
  - c Place on a PCR tube storage rack.
  - d Shake at 1800 rpm for 2 minutes.
  - e Centrifuge at 1000  $\times$  g for 30 seconds.
  - f Remove and discard all supernatant.
  - g Add the following to the FTA Master Mix tube.
    - ▶ PCR1 (4.7  $\mu$ l)
    - ▶ FEM (0.3  $\mu$ l)
    - ▶ DPMA or DMPB (5.0  $\mu$ l)
    - ▶ Nuclease-free water (5.0  $\mu$ l)
  - h Add 15  $\mu$ l FTA master mix to an FTA punch.
  - i Add 15  $\mu$ l FTA master mix to an empty well as a no template control.
  - j Dilute 2  $\mu$ l 2800M with 98  $\mu$ l nuclease-free water in a new 1.5 ml microcentrifuge tube. Gently flick the tube and then centrifuge briefly.
  - k Add 5  $\mu$ l diluted 2800M as a positive template control to the appropriate well according to the sample sheet.
- 4 Centrifuge at 1000  $\times$  g for 30 seconds.
- 5 Transport to the post-PCR area.
- 6 Place on the thermal cycler and run the PCR1 program.

### SAFE STOPPING POINT

If you are stopping, seal the plate and store at 2°C to 8°C for up to 2 days. Alternatively, leave on the thermal cycler overnight.

# Enrich Targets

This process amplifies the DNA and adds Index 1 (i7) adapters, Index 2 (i5) adapters, and sequences required for cluster amplification.

The index adapters tag DNA with a unique combination of index sequences, which allow data from each tagged library to be separated during later analysis.



**NOTE**

This procedure is described using a 96-well PCR plate. However, when processing 8 libraries, it can be performed with an 8-tube strip.

## Consumables

- ▶ ForenSeq Index Plate Fixture Kit
- ▶ Index 1 (i7) adapters and orange tube caps
- ▶ Index 2 (i5) adapters and white tube caps
- ▶ PCR2 (PCR2 Reaction Mix)
- ▶ 1.7 ml microcentrifuge tubes (1 per index adapter tube)
- ▶ Microseal 'A' film
- ▶ Microseal 'B' adhesive seal



**NOTE**

Use Microseal 'A' when sealing the plate before placing on the thermal cycler. Use Microseal 'B' for other steps that require a sealed plate.

## About Reagents

- ▶ If processing more than 8 libraries at the same time, evenly distribute PCR2 to each well of an 8-tube strip, and then use a multichannel pipette to dispense.
- ▶ Add PCR2 slowly to each well to avoid creating air bubbles.

## Preparation

- 1 Prepare the following consumables.


Item	Storage	Instructions
Index adapters (i5 and i7)	-25°C to -15°C	Only remove adapters being used. Thaw at room temperature for 20 minutes. Vortex each tube to mix. Centrifuge briefly using a 1.7 ml Eppendorf tube.
PCR2	-25°C to -15°C	Thaw at room temperature.

- 2 Save the following PCR2 program on the thermal cycler:
  - ▶ Choose the preheat lid option and set to 100°C
  - ▶ 98°C for 30 seconds
  - ▶ 15 cycles of:
    - ▶ 98°C for 20 seconds
    - ▶ 66°C for 30 seconds
    - ▶ 68°C for 90 seconds
  - ▶ 68°C for 10 minutes
  - ▶ Hold at 10°C

## Procedure

- 1 Centrifuge at  $1000 \times g$  for 30 seconds.
- 2 Arrange Index 1 (i7) adapters in columns 1–12 of the ForeSeq Index Plate Fixture.
- 3 Arrange Index 2 (i5) adapters in rows A–H of the ForeSeq Index Plate Fixture.
- 4 Place the plate on the ForeSeq Index Plate Fixture.

**Figure 3** ForeSeq Index Plate Fixture (96 libraries)



- A** Columns 1–12: Index 1 (i7) adapters (orange caps)
- B** Rows A–H: Index 2 (i5) adapters (white caps)
- C** FSP plate

- 5 Using a multichannel pipette, add  $4 \mu\text{l}$  of each Index 1 (i7) adapter to each column. Replace the cap on each i7 adapter tube with a new orange cap.
- 6 Using a multichannel pipette, add  $4 \mu\text{l}$  of each Index 2 (i5) adapter to each row. Replace the cap on each i5 adapter tube with a new white cap.
- 7 Vortex PCR2 and then centrifuge briefly.
- 8 Add  $27 \mu\text{l}$  PCR2 to each well.
- 9 Centrifuge at  $1000 \times g$  for 30 seconds.
- 10 Place on the preprogrammed thermal cycler and run the PCR2 program.

### SAFE STOPPING POINT

If you are stopping, seal the plate and store at  $2^{\circ}\text{C}$  to  $8^{\circ}\text{C}$  for up to 7 days. Alternatively, leave on the thermal cycler overnight.

# Purify Libraries

This process uses SPB (Sample Purification Beads) to purify the amplified libraries from the other reaction components.

## Consumables

- ▶ RSB (Resuspension Buffer)
- ▶ SPB (Sample Purification Beads)
- ▶ 96-well 0.3 ml PCR plate, skirted or semiskirted
- ▶ 96-well midi plates (1 + 1 if processing 16–96 libraries)
- ▶ Freshly prepared 80% ethanol (EtOH)
- ▶ Microseal 'B' adhesive seals
- ▶ RNase/DNase-free reagent reservoirs (1 + 1 if processing more than 96 libraries)

## About Reagents

- ▶ Vortex SPB before each use.
- ▶ Vortex SPB frequently to make sure that beads are evenly distributed.
- ▶ Aspirate and dispense SPB slowly due to the viscosity of the solution.

## Preparation

- 1 Prepare the following consumables.

Item	Storage	Instructions
RSB	2°C to 8°C	Let stand for 30 minutes to bring to room temperature.
SPB	2°C to 8°C	Let stand for 30 minutes to bring to room temperature.

- 2 Label plates with a marker as follows.
  - ▶ PBP - midi plate
  - ▶ PLP - PCR plate

## Procedure

- 1 Do the following depending on the number of libraries your are preparing:

Number of Libraries	Procedure
< 16	Add 50 µl SPB × the number of libraries to a 1.7 ml microcentrifuge tube.
16–96	Add [50 µl SPB × (the number of libraries/8)] + 5 µl SPB to each well of a column of a new midi plate.
> 96	Add (50 µl SPB × the number of libraries) + 200 µl SPB to a multichannel reagent reservoir.

- 2 Add 45 µl SPB to each well of the PBP plate according to the sample sheet.
- 3 Centrifuge the FSP plate at 1000 × g for 30 seconds.
- 4 Transfer 45 µl to the corresponding well of the PBP plate, according to the sample sheet.
- 5 Shake at 1800 rpm for 2 minutes.
- 6 Incubate at room temperature for 5 minutes.



- 7 Place on the magnetic stand and wait until the liquid is clear (~2 minutes).
- 8 Remove and discard all supernatant from each well.
- 9 Wash 2 times as follows.
  - a Add 200  $\mu$ l freshly prepared 80% EtOH to each well.
  - b Incubate on the magnetic stand for 30 seconds.
  - c Remove and discard all supernatant from each well.
- 10 Centrifuge at 1000  $\times$  g for 30 seconds.
- 11 Place on the magnetic stand.
- 12 Use a 20  $\mu$ l pipette to remove residual EtOH from each well.
- 13 Add 52.5  $\mu$ l RSB to each well.
- 14 Shake at 1800 rpm for 2 minutes. If the beads are not resuspended, pipette to mix or repeat shake at 1800 rpm for 2 minutes.
- 15 Incubate at room temperature for 2 minutes.
- 16 Place on the magnetic stand and wait until the liquid is clear (~2 minutes).
- 17 Transfer 50  $\mu$ l to the corresponding well of the PLP plate.
- 18 Centrifuge at 1000  $\times$  g for 30 seconds.

#### **SAFE STOPPING POINT**

If you are stopping, seal the plate and store at -25°C to -15°C for up to one year.

# Normalize Libraries

This process prepares DNA libraries for cluster generation to make sure that libraries of varying yields are equally represented within the sequencing run. This process assures that samples with varying input amounts or sample types achieve consistent cluster density to optimize the resolution of individual samples when pooled together. By normalizing the concentration of the libraries, while preserving the content of each library, post-PCR quantification and individual PCR product normalization are not necessary.

## Consumables

- ▶ HP3 (2N-NaOH)
- ▶ LNA1 (Library Normalization Additives 1)
- ▶ LNB1 (Library Normalization Beads 1)
- ▶ LNS2 (Library Normalization Storage Buffer 2)
- ▶ LNW1 (Library Normalization Wash 1)
- ▶ 1.5 ml microcentrifuge tube
- ▶ One of the following:
  - ▶ 1.5 ml microcentrifuge tube
  - ▶ 15 ml conical tube
- ▶ 96-well 0.3 ml PCR plate, skirted or semiskirted
- ▶ 96-well midi plate
- ▶ Microseal 'B' adhesive seals
- ▶ Nuclease-free water
- ▶ RNase/DNase-free Reagent Reservoir



### WARNING

This set of reagents contains formamide, an aliphatic amide that is a probable reproductive toxin.

LNA1 and LNW1 contain  $\beta$ -mercaptoethanol and prolonged exposure can be toxic to the nervous system and cause organ damage.

Perform this procedure in a hood or well-ventilated area if desired. Personal injury can occur through inhalation, ingestion, skin contact, and eye contact. Dispose of containers and any unused contents in accordance with the governmental safety standards for your region. For more information, see how to access safety data sheets (SDSs) in *Technical Assistance* on page 35.

Supernatant, excess LNA1/LNB1 Master Mix, and tips used to pipette LNA1 and LNB1 are hazardous waste. Discard in accordance with the governmental safety standards for your region.

## About Reagents

- ▶ After vortexing, hold LNA1 in front of a light and make sure that no crystals are present and all precipitate has dissolved.
- ▶ After vortexing, make sure that LNB1 beads are well-resuspended and no pellet remains at the bottom of the tube.
- ▶ It is critical to resuspend the LNB1 bead pellet at the bottom of the tube. Use a 1000  $\mu$ l pipette to make sure that the beads are homogeneously resuspended and that there is no bead mass at the bottom of the tube. Resuspension is essential for achieving consistent cluster density to optimize the resolution of individual libraries when pooled together.
- ▶ The library that remains in the PLP plate can be stored. Seal the PLP plate and store at  $-25^{\circ}\text{C}$  to  $-15^{\circ}\text{C}$  for up to 1 year.

## Preparation

- 1 Prepare the following consumables.

Item	Storage	Instructions
HP3	-25°C to -15°C	Thaw at room temperature.
LNA1	-25°C to -15°C	Thaw at room temperature. Vortex with intermittent inversion.
LNB1	2°C to 8°C	Let stand for 30 minutes to bring to room temperature. Vortex for at least 1 minute, inverting 5 times every 15 seconds. Pipette to mix until the bead pellet at the bottom is resuspended.
LNW1	2°C to 8°C	Let stand for 30 minutes to bring to room temperature.
LNS2	15°C to 30°C	Remove from storage.

- 2 Label tubes and plates with a marker as follows.
  - ▶ LNA1/LNB1 Master Mix - 1.5 ml microcentrifuge tube or 15 ml conical tube
  - ▶ NWP - midi plate
  - ▶ NLP - PCR plate
- 3 Dedicate separate hazardous waste disposal containers for liquids and solids.

## Procedure

- 1 Create a master mix in the LNA1/LNB1 Master Mix tube.
  - ▶ LNA1 (46.8 µl) (Example: 8 reactions 374 µl)
  - ▶ LNB1 (8.5 µl) (Example: 8 reactions 68 µl)
- 2 Vortex and then invert the tube several times to mix.
- 3 Pour into a reagent reservoir.
- 4 Transfer 45 µl to each well of the NWP plate that will contain a library according to the sample sheet.
- 5 Place the PLP plate on the magnetic stand and wait until the liquid is clear (~2 minutes).
- 6 Transfer 20 µl from each well of the PLP plate to the corresponding well of the NWP plate.
- 7 Shake at 1800 rpm for 30 minutes.
- 8 Prepare 0.1 N HP3 in a new 1.5 ml microcentrifuge tube.
  - ▶ Nuclease-free water (33.3 µl) (Example: 8 reactions 266 µl)
  - ▶ HP3 (1.8 µl) (Example: 8 reactions 14 µl)
- 9 Invert the tube several times to mix.
- 10 Add 30 µl LNS2 to each well of the NLP plate that will contain a library according to the sample sheet
- 11 Place the NWP plate on the magnetic stand and wait until the liquid is clear (~2 minutes).
- 12 Remove and discard all supernatant from each well.
- 13 Remove from the magnetic stand.

- 14 Add 45  $\mu$ l LNW1 to each well.
- 15 Shake at 1800 rpm for 5 minutes.
- 16 Place on the magnetic stand and wait until the liquid is clear (~2 minutes).
- 17 Remove and discard all supernatant from each well.
- 18 Repeat steps 14–17 for a total of 2 washes.
- 19 Remove from the magnetic stand.
- 20 Centrifuge at  $1000 \times g$  for 30 seconds.
- 21 Place on the magnetic stand and wait until the liquid is clear (~2 minutes).
- 22 Use a 20  $\mu$ l pipette to remove residual supernatant from each well.
- 23 Remove from the magnetic stand.
- 24 Add 32  $\mu$ l freshly prepared 0.1 N HP3 to each well.
- 25 Shake at 1800 rpm for 5 minutes. If the beads are not resuspended, pipette to mix or repeat shake at 1800 rpm for 5 minutes.
- 26 Place on the magnetic stand and wait until the liquid is clear (~2 minutes).
- 27 Transfer 30  $\mu$ l to the corresponding well of the NLP plate. Pipette to mix.
- 28 Centrifuge at  $1000 \times g$  for 30 seconds.

#### **SAFE STOPPING POINT**

If you are stopping, seal the plate and store at  $-25^{\circ}\text{C}$  to  $-15^{\circ}\text{C}$  for up to 30 days.

## Pool Libraries

This process combines equal volumes of normalized library to create a pool of libraries that are sequenced together on the same flow cell.

### Consumables

- ▶ 1.5 ml microcentrifuge tube
- ▶ Microseal 'B' adhesive seal
- ▶ RNase/DNase-free 8-tube strip and caps

### Preparation

- 1 Determine which libraries to pool for sequencing.
  - ▶ For casework samples, do not pool more than 32 libraries, including the positive and negative controls.
  - ▶ For database samples, do not pool more than 96 libraries.
- 2 Label the tube PNL with a marker.

### Procedure

- 1 Transfer 5  $\mu$ l of each library to a new 8-tube strip.
- 2 Store the plate in the post-PCR area at -25°C to -15°C for up to 30 days.
- 3 Transfer the contents of each well of the 8-tube strip to the PNL tube.
- 4 Vortex and then centrifuge briefly.

### SAFE STOPPING POINT

If you are stopping, cap the tube and store at -25°C to -15°C for up to 30 days.

## Denature and Dilute Libraries

This process dilutes the libraries in HT1 (Hybridization Buffer), adds HSC (Human Sequencing Control), and heat denatures the libraries in preparation for sequencing.



### NOTE

Perform this process immediately before loading the library onto the reagent cartridge to ensure efficient template loading on the flow cell.

### Consumables

- ▶ HP3 (2N-NaOH)
- ▶ HSC (Human Sequencing Control)
- ▶ 1.5 ml microcentrifuge tubes (2)
- ▶ MiSeq FGx Reagent Kit contents:
  - ▶ HT1 (Hybridization Buffer)
  - ▶ Reagent cartridge
- ▶ Nuclease-free water

### About Reagents

- ▶ Follow *Prepare the Reagent Cartridge* instructions in the *MiSeq FGx Instrument Reference Guide* (document # 15050524).

## Preparation

- 1 Prepare the following consumables.

Item	Storage	Instructions
HP3	-25°C to -15°C	Thaw at room temperature.
HSC	-25°C to -15°C	Thaw at room temperature.
HT1	-25°C to -15°C	Thaw at room temperature.
Reagent cartridge	-25°C to -15°C	Thaw at room temperature.

- 2 Preheat the microheating system to 96°C.
- 3 Prepare either of the following:
  - ▶ Remove a tube benchtop cooler from -25°C to -15°C storage or ice bucket.
  - ▶ Prepare an ice-water bath by combining 3 parts ice and 1 part nuclease-free water.
- 4 Label tubes with a marker as follows.
  - ▶ HSC mixture
  - ▶ DNL

## Procedure

- 1 Create an HSC denaturation reaction in the HSC mixture tube.
  - ▶ HSC (2 µl)
  - ▶ HP3 (2 µl)
  - ▶ Nuclease-free water (36 µl)
- 2 Vortex and then centrifuge briefly.
- 3 Incubate at room temperature for 5 minutes.
- 4 Add 591 µl HT1 to the DNL tube.

- 5 Transfer 7  $\mu$ l from the PNL tube to the DNL tube. Pipette to mix.
- 6 Cap the PNL tube and store at -25°C to -15°C for up to 30 days. Exceeding 30 days in storage results in a significant reduction of cluster density.
- 7 Transfer 2  $\mu$ l HSC mixture to the DNL tube. Pipette to mix. Do not store HSC mixture long term, which results in a significant reduction of cluster density.
- 8 Vortex and then centrifuge briefly.
- 9 Place on the 96°C microheating system for 2 minutes.
- 10 Invert the tube several times to mix.
- 11 Immediately place in the ice-water bath or on the -25°C to -15°C benchtop cooler for 5 minutes.
- 12 Immediately load the entire contents onto the reagent cartridge. For more information, see the *Load Sample Libraries Onto Cartridge* instructions in the *MiSeq FGx Instrument Reference Guide* (document # 15050524).

## Supporting Information

The protocols provided in this guide assume that you are familiar with the contents of this section and that you have the required equipment and consumables.

### Kit Contents

Make sure that you have all the ForenSeq DNA Signature Prep Kit (catalog # TG-450-1001) reagents identified in this section before starting the protocol. Each kit supports 384 reactions and contains 3 boxes.

#### Pre-PCR Box 1, Store as specified

Quantity	Reagent	Description	Storage Temperature
2	2800M	Control DNA 2800M	2°C to 8°C
8	PCR1	PCR1 Reaction Mix	-25°C to -15°C
8	FEM	Enzyme Mix	-25°C to -15°C
8	DPMA	DNA Primer Mix A	-25°C to -15°C
8	DPMB	DNA Primer Mix B	-25°C to -15°C

#### Post-PCR Box 2, Store at -25°C to -15°C

Quantity	Reagent	Description
4	LNA1	Library Normalization Additives 1
4	LNS2	Library Normalization Storage Buffer 2
8	LNW1	Library Normalization Wash 1
3	HP3	HP3 2N-NaOH
8	PCR2	PCR2 Reaction Mix
1	HSC	Human Seq Control
1	A501	A501 Index Adapter
1	A502	A502 Index Adapter
1	A503	A503 Index Adapter
1	A504	A504 Index Adapter
1	A505	A505 Index Adapter
1	A506	A506 Index Adapter
1	A507	A507 Index Adapter
1	A508	A508 Index Adapter
1	R701	R701 Index Adapter
1	R702	R702 Index Adapter
1	R703	R703 Index Adapter
2	R704	R704 Index Adapter
1	R705	R705 Index Adapter
1	R706	R706 Index Adapter
1	R707	R707 Index Adapter
1	R708	R708 Index Adapter
1	R709	R709 Index Adapter
1	R710	R710 Index Adapter
1	R711	R711 Index Adapter
1	R712	R712 Index Adapter
1	—	i7 Index Tube Caps, Orange
1	—	i5 Index Tube Caps, White



## Post-PCR Box 3, Store at 2°C to 8°C

Quantity	Reagent	Description
4	LNB1	Library Normalization Beads 1
1	RSB	Resuspension Buffer
2	SPB	Sample Purification Beads

## Consumables and Equipment

Make sure that you have the required user-supplied consumables and equipment before starting the protocol.

The protocol has been optimized and validated using the items listed. Comparable performance is not guaranteed when using alternate consumables and equipment.

### Consumables

Consumable	Supplier
1.5 ml microcentrifuge tubes	General lab supplier
1.7 ml microcentrifuge tubes	General lab supplier
15 ml conical tube	General lab supplier
20 µl barrier pipette tips	General lab supplier
200 µl barrier pipette tips	General lab supplier
96-well 0.3 ml semiskirted PCR plates	Eppendorf Twin-Tec, part # 951020303 or VWR, part # 89136-706
96-well storage plates, round well, 0.8 ml ('midi' plate)	Fisher Scientific, part # AB-0859
Ethanol 200 proof (absolute) for molecular biology (500 ml)	Sigma-Aldrich, part # E7023
Microseal 'A' film	Bio-Rad, part # MSA-5001
Microseal 'B' adhesive seals	Bio-Rad, part # MSB-1001
MiSeq FGx Reagent Kit	Illumina, catalog # TG-143-1001
Nuclease-free water	General lab supplier
PCR tube storage rack (If using an FTA card as input material)	VWR, part # 80086-074
One of the following, if using crude lysates as input material: <ul style="list-style-type: none"> <li>• QuickExtract DNA Extraction Solution</li> <li>• SwabSolution Kit</li> </ul>	<ul style="list-style-type: none"> <li>• Epicentre, catalog # QE09050</li> <li>• Promega, catalog # DC8271</li> </ul>
RNase/DNase-free 8-tube strips and caps	General lab supplier
RNase/DNase-free multichannel reagent reservoirs, disposable	Labcor, part # 730-001

## Equipment

Equipment	Supplier/Description	Pre-PCR	Post-PCR
1.5 ml tube benchtop cooler	VWR, catalog # 414004-286		X
96-well thermal cycler (with heated lid)	See <i>Thermal Cyclers</i> .		X
Benchtop microcentrifuge	General lab supplier	X	X
ForenSeq Index Plate Fixture	Illumina, catalog # FC-451-1001		X
Magnetic stand-96	Life Technologies, part # AM10027 or general lab supplier		X
Microplate centrifuge	General lab supplier	X	X
1.5 ml 96-well heating system	General lab supplier		X
One of the following high-speed thermal mixers: <ul style="list-style-type: none"> <li>• BioShake iQ</li> <li>• BioShake XP</li> </ul>	Q instruments, catalog #: <ul style="list-style-type: none"> <li>• 1808-0506</li> <li>• 1808-0505</li> </ul>		X
Vortexer	General lab supplier	X	X

## Thermal Cyclers

The following table lists the recommended settings for the thermal cycler. If your lab has a thermal cycler that is not listed, validate the thermal cycler before performing the protocol.

Thermal Cycler	Temp Mode	Lid Temp	Vessel Type	Ramp Mode
ABI LTI thermal cycler 9700	9600 emulation	Heated	Polypropylene plates and tubes	8%
Bio-Rad	Calculated	Heated, Constant at 100°C	Polypropylene plates and tubes	0.2°C per second
Eppendorf Mastercycler Pro S	Gradient S, Simulated Tube	Heated	Plate	2%

## Index Sequences

The kit contains the following the index adapter sequences.

Index 1 (i7)	Sequence	Index 1 (i7)	Sequence
R701	ATCACGAT	R707	CAGATCAT
R702	CGATGTAT	R708	ACTTGAAT
R703	TTAGGCAT	R709	GATCAGAT
R704	TGACCAAT	R710	TAGCTTAT
R705	ACAGTGAT	R711	GGCTACAT
R706	GCCAATAT	R712	CTTGTAAT

Index 2 (i5)	Sequence
A501	TGAACCTT
A502	TGCTAAGT
A503	TGTTCTCT
A504	TAAGACAC
A505	CTAATCGA
A506	CTAGAACA
A507	TAAGTICC
A508	TAGACCTA

## Loci



### NOTE

- The amplicon length does not include 120 bp for Illumina adapter sequences. The amplicon start and end positions are the 1-based endpoints of the entire amplicon including the sequence matching primers on the hg19 human reference genome.
- All the loci in DNA Primer Mix A are also included in DNA Primer Mix B.

The following loci are detected using DNA Primer Mix A or DNA Primer Mix B.

Amelogenin—A genetic marker that confirms the gender of the donor of the biological sample. Its size range is 106–112 bp and the control DNA is male.

**Table 1** Identity Informative SNPs

Locus	Amplicon Length (bp)	Chromosome	Amplicon Start Position	Amplicon End Position	2800M Control Alleles
rs10495407	109	1	238439234	238439342	G
rs1294331	85	1	233448359	233448443	GA
rs1413212	64	1	242806767	242806830	G
rs1490413	98	1	4367256	4367353	A
rs560681	90	1	160786641	160786730	AG
rs891700	115	1	239881850	239881964	AG
rs1109037	118	2	10085691	10085808	G
rs12997453	100	2	182413195	182413294	A
rs876724	119	2	114945	115063	C
rs907100	115	2	239563542	239563656	CG
rs993934	120	2	124109120	124109239	C
rs1355366	119	3	190806041	190806159	AG
rs1357617	120	3	961696	961815	AT
rs2399332	157	3	110300999	110301155	AC

Locus	Amplicon Length (bp)	Chromosome	Amplicon Start Position	Amplicon End Position	2800M Control Alleles
rs4364205	98	3	32417576	32417673	G
rs6444724	120	3	193207306	193207425	T
rs1979255	102	4	190318007	190318108	G
rs2046361	120	4	10968994	10969113	A
rs279844	167	4	46329584	46329750	AT
rs6811238	120	4	169663541	169663660	G
rs13182883	169	5	136633252	136633420	AG
rs159606	104	5	17374845	17374948	A
rs251934	97	5	174778619	174778715	T
rs338882	157	5	178690599	178690755	C
rs717302	110	5	2879333	2879442	G
rs13218440	170	6	12059928	12060097	AG
rs1336071	120	6	94537182	94537301	G
rs214955	120	6	152697629	152697748	G
rs727811	115	6	165045254	165045368	A
rs321198	165	7	137029715	137029879	T
rs6955448	120	7	4310285	4310404	CT
rs737681	120	7	155990742	155990861	T
rs917118	109	7	4456953	4457061	C
rs10092491	116	8	28411037	28411152	CT
rs2056277	104	8	139399038	139399141	C
rs4606077	151	8	144656710	144656860	CT
rs763869	85	8	1375576	1375660	CT
rs1015250	117	9	1823702	1823818	G
rs10776839	103	9	137417271	137417373	G
rs1360288	119	9	128967994	128968112	C
rs1463729	99	9	126881396	126881494	GA
rs7041158	115	9	27985907	27986021	C
rs3780962	94	10	17193284	17193377	T
rs735155	170	10	3374133	3374302	A
rs740598	120	10	118506839	118506958	AG
rs826472	153	10	2406511	2406663	T
rs964681	105	10	132698394	132698498	CT
rs10488710	118	11	115207134	115207251	CG
rs1498553	111	11	5708981	5709091	CT
rs2076848	118	11	134667502	134667619	AT
rs901398	90	11	11096173	11096262	T
rs10773760	99	12	130761623	130761721	AG
rs2107612	103	12	888262	888364	AG
rs2111980	94	12	106328186	106328279	G
rs2269355	65	12	6945881	6945945	C
rs2920816	157	12	40862976	40863132	T
rs1058083	76	13	100038193	100038268	AG
rs1335873	109	13	20901665	20901773	T
rs1886510	116	13	22374646	22374761	CT
rs354439	170	13	106938320	106938489	T
rs1454361	118	14	25850765	25850882	AT

Locus	Amplicon Length (bp)	Chromosome	Amplicon Start Position	Amplicon End Position	2800M Control Alleles
rs4530059	170	14	104769099	104769268	G
rs722290	101	14	53216686	53216786	G
rs873196	114	14	98845506	98845619	CT
rs1528460	115	15	55210664	55210778	T
rs1821380	118	15	39313343	39313460	G
rs8037429	63	15	53616876	53616938	T
rs1382387	89	16	80106318	80106406	GT
rs2342747	104	16	5868645	5868748	AG
rs430046	119	16	78016980	78017098	C
rs729172	104	16	5606153	5606256	C
rs740910	113	17	5706552	5706664	A
rs8078417	143	17	80461847	80461989	CT
rs938283	98	17	77468433	77468530	T
rs9905977	170	17	2919324	2919493	G
rs1024116	98	18	75432317	75432414	A
rs1493232	75	18	1127945	1128019	A
rs1736442	153	18	55225698	55225850	G
rs9951171	119	18	9749789	9749907	G
rs576261	76	19	39559780	39559855	AC
rs719366	170	19	28463281	28463450	T
rs1005533	158	20	39487066	39487223	A
rs1031825	126	20	4447416	4447541	C
rs1523537	117	20	51296076	51296192	C
rs445251	119	20	15124865	15124983	CG
rs221956	97	21	43606933	43607029	C
rs2830795	114	21	28608089	28608202	A
rs2831700	79	21	29679639	29679717	A
rs722098	101	21	16685561	16685661	AG
rs914165	156	21	42415865	42416020	AG
rs1028528	78	22	48362256	48362333	AG
rs2040411	68	22	47836378	47836445	A
rs733164	120	22	27816711	27816830	AG
rs987640	120	22	33559450	33559569	AT

Table 2 Autosomal STRs

Locus	Minimum Amplicon Length (bp)	Maximum Amplicon Length (bp)	Chromosome	2800M Control Alleles
D1S1656	141	189	1	12,13
TPOX	85	145	2	11,11
D2S441	144	180	2	10,14
D2S1338	114	182	3	22,25
D3S1358	138	186	3	17,18
D4S2408	93	117	4	9,9
FGA	150	306	22	20,23
D5S818	102	150	5	12,12
CSF1PO	85	129	5	12,12

Locus	Minimum Amplicon Length (bp)	Maximum Amplicon Length (bp)	Chromosome	2800M Control Alleles
D6S1043	163	227	6	12,20
D7S820	135	179	7	8,11
D8S1179	86	138	8	14,15
D9S1122	108	140	9	12,12
D10S1248	128	172	10	13,15
TH01	100	148	11	6,9,3
vWA	132	192	12	16,19
D12S391	237	281	12	18,23
D13S317	138	186	13	9,11
PentaE	362	467	15	7,14
D16S539	132	180	16	9,13
D17S1301	114	142	17	11,12
D18S51	140	227	18	16,18
D19S433	154	212	19	13,14
D20S482	125	165	20	14,15
D21S11	158	276	20	29,31.2
PentaD	209	293	21	12,13
D22S1045*	193	229	21	16,16

\* Interpret the locus D22S1045 with caution. More imbalance in read counts may be seen between alleles of a heterozygote than observed at other loci. Consider multilocus genotype when determining the presence of a DNA mixture.

Table 3 Y Haplotype Markers

Locus	Minimum Amplicon Length (bp)	Maximum Amplicon Length (bp)	Chromosome	2800M Control Alleles
DYF387S1	123	155	Y	37,38
DYS19	261	345	Y	14
DYS385a-b	316	354	Y	13,16
DYS389I	231	275	Y	14
DYS389II	255	299	Y	31
DYS390	242	286	Y	24
DYS391	123	167	Y	10
DYS392	346	358	Y	13
DYS437	178	210	Y	14
DYS438	144	169	Y	9
DYS439	199	239	Y	12
DYS448	288	324	Y	19
DYS460	356	380	Y	11
DYS481	102	129	Y	22
DYS505	154	194	Y	11
DYS522	294	334	Y	12
DYS533	198	258	Y	12
DYS549	214	262	Y	13
DYS570	162	214	Y	17
DYS576	183	235	Y	18
DYS612	215	248	Y	29

Locus	Minimum Amplicon Length (bp)	Maximum Amplicon Length (bp)	Chromosome	2800M Control Alleles
DYS635	214	306	Y	21
DYS643	115	115	Y	10
Y-GATA-H4	151	203	Y	11

Table 4 X Haplotype Markers

Locus	Minimum Amplicon Length (bp)	Maximum Amplicon Length (bp)	Chromosome	2800M Control Alleles
DXS10074	211	309	X	21
DXS10103	161	185	X	18
DXS10135	228	334	X	28
DXS7132	176	208	X	13
DXS7423	147	215	X	15
DXS8378	430	462	X	12
HPR1B	193	237	X	12

In addition to the SNPs and haplotype markers above, the following loci are also detected when using DNA Primer Mix B. These are not present in DNA Primer Mix A.

Table 5 Phenotypic Informative SNPs

Locus	Amplicon Length (bp)	Chromosome	Amplicon Start Position	Amplicon End Position	2800M Control Alleles
rs28777	92	5	33958916	33959007	A
rs12203592	110	6	396273	396382	C
rs4959270	161	6	457655	457815	AC
rs683	120	9	12709246	12709365	AC
rs1042602	113	11	88911659	88911771	AC
rs1393350	99	11	89010977	89011075	G
rs12821256	119	12	89328278	89328396	CT
rs12896399	73	14	92773627	92773699	G
rs2402130	120	14	92801169	92801288	A
rs1800407	119	15	28230246	28230364	G
N29insA	112	16	89985688	89985799	C
rs1110400	213	16	89985774	89985986	T
rs11547464	213	16	89985774	89985986	G
rs1805005	173	16	89986044	89986216	G
rs1805006	173	16	89986044	89986216	C
rs1805007	213	16	89985774	89985986	C
rs1805008	213	16	89985774	89985986	C
rs1805009	227	16	89986484	89986710	G
rs201326893_Y152OCH	213	16	89985774	89985986	C
rs2228479	173	16	89986044	89986216	G
rs885479	213	16	89985774	89985986	G
rs2378249	118	20	33218028	33218145	A

Table 6 Ancestry Informative SNPs

Locus	Amplicon Length (bp)	Chromosome	Amplicon Start Position	Amplicon End Position	2800M Control Alleles
rs2814778	120	1	159174650	159174769	A
rs3737576	98	1	101709521	101709618	A
rs7554936	106	1	151122413	151122518	CT
rs10497191	101	2	158667153	158667253	C
rs1834619	84	2	17901444	17901527	G
rs1876482	120	2	17362526	17362645	C
rs260690	115	2	109579681	109579795	A
rs3827760	108	2	109513546	109513653	T
rs6754311	98	2	136707920	136708017	CT
rs798443	84	2	7968221	7968304	A
rs12498138	119	3	121459545	121459663	G
rs1919550	117	3	121364112	121364228	A
rs1229984	120	4	100239288	100239407	G
rs3811801	114	4	100244261	100244374	C
rs4833103	95	4	38815462	38815556	AC
rs7657799	116	4	105375396	105375511	T
rs7722456	114	5	170202901	170203014	T
rs870347	119	5	6844995	6845113	T
rs16891982*	108	5	33951621	33951728	G
rs192655	70	6	90518235	90518304	AG
rs3823159	119	6	136482701	136482819	A
rs917115	71	7	28172543	28172613	T
rs1462906	84	8	31896545	31896628	C
rs1871534	71	8	145639652	145639722	C
rs2196051	120	8	122124216	122124335	T
rs6990312	111	8	110602270	110602380	G
rs3814134	104	9	127267664	127267767	T
rs4918664	168	10	94920962	94921129	A
rs1079597	167	11	113296227	113296393	G
rs174570	120	11	61597179	61597298	C
rs2238151	113	12	112211753	112211865	CT
rs671	136	12	112241658	112241793	G
rs1572018	116	13	41715225	41715340	AG
rs2166624	71	13	42579949	42580019	AG
rs7326934	96	13	49070482	49070577	G
rs7997709	85	13	34847693	34847777	T
rs9522149	119	13	111827125	111827243	C
rs200354	165	14	99375246	99375410	G
rs12439433	100	15	36219979	36220078	G
rs1426654	92	15	48426457	48426548	A
rs1800414	116	15	28196969	28197084	A
rs735480	108	15	45152321	45152428	T
rs12913832*	119	15	28365523	28365641	AG
rs459920	78	16	89730800	89730877	T
rs11652805	119	17	62987113	62987231	T
rs17642714	118	17	48726060	48726177	AT



Locus	Amplicon Length (bp)	Chromosome	Amplicon Start Position	Amplicon End Position	2800M Control Alleles
rs2593595	102	17	41056210	41056311	TC
rs4411548	158	17	40658440	40658597	G
rs4471745	67	17	53568849	53568915	G
rs2042762	83	18	35277568	35277650	A
rs3916235	120	18	67578894	67579013	AG
rs4891825	106	18	67867615	67867720	AG
rs7226659	149	18	40488180	40488328	G
rs7251928	200	19	4077044	4077243	A
rs310644	89	20	62159472	62159560	A
rs2024566	88	22	41697312	41697399	A

\* Also used for phenotype prediction.

## Acronyms

Acronym	Definition
2800M	Control DNA 2800M
A7XX	i7 Index Adapter
A50X	i5 Index Adapter
DNL	Diluted Normalized Libraries
DPMA	DNA Primer Mix A
DPMB	DNA Primer Mix B
FEM	Enzyme Mix
FSP	ForenSeq Sample Plate
HP3	2N NaOH
HSC	Human Sequencing Control
HT1	Hybridization Buffer
LNA1	Library Normalization Additives 1
LNB1	Library Normalization Beads 1
LNS2	Library Normalization Storage Buffer 2
LNW1	Library Normalization Wash 1
NLP	Normalized Library Plate
NWP	Normalization Working Plate
PBP	Purification Bead Plate

Acronym	Definition
PCR1	PCR1 Reaction Mix
PCR2	PCR2 Reaction Mix
PLP	Purified Library Plate
PNL	Pooled Normalized Libraries
RSB	Resuspension Buffer
SPB	Sample Purification Beads

## Technical Assistance

For technical assistance, contact Illumina Technical Support.

**Table 7** Illumina General Contact Information

<b>Website</b>	www.illumina.com
<b>Email</b>	techsupport@illumina.com

**Table 8** Illumina Customer Support Telephone Numbers

Region	Contact Number	Region	Contact Number
North America	1.800.809.4566	Italy	800.874909
Australia	1.800.775.688	Netherlands	0800.0223859
Austria	0800.296575	New Zealand	0800.451.650
Belgium	0800.81102	Norway	800.16836
Denmark	80882346	Spain	900.812168
Finland	0800.918363	Sweden	020790181
France	0800.911850	Switzerland	0800.563118
Germany	0800.180.8994	United Kingdom	0800.917.0041
Ireland	1.800.812949	Other countries	+44.1799.534000

**Safety data sheets (SDSs)**—Available on the Illumina website at [support.illumina.com/sds.html](http://support.illumina.com/sds.html).

**Product documentation**—Available for download in PDF from the Illumina website. Go to [support.illumina.com](http://support.illumina.com), select a product, then select **Documentation & Literature**.



Material # 20000923



Illumina  
San Diego, California 92122 U.S.A.  
+1.800.809.ILMN (4566)  
+1.858.202.4566 (outside North America)  
techsupport@illumina.com  
[www.illumina.com](http://www.illumina.com)