

# ForenSeq™ Universal Analysis Software

Interrogate the broadest range of forensically relevant loci with a simple user interface and powerful analysis algorithms.

## Highlights

- Complete, Sample-to-Answer Workflow**  
 Manage the sequencing workflow from run setup to data analysis, data portability, and report generation
- Perform Automated Data Visualization and Reporting with Simple Graphical User Interface**  
 Compare samples automatically, generate population statistics, and view data through an intuitive user interface
- Generate Investigative Leads Through Analysis Software Advances**  
 Estimate visible traits such as hair color and eye color as well as biogeographical ancestry
- Discriminate Further with Comprehensive Sequence Variant Analysis**  
 Increase resolution with tools for intra-STR and flanking region variant detection

## Introduction

The ForenSeq Universal Analysis Software is part of a fully validated<sup>1</sup> sequencing workflow specifically designed for use in forensic genomics applications. The workflow includes the MiSeq FGx™ instrument and ForenSeq DNA Signature Prep Kit and enables simultaneous analysis of approximately 200 forensically relevant single nucleotide polymorphisms (SNPs) and short tandem repeats (STRs)—including marker sets not routinely available with traditional methods (Table 1). From sequencing library<sup>2</sup> and file management to complete data analysis and reporting, the analysis software provides the full functionality needed to analyze and interpret the widest range of casework and database samples.

For remote labs or criminal justice partners that require access to data independent of the MiSeq FGx instrument, the ForenSeq Universal Analysis Software is sold separately on a standalone, dedicated server.

## Complete, Sample-to-Answer Workflow

The Verogen MiSeq FGx Forensic Genomics Solution uses a pipeline of software applications to perform sequencing runs and complete data analysis (Figure 1). ForenSeq Universal Analysis Software guides the sequencing workflow from run-setup and sample data entry to the final stages of data analysis and report generation. MiSeq FGx Control Software (MCS), installed on the MiSeq FGx instrument control computer, captures flow cell images, operates the flow cell stage, and controls reagent delivery and temperature. During the run, Real-Time Analysis (RTA) software performs image analysis, base calling, and




Figure 1: ForenSeq Universal Analysis Software Workflow

Table 1: Simultaneous Analysis of Forensically Relevant Loci with ForenSeq Universal Analysis Software

Feature	Markers <sup>a</sup>
Global Autosomal STRs	27
Y-STRs	24
X-STRs	7
Identity SNPs	94
Phenotypic SNPs	22
Biogeographical Ancestry SNPs	56
Total Number of Loci <sup>b</sup>	> 200

- a. SNP and STR chromosome locations can be found in the ForenSeq DNA Signature Prep Reference Guide.
- b. Over 200 markers analyzed when running primer set B.

assigns base-by-base quality scores. ForenSeq Universal Analysis Software then initiates the final stages of analysis, including demultiplexing, sequence alignment, allele calling, genotyping, and reporting.

## Manage Data with Data Portability Features and Sample History

The ForenSeq Universal Analysis Software allows creation of data packages for easy management and portability of sample information. Data can be packaged specifically for archive to a separate data storage location, or copied for discovery requests and data sharing between labs. Document and view both system and user actions within the ForenSeq Universal Analysis Software by enabling visibility of the Sample History.

\* A sequencing "library" is a collection of amplified DNA fragments from a single DNA sample.

The screenshot displays the ForenSeq™ software interface. At the top, the logo and user information 'user@yourlab.com' and 'VEROGEN' are visible. The main area is titled 'Run 13 Project' and includes a search bar for 'Project Name'. On the left, there are two summary cards for runs: '2.0 - Run 16' (Completed 11 Apr 2018, 3 Samples) and '1.0 - Run 16' (Completed 07 Mar 2018, 4 Samples). The main content area shows a 'GENERATE REPORT(S)' button and a table of reports organized by project level, sample genotype, comparisons, and phenotypes. The table lists reports for Sample 1, Sample 6, and Positive Control, with columns for Sample Name, Analysis Name, Report Type, and Date Created. A '1 - 6 of 6 Items' indicator is at the bottom right of the table.

Sample Name	Analysis Name	Report Type	Date Created
Sample 1	1.0 - Run 16	Sample Details	05 Apr 2018
Sample 1	1.0 - Run 16	Sample Summary	05 Apr 2018
Sample 6	1.0 - Run 16	Sample Details	05 Apr 2018
Positive Control	1.0 - Run 16	CODIS	08 Mar 2018
Positive Control	1.0 - Run 16	CODIS	08 Mar 2018
Sample 1	1.0 - Run 16	CODIS	08 Mar 2018

**Figure 2: ForenSeq Universal Analysis Project Reports**—Reports are organized by project, sample, comparisons, and phenotypes. Project level reports include CODIS, Project Details, and Flanking Region reports. Sample Genotypes provide all allele information for a specific sample. Comparisons and Phenotypic reports provide their relevant information by sample.

## Perform Automated Data Visualization and Reporting with Simple Graphical User Interface

The Verogen ForenSeq Universal Analysis Software provides an easy data analysis workflow and user-friendly viewing and reporting features. The graphical user interface has a simple, intuitive design that can be viewed from any laboratory computer through the Google Chrome web browser. The Run Overview screen features easy run-setup, sample information, and index tracking options. Within the analysis software, sequencing run data can be associated with specific, user-defined projects and results can be viewed by run or by project in the Projects screen.

At run completion, ForenSeq Universal Analysis Software offers data visualization through several types of display screens. The Sample Details screen provides a summarized view of SNP and STR locus intensity scatter plots. Below the intensity charts, individual allele calls are displayed, along with convenient icons displayed when specific quality flags are triggered (Figure 3). To view more in-depth reports on specific loci, drill down using the Locus Detail pop-up screens. The Locus Detail screens display intensity charts, the base-by-base DNA target sequence, stutter, and more (Figure 4). The analysis software also offers a full suite of run metrics and sample quality controls that can be evaluated after the run. Quality metrics and sample quality controls are displayed with color indicators allowing forensic analysts to quickly scan the results and determine whether controls fall within the recommended ranges.

Project Detail reports and Sample Detail reports can be generated automatically within the software (Figure 2). Optional offline flanking region analysis provides additional sequence variant information from select amplicons for even more discrimination power. CODIS-compatible sample-level and project-level reports can be generated for upload to criminal databases. Each report can be easily printed or

exported in .xlsx format. Password protected user accounts ensure secure access to the MiSeq FGx sequencing instrument and to the analysis server.

## Generate Investigative Leads through Analysis Software Advances

Compared to current capillary electrophoresis-based methods, ForenSeq Universal Analysis Software provides a number of data analysis advantages. These include, the capacity to interrogate a greater number of markers (Table 1) and the ability to recover the maximum amount of useful genomic information from degraded DNA, low quality DNA, or complex mixtures<sup>2</sup>. Furthermore, all markers, including SNPs, autosomal, X-, and Y-STRs, can be analyzed simultaneously using a single, streamlined workflow.

Beyond analysis of autosomal, X-, and Y-STRs, the software enables the analysis of marker sets not routinely available with traditional capillary electrophoresis methods. These include a dense set of identity-informative SNPs (iiSNPs)<sup>3,4</sup>, which are informative for source attribution, phenotypic-informative SNPs (piSNPs)<sup>5</sup>, which provide estimates of eye color (blue, intermediate, brown) and hair color (brown, red, black, blond), and biogeographical ancestry-informative SNPs (aiSNPs)<sup>6</sup>. Biogeographical ancestry estimation is presented as a principal component analysis (PCA) plot relative to major population groups<sup>7</sup> (Figure 5). Centroids within the plots then provide logical groupings that give meaningful context to the unknown sample.

Additionally, ForenSeq Universal Analysis Software can detect intra-STR sequence variation that discloses additional identity informative data beyond length variation (Figure 4). These data can be important in kinship analyses and mixture analysis, where DNA contributors share equi-length alleles that can only be distinguished by comparing their base-by-base sequences.




Figure 3: ForenSeq Universal Analysis Software Sample Details Screen—The Sample Details screen displays overall sample results as an intensity plot, along with detailed summary tables for all loci included in the sequencing library.




Figure 4: ForenSeq Universal Analysis Software Locus Details Screen—The pop-up Locus Detail screen shows the allele call, intensity bar graph, and stutter. View the Repeat Sequence column to see the full, base-by-base target sequence, and any intra-STR variation present in the sample.




Figure 5: ForenSeq Universal Analysis Software Estimation of Visible Traits—When the appropriate primer mix is used, the software offers an optional estimation of visible traits including eye and hair color. Markers for biogeographical ancestry can also be tested and reviewed.

**Table 2: ForenSeq Server**

Components and Specifications
4 × Seagate 1 TB SATA hard drives
RAID controller allowing for data redundancy and speed
Intel 2.4 GHz × 64 processor with 6 cores/12 threads
Intel socket R3 server board
32 GB DDR4 RAM
550 W power supply
ULL, FCC, CE certified
Windows Server 2012 R2 Standard with 5 CALs

## Ordering Information

Product	Catalog No.
ForenSeq Universal Analysis Software and Server	SE-550-1001

## Summary

ForenSeq Universal Analysis Software contains comprehensive sample management and analytical capabilities, including sample and index management, application-specific workflows, data visualization at sample and locus levels, quality flags to simplify data analysis, and easily exportable reports.

## Learn More

To learn more about ForenSeq Universal Analysis Software, visit: [www.verogen.com/products](http://www.verogen.com/products)

## References

- Jäger AC, Alvarez ML, Davis CP, et al. Developmental validation of the MiSeq FGx Forensic Genomics System for targeted next generation sequencing in forensic DNA casework and database laboratories. *Forensic Sci Int.* 2017;28: 52-70.
- Verogen (2018) ForenSeq DNA Signature Prep Kit Data Sheet ([www.verogen.com/products](http://www.verogen.com/products)).
- Kidd KK, Pakstis AJ, Speed WC, Grigorenko EL, Kajuna SL, et al. (2006) Developing a SNP panel for forensic identification of individuals. *Forensic Sci Int* 164(1): 20–32.
- Sanchez JJ, Phillips C, Børsting C, Balogh K, Bogus M, et al. (2006) A multiplex assay with 52 single nucleotide polymorphisms for human identification. *Electrophoresis* 27(9): 1713–1724.
- Walsh S, Liu F, Wollstein A, Kovatsi L, Ralf A, et al. (2013) The HirisPlex system for simultaneous prediction of hair and eye colour from DNA. *Forensic Sci Int Genet* 7(1): 98–115.
- Kidd KK, Speed WC, Pakstis AJ, Furtado MR, Fang R, et al. (2013) Progress toward an efficient panel of SNPs for ancestry inference. *Forensic Sci Int Genet* 10: 23–32.
- Phillips C, Prieto L, Fondevila M, Salas A, Gómez-Tato A, et al. (2009) Ancestry analysis in the 11-M Madrid bomb attack investigation. *PLoS One* 4(8): e6583.