

## Amplify and Tag Targets

### Purified DNA

1. Using nuclease-free water, dilute 1 ng purified DNA input to 0.2 ng/μl.
2. In the Master Mix tube, combine the following volumes per sample:
  - PCR1 (4.7 μl)
  - FEM (0.3 μl)
  - DPMA or DPMB (5 μl)
3. Pipette and centrifuge briefly.
4. [Optional] Distribute master mix among an 8-tube strip.
5. Add 10 μl master mix to the FSP.
6. In a 1.7 ml tube, combine:
  - 2800M (2 μl)
  - Nuclease-free water (98 μl)
7. Invert three times and centrifuge briefly.
8. Add 5 μl diluted 2800M to the FSP.
9. Add 5 μl nuclease-free water to the FSP.
10. Add 5 μl 0.2 ng/μl DNA to the FSP and pipette.
11. Centrifuge at 1000 × g for 30 seconds.
12. Place on the thermal cycler and run the PCR1 program.

### SAFE STOPPING POINT

If you are stopping, store the plate at 2°C to 8°C for up to 2 days. Alternatively, leave on the thermal cycler overnight.

### Crude Lysate

1. In the Master Mix tube, combine the following volumes per sample:
  - PCR1 (4.7 μl)

- FEM (0.3 μl)
  - DPMA or DPMB (5 μl)
  - Nuclease-free water (3 μl)
2. Pipette and centrifuge briefly.
  3. [Optional] Distribute master mix among an 8-tube strip.
  4. Add 13 μl master mix to the FSP.
  5. In a 1.7 ml tube, combine:
    - 2800M (2 μl)
    - Nuclease-free water (38 μl)
  6. Invert three times and centrifuge briefly.
  7. Add 2 μl diluted 2800M to the FSP.
  8. Add 2 μl nuclease-free water to the FSP.
  9. Add 2 μl diluted crude lysate to the FSP.
  10. Centrifuge at 1000 × g for 30 seconds.
  11. Place on the thermal cycler and run the PCR1 program.

### SAFE STOPPING POINT

If you are stopping, store the plate at 2°C to 8°C for up to 2 days. Alternatively, leave on the thermal cycler overnight.

### FTA Card

1. Place a 1.2 mm FTA card punch into each FSP well.
2. Add 100 μl 1X TBE Buffer.
3. Place on a PCR tube storage rack.
4. Shake at 1800 rpm for 2 minutes.
5. Centrifuge at 1000 × g for 30 seconds.
6. Remove and discard all supernatant.
7. Add the following volumes to each control well:
  - PCR1 (4.7 μl)
  - FEM (0.3 μl)
  - DPMA or DPMB (5 μl)
8. In a 1.7 ml tube, combine:
  - 2800M (2 μl)

- Nuclease-free water (98 μl)
9. Invert three times and centrifuge briefly.
  10. Add 5 μl diluted 2800M to the positive control wells.
  11. Pipette to mix.
  12. Add 5 μl nuclease-free water to the negative control wells.
  13. Pipette to mix.
  14. In the FTA Master Mix tube, combine the following volumes per sample:
    - PCR1 (4.7 μl)
    - FEM (0.3 μl)
    - DPMA or DPMB (5 μl)
    - Nuclease-free water (5 μl)
  15. Pipette, and then centrifuge briefly.
  16. [Optional] Distribute FTA Master Mix among an 8-tube strip.
  17. Add 15 μl FTA Master Mix to the FSP.
  18. Centrifuge at 1000 × g for 30 seconds.
  19. Place on the thermal cycler and run the PCR1 program.

### SAFE STOPPING POINT

If you are stopping, store the plate at 2°C to 8°C for up to 2 days. Alternatively, leave on the thermal cycler overnight.

## Enrich Targets

1. Centrifuge the FSP at 1000 × g for 30 seconds.
2. Place the Index 1 Adapter tubes in the ForenSeq Index Plate Fixture.
3. Place the Index 2 Adapter tubes in the ForenSeq Index Plate Fixture.
4. Place the FSP on the ForenSeq Index Plate Fixture.
5. Add index adapters:
  - a. Add 4 μl R7XX down each column.
  - b. Replace the cap on each Index 1 Adapter tube.

# ForenSeq DNA Signature Prep Kit Checklist



- c. Add 4  $\mu$ l A50X across each row.
- d. Replace the cap on each Index 2 Adapter tube.
- 6. Vortex PCR2, and then centrifuge briefly.
- 7. [Optional] Distribute PCR2 among an 8-tube strip.
- 8. Add 27  $\mu$ l PCR2.
- 9. Centrifuge at 1000  $\times$  g for 30 seconds.
- 10. Place on the thermal cycler and run the PCR2 program.

## SAFE STOPPING POINT

If you are stopping, store the plate at 2°C to 8°C for up to 7 days. Alternatively, leave on the thermal cycler overnight.

## Purify Libraries

- 1. Add 45  $\mu$ l SPB to the PBP.
- 2. Centrifuge the FSP at 1000  $\times$  g for 30 seconds.
- 3. Transfer 45  $\mu$ l reaction from the FSP to the PBP.
- 4. Shake at 1800 rpm for 2 minutes.
- 5. Incubate at room temperature for 5 minutes.
- 6. Place on the magnetic stand until clear.
- 7. Remove and discard all supernatant.
- 8. Wash as follows.
  - a. Add 200  $\mu$ l fresh 80% EtOH.
  - b. Incubate for 30 seconds.
  - c. Remove and discard all supernatant.
- 9. Wash a **second** time.
- 10. Centrifuge at 1000  $\times$  g for 30 seconds.
- 11. Place on the magnetic stand.
- 12. Remove residual EtOH.
- 13. Remove from the magnetic stand.
- 14. Add 52.5  $\mu$ l RSB.
- 15. Shake at 1800 rpm for 2 minutes.
- 16. If necessary, pipette or reshake.
- 17. Incubate at room temperature for 2 minutes.

- 18. Place on the magnetic stand until clear.
- 19. Transfer 50  $\mu$ l supernatant from the PBP to the PLP.
- 20. Centrifuge at 1000  $\times$  g for 30 seconds.

## SAFE STOPPING POINT

If you are stopping, store the plate at -25°C to -15°C for up to 1 year.

## Normalize Libraries

- 1. In the LNA1/LNB1 Master Mix tube, combine the following volumes per sample:
  - LNA1 (46.8  $\mu$ l)
  - LNB1 (8.5  $\mu$ l)
- 2. Vortex and invert several times.
- 3. Transfer to a reagent reservoir.
- 4. Add 45  $\mu$ l LNA1/LNB1 Master Mix to the NWP.
- 5. To clear any aspirated beads, place the PLP on the magnetic stand until clear.
- 6. Transfer 20  $\mu$ l supernatant from the PLP to the NWP.
- 7. Shake at 1800 rpm for 30 minutes.
- 8. While shaking, perform steps 9–11.
- 9. In the 0.1 N HP3 tube, combine the following volumes per sample:
  - Nuclease-free water (33.3  $\mu$ l)
  - HP3 (1.8  $\mu$ l)
- 10. Invert several times and set aside.
- 11. Add 30  $\mu$ l LNS2 to the NLP.
- 12. Immediately after shaking, place the NWP on the magnetic stand until clear.
- 13. Remove and discard all supernatant.
- 14. Remove from the magnetic stand.
- 15. Wash as follows.
  - a. Add 45  $\mu$ l LNW1 to each well.
  - b. Shake at 1800 rpm for 5 minutes.
  - c. Place on the magnetic stand until clear.

- d. Remove and discard all supernatant.
- e. Remove from the magnetic stand.

- 16. Wash a **second** time.
- 17. Centrifuge at 1000  $\times$  g for 30 seconds.
- 18. Place on the magnetic stand until clear.
- 19. Remove residual LNW1.
- 20. Remove from the magnetic stand.
- 21. Add 32  $\mu$ l 0.1 N HP3.
- 22. Shake at 1800 rpm for 5 minutes.
- 23. Place on the magnetic stand until clear.
- 24. Transfer 30  $\mu$ l supernatant from the NWP to the NLP.
- 25. Pipette to mix.
- 26. Centrifuge at 1000  $\times$  g for 30 seconds.

## SAFE STOPPING POINT

If you are stopping, store the plate at -25°C to -15°C for up to 30 days.

## Pool Libraries

- 1. Transfer 5  $\mu$ l of each library to an 8-tube strip.
- 2. Store the NLP at -25°C to -15°C for up to 30 days.
- 3. Transfer libraries from the 8-tube strip to the PNL tube.
- 4. Vortex and centrifuge briefly.

## SAFE STOPPING POINT

If you are stopping, store the tube at -25°C to -15°C for up to 30 days.

## Denature and Dilute Libraries

- 1. In the Denatured HSC tube, combine:
  - HSC (2  $\mu$ l)
  - HP3 (2  $\mu$ l)
  - Nuclease-free water (36  $\mu$ l)

# ForenSeq DNA Signature Prep Kit Checklist

2. Pipette gently to mix. Cap and centrifuge briefly to mix.
3. Incubate at room temperature for 5 minutes.
4. Add 591 µl HT1 to the DNL tube.
5. Transfer 7 µl library from the PNL tube to the DNL tube.
6. Pipette to mix.
7. Store the PNL tube at -25°C to -15°C for ≤ 30 days.
8. Add 4 µl denatured HSC to the DNL tube.
9. Pipette to mix.
10. Vortex and centrifuge briefly.
11. Place on the microheating system for 2 minutes.
12. Invert several times to mix.
13. Immediately cool for 5 minutes.
14. Transfer entire volume to the reagent cartridge.

## Acronyms

Acronym	Definition
2800M	Control DNA 2800M
A50X	Index 2 (i5) Index Adapter
DNL	Diluted Normalized Libraries
DPMA	DNA Primer Mix A
DPMB	DNA Primer Mix B
FEM	Enzyme Mix
FSP	ForenSeq Sample Plate
HP3	2 N NaOH
HSC	Human Sequencing Control
HT1	Hybridization Buffer

Acronym	Definition
LNA1	Library Normalization Additives 1
LNB1	Library Normalization Beads 1
LNS2	Library Normalization Storage Buffer 2
LNW1	Library Normalization Wash 1
NLP	Normalized Library Plate
NWP	Normalization Working Plate
PBP	Purification Bead Plate
PCR1	PCR1 Reaction Mix
PCR2	PCR2 Reaction Mix
PLP	Purified Library Plate
PNL	Pooled Normalized Libraries
R7XX	Index 1 (i7) Index Adapter
RSB	Resuspension Buffer
SPB	Sample Purification Beads