

Universal Analysis Software mtDNA Module

Version 2

Reference Guide

VEROGEN PROPRIETARY

Document # VD2022003 Rev. A

June 2022

Revision History

Document #	Date	Description of Change
VD2022003 Rev. A	June 2022	Initial release of the Universal Analysis Software – mtDNA Module
VD2019002 Rev. E	October 2021	<ul style="list-style-type: none"> • Updated Dark Mode to include color deficiency accessibility • Updated Icons for Run Status and added new status for Analysis Status • Added instructions for creating Flanking Regions report and the report format for ForenSeq MainstAY analysis • Added ability to find SW version number in system settings • Updated calculation for total Alleles in MainstAY sample summary, for female samples • Included instructions to identify software version number
VD2019002 Rev. D	May 2021	<p>Updated software descriptions to Universal Analysis Software v2.3, which supports the ForenSeq MainstAY Kit with NA24385 as the positive control and refreshes the user interface.</p> <ul style="list-style-type: none"> • Added the ForenSeq MainstAY assay, Analysis Method, and STR reports. • Added run metrics for ForenSeq MainstAY libraries. • Added instructions for performing a status check for troubleshooting. • Updated instructions for creating CODIS reports, including the addition of options for batch IDs, case IDs, kit names, and partial profiles. • Updated ForenSeq Positive Amplification Control options for ForenSeq mtDNA Control Region and ForenSeq mtDNA Whole Genome to HL60 and ForenSeq Positive Amplification Control options for ForenSeq Kintelligence and ForenSeq MainstAY to NA24385 <p>Updated the safety data sheet (SDS) link to verogen.com/product-documentation. Renamed the software to Universal Analysis Software and this guide to <i>Universal Analysis Software v2.0 Reference Guide</i>. Added the term stutter to the glossary. Corrected target run metrics for ForenSeq Kintelligence libraries. Replaced kinship with Kintelligence.</p>

Universal Analysis Software – mtDNA Module Reference Guide

Document #	Date	Description of Change
VD2019002 Rev. C	February 2021	<p>Updated software descriptions to ForenSeq Universal Analysis Software v2.2, which supports the ForenSeq Kintelligence Kit with NA24385 as the positive control.</p> <ul style="list-style-type: none">• Added the ForenSeq Kintelligence assay, Verogen Kintelligence Analysis Method, and Kintelligence reports.• Added run metrics for ForenSeq Kintelligence libraries.• Added instructions for assigning Unique Dual Index adapters, accessing the analysis folder, switching analysis methods, and printing. <p>Updated the ForenSeq trademark to a registered trademark.</p> <p>Refreshed the format of the guide:</p> <ul style="list-style-type: none">• Updated the layout and formatting.• Updated fonts, table styles, and other design elements. <p>Added guidance for failed run metrics.</p> <p>Added a link to the Verogen Documentation web page.</p> <p>Added the glossary terms clusters, intensity, intralocus balance, positions called, reference SNP ID number, short tandem repeat, typed and untyped, and Unique Dual Index.</p> <p>Added ambient temperature and reagent chiller temperature as potential sources of high phasing or prephasing values.</p> <p>Reorganized information on primary workspaces, reports, analysis methods, and reviewing sample data.</p> <p>Renamed the Sample Details page to Results and updated other user interface terminology as needed.</p> <p>Clarified that Dark Mode is the only software setting available to non-administrators.</p> <p>Corrected flow cell compatibility for the ForenSeq mtDNA Control Region Kit.</p>

Universal Analysis Software – mtDNA Module Reference Guide

Document #	Date	Description of Change
VD2019002 Rev. B	August 2020	<p>Updated software descriptions for ForenSeq Universal Analysis Software v2.1, which provides the following enhancements:</p> <ul style="list-style-type: none">• Support for the ForenSeq mtDNA Whole Genome Kit• A Run, Project, and Sample History setting in System Settings• Project History and Sample History worksheets for Project and Sample Reports <p>Updated guide format and organization, including adding an index and consolidating workspace descriptions.</p> <p>Revised step-by-step instructions to be more succinct and include field descriptions.</p> <p>Updated telephone numbers for Verogen contact information.</p> <p>Updated glossary term definitions and added the terms coverage plot, indel, library, mtDNA navigator, plexity, position viewer, system event, and user action.</p> <p>Updated figures to include callouts highlighting user interface elements.</p> <p>Added information on strand bias and lost passwords.</p> <p>Added instructions for creating projects, assigning a sample to multiple projects, modifying user-created analysis methods, and enabling and disabling dark mode.</p> <p>Added another possible cause of low cluster density: belatedly combining bead-based normalized libraries and Hybridization Buffer (HT1).</p> <p>Noted that the analysis folder contains VCF and BAM files.</p> <p>Noted that the software follows SWGDAM nomenclature rules.</p> <p>Noted the conditions for displaying the Call or No Call filter.</p> <p>Clarified that each sample name must be unique in the software.</p> <p>Clarified that phasing and prephasing values are for Read 1 only.</p>
VD2019002 Rev. B	August 2020	<p>Distributed troubleshooting information among individual sections.</p> <p>Moved information on analysis method settings to <i>Analysis Defaults for ForenSeq Universal Analysis Software v2 (Pub. No. VD20200045)</i>.</p> <p>Corrected duplicate entries of D7S820 and D16S539 in the Human Sequencing Control.</p> <p>Corrected the list of information included in project and sample reports.</p> <p>Corrected citations that appear in the appendix.</p>
VD2019002 Rev. A	August 2019	Initial release

Table of Contents

Revision History	2
Table of Contents	5
Chapter 1 Software Overview	8
Introduction	8
Access the Software	9
Primary Workspaces	10
Software Status	16
Glossary	16
Chapter 2 Run Setup	21
Create a Run	21
Edit a Run	21
Add Samples to a Run	21
Manage Samples in a Run	23
Monitor a Run	25
Chapter 3 Data Analysis	26
Review Run Metrics	26
Review Results	27
Add Analyzed Samples to a Project	29
Reanalyze Samples	29
Chapter 4 Reports	31
Generating Reports	31
File Formats for Reports	35
Chapter 5 Run, Project, and Sample Management	38
Create a Project	38
Edit a Project	38
Obtain the Analysis Folder Path	39
Printing	39
History	39

Chapter 6 Software Settings	43
Settings Overview	43
Analysis Methods	43
CODIS Report Defaults	46
User Management	47
System Settings	47
Enable or Disable Dark Mode	48
Display Software Version	48
Chapter 7 Troubleshooting	49
Troubleshoot a Run	49
Perform a Status Check	53
Appendix A Run Metrics	54
Accessing Run Metrics	54
Quality Metrics	54
Positive and Negative Controls	58
Sample Representation	60
Appendix B Supporting Information	62
Special Handling of Sequencing Data	62
Index	66
Technical Support	69

Universal Analysis Software – mtDNA Module Reference Guide

This document and its contents are proprietary to Verogen, Inc. and its affiliates, and are intended solely for the contractual use of its customer in connection with the use of the product(s) described herein and for no other purpose. This document and its contents shall not be used or distributed for any other purpose and/or otherwise communicated, disclosed, or reproduced in any way whatsoever without the prior written consent of Verogen. Verogen does not convey any license under its patent, trademark, copyright, or common-law rights nor similar rights of any third parties by this document.

The instructions in this document must be strictly and explicitly followed by qualified and properly trained personnel in order to ensure the proper and safe use of the product(s) described herein. All of the contents of this document must be fully read and understood prior to using such product (s).

FAILURE TO COMPLETELY READ AND EXPLICITLY FOLLOW ALL OF THE INSTRUCTIONS CONTAINED HEREIN MAY RESULT IN DAMAGE TO THE PRODUCT(S), INJURY TO PERSONS, INCLUDING TO USERS OR OTHERS, AND DAMAGE TO OTHER PROPERTY.

VEROGEN DOES NOT ASSUME ANY LIABILITY ARISING OUT OF THE IMPROPER USE OF THE PRODUCT(S) DESCRIBED HEREIN (INCLUDING PARTS THEREOF OR SOFTWARE).

© 2022 Verogen, Inc. All rights reserved.

All other names, logos, and other trademarks are the property of their respective owners.

For Research, Forensic, or Paternity Use Only. Not for any clinical or therapeutic use in humans or animals.

Software Overview

Introduction	8
Access the Software	9
Primary Workspaces	10
Software Status	16
Glossary	16

Introduction


Universal Analysis Software (UAS v2) is a forensic genomics solution for run setup, sample management, data analysis, and reporting. Pre-installed on a dedicated server, the software integrates with the MiSeq FGx[®] Sequencing System to automatically generate FASTQ files and call variants from run data.

This reference guide provides details on the analysis of data generated from the ForenSeq mtDNA Whole Genome and mtDNA Control Region kits using the mtDNA analysis method, as well as data generated from third-party kits sequenced in FGx mode.

Analysis Workflow

The following diagram outlines the analysis workflow, which includes UAS and two applications that run on the MiSeq FGx System: MiSeq FGx Control Software (MFCS) and Real-Time Analysis (RTA). For more information on these applications, see the *MiSeq FGx Sequencing System Reference Guide (document # VD2018006)*.

Figure 1 Overview of data analysis steps



 UAS v2

 MiSeq FGx System

Access the Software

UAS is a browser-based application that does not require an internet connection. Open a web browser (Chrome is recommended) on a computer with access to the same network as the MiSeq FGx System. Set screen resolution to at least 992 pixels for optimal display.

To access the UAS user interface, enter the server address in the web browser address bar. If the Domain Name System (DNS) does not list the server, enter the server IP address for the network.

Create an Account

Before using UAS, you must create an account. The first user to create an account is automatically assigned administrator privileges. The administrator is responsible for approving new accounts. For instructions, see [Enable and Disable Accounts \(page 47\)](#).

Lost passwords are not recoverable, so creating a second administrator account is highly recommended. If the administrator password is lost and only one administrator account exists, the software must be uninstalled and reinstalled.

Universal Analysis Software – mtDNA Module Reference Guide

1. Access the network server through a web browser.
2. On the Login page, select **Register as a new user**.
3. On the Register page, complete the following fields:
 - **Email** — Enter an email address to be the username for the account.
 - **Password** — Enter a password for the account.
 - **Confirm password** — Reenter the password for the account.
4. Select **Register** to submit your account for approval from the system administrator. A confirmation message appears. You can access the software after approval.

Sign In or Sign Out

1. Sign in as follows.
 1. **a.** Access the network server through a web browser.
 2. **b.** Enter your email address (username) and password.
 3. **c.** Select **Login**.

The software opens to the Home page.

2. When ready to sign out, select  **Account**, and then select **Logout**.

Change a Password

1. On the navigation rail, select  **Account**.
2. Select **Change Password**.
3. In the Current Password field, enter your password.
4. In the New Password field, enter a new password.
5. In the Confirm Password field, repeat the new password.
6. Select **Save** to apply the new password.

Primary Workspaces

The following pages comprise the primary workspaces in UAS:

- **Home** — View and create runs and projects. This page functions as a system dashboard.
- **Runs** — Set up and manage runs before and during a run.
- **Projects** — View and manage projects and samples after a run and generate reports.
- **Results** — Review data on a sample-by-sample basis after a run.

The software also includes printing, commenting, and other functions that supplement run setup, data analysis, and reporting. For instructions on using these supplementary features, see [Chapter 5 Run, Project, and Sample Management](#).


Home Page and Navigation

The Home page displays information and functions organized into runs and projects. Selecting a run opens the Run Details page, and selecting a project opens the Project Details page. Color-coded Quality Metrics icons provide an overall view of run performance. Hovering over an icon displays more detailed metrics.

Universal Analysis Software – mtDNA Module Reference Guide

The navigation rail vertically arranges icons and puts key destinations within easy reach. It appears on the left side of each page so you can go to runs, projects, settings, and your account or can return to the Home page.

Figure 2 Overview of the Home page



- A** Navigation rail — Icons that go to key destinations
- B** Runs — List of runs with quality metrics and options to search and create
- C** Projects — List of projects with options to search and create
- D** Quality Metrics icon — Quality metrics for a completed run
- E** Scroll arrows — Move through all runs or projects

Runs


Runs is the sequencing workspace for creating and monitoring runs and the input and organization of sample information. Selecting Runs on the navigation rail opens the Runs navigation drawer, which provides functions for creating and searching runs.

The Runs navigation drawer also lists every run in the software with the run state, number of samples in the run, flow cell type, and color-coded Quality Metrics icon. Runs are listed from newest to oldest. Use the scroll arrows to move through the complete list.

Run Details Page

Selecting a run from the Runs navigation drawer opens the Run Details page. Each Run Details page is a run-specific workspace for managing the run and samples in the run. Throughout the software, selecting a run name hyperlink goes to the Run Details page.

Figure 3 Run Details Page on the Universal Analysis Software user interface



A table displays all samples in the run with information about each sample and options for sorting and filtering. Functions for importing, adding, searching, and deleting samples appear above the table. The header section of a Run Details page displays the run name, quality metrics, sample representation, and the following icons.

Table 1 Run Details page icons

Icon	Name	Function
	Run Details	Edit the run and view the run state, cycles per read, and other run information.
	Run Activity	See system events and user actions for the run and add comments. This icon is visible when the history setting is enabled. See System Settings (page 47) .
	Print	Print the Run Details page.

Projects

Projects is the post-sequencing workspace for creating and managing projects, accessing analysis results, and generating reports. Selecting Projects on the navigation rail opens the Projects navigation drawer, which provides functions for creating and searching projects.

The Projects navigation drawer also lists every project in the software with the number of samples and runs in the project and the date and time of project creation. Projects are listed from newest to oldest. Use the scroll arrows to move through the complete list.

Project Side Sheet

Selecting a project from the Projects navigation drawer opens the Project side sheet and the Samples page. The Project side sheet lists all samples in the project, grouped by run. Expand the runs to see which samples are included. A green checkmark indicates an analyzed sample.

Universal Analysis Software – mtDNA Module Reference Guide

The Project side sheet provides the following icons for accessing additional pages and project information.

Table 2 Project side sheet icons

Icon	Name	Function
	Edit	Change the project name and add or change the optional description.
	Samples	Open the Samples page for the project.
	Reports	Open the Reports page for the project.
	Project Activity	See system events and user actions for the project and add comments. This icon is visible when the history setting is enabled. See System Settings (page 47) .
	Print	Print the Samples page.
	QPN	Open the QPN page for each run in the project. Each icon letter is color-coded to provide high-level status. See Appendix A Run Metrics .

Samples Page

The Samples page lists information for each sample in a selected project with functions for adding, reanalyzing, and searching samples. Filters narrow results by sample or assay type and sort by sample name, run name, and analysis method.

A QC Indicators column displays the number of quality control (QC) indicators for a sample. The column displays 0 for a sample with zero QC indicators. Samples with QC indicators may require further review.

Reports Page

The Reports page lists reports in a project by creation date, starting with the most recent report. On this page, you can create reports or find and download existing reports. Sort by the Date Created and File Name columns and filter the Report Type column to organize results.

QPN Page

The QPN page is divided into three tabs: Quality Metrics, Positive and Negative Controls, and Sample Representation. Together, these tabs provide an overall view of run performance. For details, see [Appendix A Run Metrics](#).

Results Page

Selecting an mtDNA sample from a project opens the Results page. Each Results page is a sample-specific workspace for in-depth analysis of results. A header section displays the sample name, analysis method applied to the sample, buttons to create reports and filter and compare samples, and the following icons.

Table 3 Results page icons


Icon	Name	Function
	Sample Details	View the sample type and analysis state, analysis folder path, and other analysis information.
	Sample Activity	See system events and user actions for the sample and add comments. This icon is visible when the history setting is enabled. See System Settings (page 47) .

The remainder of the Results page is divided into an mtDNA navigator, position viewer, and coverage plot that examine each position. Changing the position in one of these tools simultaneously updates the position in the others.

mtDNA Navigator

The mtDNA navigator displays data spanning the whole mtGenome or control region. Data include the IT, variants, calls or no calls, and coverage at various positions. The movable pointer updates the position.

Figure 4 mtDNA navigators for whole genome (left) and control region (right) results



- A** Maximum coverage indicator
- B** Movable pointer
- C** Control region view (filter applied)

By default, the mtDNA navigator shows the entire mtGenome. To focus on the control region only, apply the Control Region filter.

Position Viewer

The position viewer compares the sample sequence to the revised Cambridge Reference Sequence (rCRS). Color-coded dashes indicate ITs, variants, and no calls. For example, the blue dash in [Figure 5](#) underlines the base at position 16,362 because the call (C) differs from the reference (T).


Individually select or enter each position number, or use the scroll arrows to move through pages of sequence segments and hotspots. For each selected position, the position viewer displays the following data:

- **Reference and call** — The base calls for the reference and sample.
- **Total read count** — The total number of reads from the forward and reverse strands.

Universal Analysis Software – mtDNA Module Reference Guide

- **Strand depth** — The strand (forward or reverse) read count with the majority of reads.
- **Percentages** — The percentage of each base compared to the AT and IT.

Figure 5 Position viewer set to position 16,362



- A Selected position number
- B Scroll arrows
- C Data for the selected position
- D Sample sequence
- E Reference sequence

Strand Bias

If you want to evaluate strand bias, which provides an additional metric for evaluating quality, use the strand depth and total read count. The software reports strand depth as the read count from the forward or reverse strand that has the majority of reads for a base position. Subtract from the total read to determine the read count for the opposite strand.

Coverage Plot

The coverage plot displays the coverage (number of reads) for each mtDNA position. Positions are plotted along the x-axis and the y-axis shows coverage. A pointer marks the selected position. Move the pointer to shift the position.

Figure 6 Coverage at each position



- A Coverage along the y-axis
- B Positions along the x-axis
- C Movable pointer

Software Status

The software displays a state for runs and analyses. The state determines which modifications are permitted and are also useful for general status updates.

Run State

After a run is created and saved, you can view the run state by selecting the Information icon on the Run Details page displays the run state. The following table describes each possible state.

Table 4 Run states

State	Description
Created	The run information is saved in UAS.
Sequencing	The run is in process on the MiSeq FGx System.
Completed	The run completed successfully.
Errored	An error prevented the run from finishing, or it was manually stopped.
Paused	A user suspended the run.
Stopped	A user ended the run due to incorrect setup, poor data quality, or other issue.

Analysis State

After a sample is added to a run, you can view the analysis state by selecting the Information icon on the Results page. The following table describes each possible state.

Table 5 Analysis states

Symbol	Retain Status	Description
	Created	UAS has not started analysis.
	Queued	The analysis is in line to start.
	Processing	Analysis or reanalysis is in progress.
	Completed	Analysis has successfully completed.
	Errored	Analysis has unsuccessfully completed.

Glossary

Term	Description
Amplicon	The product of PCR amplification of a targeted region of interest in a genome.
Amplicon end position	The last base position of an amplicon, including the PCR primer binding sequence.

Universal Analysis Software – mtDNA Module Reference Guide

Term	Description
Amplicon start position	The first base position of an amplicon, including the PCR primer binding sequence.
Analysis method	Settings and thresholds that inform the analysis of sequencing data to generate variant calls and QC indicators.
Analytical threshold (AT)	The percentage that a read count must reach for the software to call a base.
Assay type	The library prep kit used to process samples for sequencing and analysis.
Call	The reported base for a position (coordinate) based on the analysis method.
Cluster density	The number of thousands of clusters on a flow cell per square millimeter (K/mm ²).
Clusters	A clonal group of DNA bound to the surface of a flow cell that produces one read. One template DNA strand seeds each cluster and is clonally amplified using bridge amplification until the cluster has about 1000 copies.
Clusters passing filter	The percentage of clusters that pass quality filters. For more information, see Clusters Passing Filter (page 55) .
Control region	A subset of the human mtGenome (positions 16,008–594) containing three hypervariable regions.
Coverage	The number of reads for a locus or variant position in a sample that align to or cover known reference bases. Sequencing coverage level (referred to as <i>read depth</i> , <i>read count</i> , or <i>intensity</i> in the software) can help determine whether a variant call is made with confidence.
Coverage plot	A graph on a Results page that displays the coverage at each position in an mtDNA sample. For details, see Coverage Plot (page 15) .
Cycle	A chemistry step and an imaging step to call one base in a sequencing run. For cycle numbers, see Cycles per Read (page 57) .
Deletion	A base deleted from the sample sequence that is present in the reference sequence.
Differences indicator	A software feature that highlights different variant calls between two samples being compared.
Flow cell	A glass slide coated with oligos complementary to the library adapter sequences, allowing the library to adhere for a run.
Forward primer length	The number of bases in the forward PCR primer (oligonucleotide).
Forward primer start coordinate	The starting base at the 5' end of the forward PCR primer. The base depends on the amplicon and analysis method.
Hotspot	An mtDNA position with a variant, interpretation threshold, differences, no call, call indicator, or user modification.
Human Sequencing Control (HSC)	A pre-amplified and tagged positive control added to a library pool to facilitate troubleshooting of instrument and sequencing reagent issues.
Index adapter sequence	A short stretch of oligonucleotides added to each sample during library prep for identification purposes. The software uses the oligonucleotide sequences to demultiplex pooled libraries.

Universal Analysis Software – mtDNA Module Reference Guide

Term	Description
Index CV	The percentage of reads assigned to each library, representing the distribution of read counts of the libraries in the run. CV is the coefficient of variation for the number of read counts across all indexes.
Indel	The insertion or deletion of bases in a genome.
Insertion	A base in the sample sequence that is not present in the reference sequence.
Intensity	The software reports signal intensity as the number of reads.
Interpretation threshold (IT)	The percentage that determines the read count that must be reached for the software to call a base. Signal below the IT but greater than or equal to the AT can be manually typed.
Library	A DNA sample that uses primers and adapters containing specific regions of the genome for sequencing.
Library type	The function of a library, either primer-directed sequencing or non-directed sequencing, for custom assays and analysis.
Minimum quality score	The Q-score reads must meet to contribute to the variant call at each DNA coordinate.
Minimum read count	The minimum read count a position must reach for a base to be called.
mtDNA navigator	A section of a Results page that displays variants and other sequencing results for an mtDNA sample. For details, see mtDNA Navigator (page 14) .
N indicator	The quality indicator for reagent blanks and negative amplification controls.
Negative amplification control	A sample that does not contain DNA as a template for amplification.
No call	The call assigned when a position lacks read counts for a specific base.
Original call	An unmodified variant call based on the run analysis method.
P indicator	The quality indicator for positive amplification control DNA, such as HL60.
Percent A, C, G, T	The percentage of forward and reverse reads out of the total reads for the A, C, G, or T bases at the call position.
Percent del	The percentage of forward and reverse reads out of the total reads for a deletion.
Percent ref	The percentage of forward and reverse reads out of the total reads that do not support an insertion at an insertion point in mtDNA.
Percentages	A position viewer feature that shows the percentage of each nucleotide called at a position in relation to the thresholds.
Phasing	The percentage of bases in a cluster that fall behind the current cycle in Read 1. For details, see Phasing and Prephasing (page 55) .
Plexity	The number of libraries combined in a pool for sequencing. For example, if eight libraries are combined in a pool, the plexity is eight.
Position	A base coordinate of mtDNA using the numbering convention of the rCRS, comprised of positions 1–16,569.

Universal Analysis Software – mtDNA Module Reference Guide

Term	Description
Position viewer	A section of a Results page that displays reference and sample sequences and data for the selected position in an mtDNA sample. For details, see Position Viewer (page 14) .
Positions called	The number of positions out of the total number of interrogated mtDNA positions, including indels, for which the software called a base.
Positive amplification control	A known sequence of DNA provided in the library prep kit and prepared with forensic and reference samples to ensure that library prep and sequencing reagents are functional.
Prephasing	The percentage of bases in a cluster that jump ahead of the current cycle in Read 1. For details, see Phasing and Prephasing (page 55) .
Project	A collection of analyzed results for at least one sample, and the primary workspace for viewing and modifying results and creating reports.
Revised Cambridge Reference Sequence (rCRS)	The reference genome for calling variants in the mtGenome.
Read	The sequence string of A, T, C, and G bases corresponding to an amplicon. A run can simultaneously generate millions of reads from multiple loci and multiple samples.
Reads A, C, G, T	The number of forward and reverse reads for the A, C, G, or T base at the DNA call position.
Reads del	The number of forward and reverse reads for a deletion at the call position.
Reads ref	The number of forward and reverse reads that do not support an insertion at an insertion point.
Reagent blank	An extraction control without DNA that is carried through the entire sample processing.
Read 1 (i5) primer sequence	The sequence of the PCR primer that is tagged with the Index 2 (i5) adapter sequence. Index 2 is part of the adapter sequence, as is the binding site for the Read 1 sequencing primer.
Read 2 (i7) primer sequence	The sequence of the PCR primer that is tagged with the Index 1 (i7) adapter sequence. Index 1 is part of the adapter sequence, as is the binding site for the Read 2 sequencing primer.
Reagent cartridge kit lot	The batch code to identify the batch or lot that a component of the MiSeq FGx Reagent Kit or MiSeq FGx Reagent Micro Kit was manufactured in.
Reanalysis	The process of reanalyzing a sample with another analysis method, after the original run analysis method.
Reference call	The base call of the rCRS at a given position.
Reference SNP ID number (rsID)	An identification tag the National Center for Biotechnology Information (NCBI) assigns to a group of SNPs that map to the same location.
Reverse primer length	The number of bases in the reverse PCR primer (oligonucleotide).
Reverse primer start coordinate	The starting base at the 5' end of the reverse PCR primer. The base depends on the amplicon and analysis method.
Run	The process of clustering and sequencing libraries to generate base calls.

Universal Analysis Software – mtDNA Module Reference Guide

Term	Description
Run analysis method	The analysis method assigned to a sample for the initial run analysis.
Sample	The DNA input for a library prep kit, before the addition of primers and adapters.
Sample compare	A software feature that compares an mtDNA sample of interest to other mtDNA samples.
Sample of interest	A user-designated sample that is compared to other selected samples.
Sample plate	The 96-well plate that contains libraries for sequencing.
Sample representation	The number of reads per sample for a run. These data are available after sequencing and analysis are complete.
Sample type	The function of a sample: forensic or reference sample, positive amplification control, negative amplification control, and reagent blank.
Sample well position	The column and row position in a 96-well plate containing samples or libraries. For example, position A1 indicates a well at the intersection of row A and column 1.
Single-nucleotide polymorphism (SNP)	Variation of one base pair in a DNA sequence relative to a known DNA reference sequence.
Strand depth	The read count from the forward or reverse strand that has the majority of reads for a base position.
Strand direction	The selection of a forward or reverse strand for a primer-directed library type.
System event	An action the software automatically performs and records.
Total read count	The total number of reads from the forward and reverse DNA strands for a base position.
User action	A user interaction with the software, such as adding a sample to a run.
User-modified call	A base call that a user changed.
Variant	A base at a position in an mtDNA sequence that differs from the rCRS base and a synonym for SNP.
Whole mitochondrial genome	The entire 16,569 bp mtGenome.

Run Setup

Create a Run	21
Edit a Run	21
Add Samples to a Run	21
Manage Samples in a Run	23
Monitor a Run	25

Create a Run

1. On the Home page, select **Create Run**.
2. In the Create Run dialog box, complete the following fields. (Required fields are indicated on-screen with an asterisk *.)
 - **Run Name** — Enter a preferred name to identify the run.
 - **Description** — [Optional] Enter a description of the run.
 - **Flow Cell Type** — Select **Standard** or **Micro**, depending on your reagent kit.
 - **Reagent Cartridge Kit Lot #** — [Optional] Enter the LOT number printed on the reagent cartridge packaging.
3. Select **Save** to create the run.
The Run Details page opens.

Edit a Run

1. On the navigation rail, select  **Runs**. Then select a run.
2. Select  **Information**, and then select  **Edit**.
3. In the Edit Run dialog box:
 - Edit the run name.
 - Add or modify the optional run description.
 - Select a new flow cell type, **Standard** or **Micro**.
 - Add or modify the optional reagent cartridge kit lot number.
 Editing the flow cell type requires a run state of Created.
4. Select **Save**.

Add Samples to a Run

Add samples and sample information to a run either by importing a sample sheet or by adding individual samples. After samples are added, the Run Details page lists all samples included in the run.

All samples in a run must have the same assay type.

When adding samples to a run, you associate each sample with at least one project. If you enter the name of a new project, the software creates the project. Because a sample can be assigned to multiple projects, a run can include samples from multiple projects.

Import a Sample Sheet

1. On the navigation rail, select  **Runs**. Then select a run.
2. Select **Import Samples**.
3. [Optional] Select **Sample Sheet Template** to download a .txt file that provides examples of Sample Sheet formats for each Assay format.
4. Select **Drop Files Here**, and then navigate to the sample sheet location.
5. Select the sample sheet (.txt file), and then select **Open**.
The sample sheet appears under the file upload area. Uploading another .txt file overrides the current file because the software accepts only one sample sheet per run.
6. Select **Next**.
7. Review the information from the uploaded sample sheet to ensure accuracy. Scroll to see all entries. If utilizing the template, make sure to delete the sections pertaining to assay types that are not utilized in the current analysis.
8. Select **Import** to add the samples to the run.

Assigning Indexes

The sample sheet template includes two columns for entering index adapter names: Index 1 (i7 index) and Index 2 (i5 index).

Assigning Projects

The Initial Project List column of a sample sheet specifies which project a sample is assigned to. To assign multiple projects to a sample, enter each project name separated by a comma followed by a space and put quotes around the entry. For example:

```
"Project 1, Project 2, Project 3"
```

Add Individual Samples

1. On the navigation rail, select  **Runs**. Then select a run.
2. Select **Add Sample**.
3. In the Add Sample dialog box, complete the following fields (Required fields are indicated on-screen with an asterisk *.)

Field	Instruction
Sample Name	Enter a unique identifier for the sample. Each sample name must be unique in UAS, regardless of sample type.
Description	Enter a preferred description of the sample.

Universal Analysis Software – mtDNA Module Reference Guide

Field	Instruction
Assay Type*	Select the library prep method: <ul style="list-style-type: none">• ForenSeq mtDNA Whole Genome — Libraries prepared with the ForenSeq mtDNA Whole Genome kit.• ForenSeq mtDNA Control Region — Libraries prepared with the ForenSeq mtDNA Control Region kit.• mtDNA Custom — Libraries prepared with a custom or third-party assay.
Sample Type	Select a sample type: <ul style="list-style-type: none">• Sample — A forensic or reference sample.• For the positive amplification, select HL60.• Negative Amplification Control — The negative control (water).• Reagent Blank — The extraction control.
Index 1	Select the name of an i7 index adapter.
Index 2	Select the name of an i5 index adapter.
Run Analysis Method	Select a method for analyzing samples in the run.
Initial Projects	Enter the name of a project to assign the sample to. To assign the sample to multiple projects, enter multiple project names.
Sample Plate	Enter a preferred name for the 96-well plate.
Sample Well Position	Enter the position of the plate well the sample occupies.

* Because the software allows one assay type per run, the assay type assigned to the first sample is automatically assigned to the others.

4. Select **Add** to add the sample to the run.

Manage Samples in a Run

After adding samples to a run, you can edit sample information, change project assignments, or remove samples from the run. Use filters, sorting, and the search function to quickly find samples.

Find a Sample

1. On the navigation rail, select  **Runs**. Then select a run.
2. To search for samples, enter any of the following information in the Search field:
 - Sample name
 - Description
 - Index adapter name
 - Run analysis method
3. To filter the samples, do as follows.
 1. **a.** In the Sample Type or Assay Type column, select **Filter**.
 2. **b.** Select at least one checkbox.
 3. **c.** Reselect **Filter** to apply the selected filters.
4. To sort the samples, select a column heading.
The Name, Description, and Index Combination columns are sortable.

Edit Sample Information

1. On the navigation rail, select  **Runs**. Then select a run.
2. In the Actions column, select  **Edit**.
3. In the Edit Sample dialog box, modify the desired fields.

Field	Instruction
Sample Name	Edit the unique identifier for the sample.
Description	Enter or edit a preferred description of the sample.
Assay Type*	Select a different library prep method: <ul style="list-style-type: none"> • ForenSeq mtDNA Whole Genome — Libraries prepared with the ForenSeq mtDNA Whole Genome kit. • ForenSeq mtDNA Control Region — Libraries prepared with the ForenSeq mtDNA Control Region kit. • mtDNA Custom — Libraries prepared with a custom or third-party assay
Sample Type	Select a different sample type: <ul style="list-style-type: none"> • Sample — A forensic or reference sample. • For the positive amplification, select HL60. • Negative Amplification Control — The negative control (water). • Reagent Blank — The extraction control.
Index 1*	Select the name of a different i7 index adapter.
Index 2*	Select the name of a different i5 index adapter.
Run Analysis Method*	Select a different run analysis method.
Initial Projects*	Enter additional projects or remove current projects. For detailed instructions, see Modify Project Assignments (below) .
Sample Plate	Edit or enter a preferred name for the 96-well plate.
Sample Well Position	Edit or enter the position of the well the sample occupies.

* Editing these fields requires a run state of Created. To edit them after run completion, remove the sample from the run, make the necessary edits, and add the sample to the run again.

4. Select **Save** to apply the changes.

Modify Project Assignments

When the run state is Created, you can change which projects a sample is assigned to or assign the sample to additional projects.

1. On the navigation rail, select  **Runs**. Then select a run.
2. In the Actions column, select **Edit Sample Details** for the applicable sample.
3. In the Edit Sample dialog box, edit the Initial Projects field as follows.
 - To remove a sample from a project, select **X** next to the project name.
 - To add a sample to a project, enter the project name.
 - To add a sample to multiple projects, enter multiple project names.
 - To create a project for the sample, enter a new project name.

4. Select **Save**.

Remove Samples from a Run

1. On the navigation rail, select  **Runs**. Then select a run.
2. Select the checkbox of each sample you want to remove from the run.
3. Select **Remove Samples**.
4. In the Remove Samples dialog box, select **Confirm**.

Monitor a Run

During a run, the Run Details page indicates which cycle the run is on and displays metrics so you can monitor run status and performance. Green indicates passing metrics and orange indicates metrics that require further evaluation. These metrics are a subset of the run metrics used to evaluate performance after sequencing is complete.

1. On the navigation rail, select  **Runs**. Then select a run.
2. Select **Quality Metrics**.
 - Review the values for cluster density and clusters passing filter.
 - Review the phasing and prephasing values for Read 1.
 - Check Read 1, Index 1, Index 2, and Read 2 to confirm that all four reads passed.
 - Check Overall Intensity and 0 Discordant Loci to confirm that the HSC passed.

For more information, see [Quality Metrics \(page 54\)](#).

3. Select  **Run Activity**.
 - Review the run activity.
 - Select **Close** to return to the run.
4. Select **Sample Representation**.
 - Compare the samples to the sample read count guideline to confirm performance.
 - Review the number of reads in parentheses next to the sample name or hover over each sample.

For more information, see [Sample Representation \(page 60\)](#).

5. Select **Close** to return to the run.

Data Analysis

Review Run Metrics	26
Review Results	27
Add Analyzed Samples to a Project	29
Reanalyze Samples	29

Review Run Metrics

Run metrics provide run quality, results for positive and negative controls, and sample representation for an overall view of run performance. Green indicates passing metrics and orange indicates metrics that require further evaluation.

For help with metrics requiring further evaluation, see [Troubleshoot a Run \(page 49\)](#). For detailed descriptions of each metric, see [Appendix A Run Metrics](#).

1. On the navigation rail, select  **Projects**.
2. Select a project, and then select **QPN** for a run in the project.
3. On the Quality Metrics page, do as follows.
 - Review the values for cluster density and clusters passing filter.
 - Review the phasing and prephasing values for Read 1.
 - Check Read 1, Index 1, Index 2, and Read 2 to confirm that all four reads passed.
 - Check Overall Intensity and 0 Discordant Loci to confirm that the HSC passed.
4. Select **Positive & Negative Controls**.
 - Under HL60:
 - Expand **Discordance** to see the Discordance table.
 - Use the scroll arrows to move through all positions in the table.
 - Under Reagent Blanks, review the number of positions called for each reagent blank.
 - Under Negative Amplification Controls, review the number of positions called for each negative control.
 - Select a sample name to go to the Results page.
5. Select **Sample Representation**.
 - Compare the samples to the sample read count guideline to confirm performance.
 - Review the number of reads in parentheses next to the sample name or hover over each sample.
6. If any metrics do not meet the minimum requirements, see [Troubleshoot a Run \(page 49\)](#) for troubleshooting guidance. For additional assistance, contact Verogen Technical Support. Collect the following items from D:\Illumina\MiSeqFGxOutput to help Verogen with troubleshooting:
 - InterOp folder
 - runInfo.xml file
 - runParameters.xml file

Review Results

The Results page for each mtDNA sample includes filters, fields, and buttons for a dynamic view of results for the sample. The following sections describe how to use these features to review the results. For feature descriptions, see [Results Page \(page 13\)](#).

Filter Results

1. On the navigation rail, select  **Projects**.
2. Select a project, and then select a sample to open the Results page.
3. Select **Filter**, and then select at least one of the following checkboxes:
 - **Control Region** — Zooms in on the control region of the mtGenome.
 - **Interpretation Threshold** — Indicates positions with reads between the AT and IT.
 - **Variants** — Indicates positions with base calls that differ from the reference.
 - **Call** — For reagent blanks and negative amplification controls, indicates positions with calls to identify possible contamination.
 - **No call** — For samples and HL60, indicates positions with no calls to identify areas without coverage.

The sample type determines whether Call or No Call filter is visible.

4. Reselect **Filter** to apply the selected filters.
The mtDNA navigator, position viewer, and coverage plot refresh to reflect the filtered data.
5. To clear the filters:
 1. **a.** Select **Active Filter**.
 2. **b.** Clear the checkboxes.
 3. **c.** Reselect **Active Filter**.

Set Maximum Coverage

Setting a maximum coverage filters sample data by total read count.

1. On the navigation rail, select  **Projects**.
2. Select a project, and then select a sample to open the Results page.
3. In the Max Coverage field, enter a read number without commas.
The mtDNA navigator and coverage plot adjust to the specified value.

Compare Samples

You can compare and filter results among 2–9 mtDNA samples in a project. During a comparison, the mtDNA navigator and coverage plot show data for the sample of interest in gray and a selected comparison sample in blue. The position viewer displays sequences and data for the sample of interest and selected comparison sample. The name of the sample of interest appears blue.

1. Navigate to the Results page:

Universal Analysis Software – mtDNA Module Reference Guide

- On the navigation rail, select  **Projects**.
 - Select a project, and then select a sample to be the sample of interest.
 - On the Samples page, select a sample.
This sample is the sample of interest.
 - On the Results page, select **Compare**.
 - In the Sample Compare dialog box, select the checkboxes of up to eight samples.
 - Select **Submit**.
The position viewer lists the selected comparison samples. The mtDNA navigator and coverage plot refresh to show data for the sample of interest and the first selected comparison sample. A green checkmark next to a selected comparison sample indicates no differences with the sample of interest.
 - If you are comparing more than two samples, select a comparison sample in the position viewer.
The mtDNA navigator and coverage plot refresh to compare the selected comparison sample to the sample of interest.
 - [Optional] Select **Filter**, and then select at least one of the following checkboxes:
 - Control Region** — Zooms in on the control region of the mtGenome.
 - Interpretation Threshold** — Indicates positions with reads between the AT and IT.
 - Variants** — Indicates positions with base calls that differ from the reference.
 - Call** — For reagent blanks and negative amplification controls, indicates positions with calls to identify possible contamination.
 - No call** — For samples and HL60, indicates positions with no calls to identify areas without coverage.
 - Differences** — Indicates call differences between the sample of interest and the selected comparison samples.
- The sample of interest determines whether the Call or No Call filter is displayed. For example, when the sample of interest is a negative amplification control or reagent blank, the software displays the Call filter. When the sample of interest is a sample or positive amplification control, the software displays the No Call filter.
- When finished comparing samples, select **Stop Compare**.

Edit Base Calls

- On the navigation rail, select  **Projects**.
- Select a project, and then select a sample to open the Results page.
- In the position viewer, select a position by using the scroll arrows or entering the position.
- Select the **Call** arrow, and then select a call. Scroll to see all available calls.
The position viewer refreshes to reflect the modified call. An asterisk (*) marks the modified call and the original call appears below the Call list.

Add Analyzed Samples to a Project

You can add analyzed samples in a project to other projects. After adding a sample to another project, sample modifications are **not** applied in the newly assigned project.

1. On the navigation rail, select  **Projects**. Then select a project.
2. On the Samples page, select **Add Samples**.
3. In the Add Samples dialog box, in the **Search** field, enter a sample name.
4. In the search results, select the checkbox of each sample you want to add to the project.
Each selection appears in a summary list with the corresponding analysis method.
5. To remove a sample, clear the checkbox or select **X** in the summary.
6. Select **Add** to add the samples.

Reanalyze Samples

Reanalyze samples to apply an analysis method with user-defined parameters. The software preserves results from all analysis methods.

1. On the navigation rail, select  **Projects**.
2. Select a project.
3. On the Samples page, select the checkbox of each sample you want to reanalyze.
4. Select **Reanalyze**.
5. In the Reanalyze Samples dialog box, choose one of the following options:

Reanalysis Option	Instruction
Assign the same analysis method to all samples selected for reanalysis.	Select a method from the Analysis Methods list, and then select Apply to All .
Assign an analysis method to each sample selected for reanalysis.	Select a method from each list in the New Analysis Method column.

The software lists only analysis methods that are compatible with the assay used to prepare the sample.

6. Select **Next**.
7. Review the sample names, previous analysis methods, and new analysis methods.
8. To make changes, select **Previous**.
9. Select **Save** to initiate reanalysis.
In the Projects sidebar, an icon appears next to each sample undergoing reanalysis.
10. When reanalysis is complete, select a sample to view results.
The Results page displays results for the most recent analysis.

Switch Analysis Methods

A Results page shows results for the analysis method displayed in the Select an Analysis Method list. For reanalyzed samples, the list includes all analysis methods applied to the sample with the most recent first. Select different analysis methods to see other results.

Universal Analysis Software - mtDNA Module Reference Guide

1. On the navigation rail, select  Projects.
2. Select a project, and then select a sample to open the Results page.
3. In the Select an Analysis Method list, select an analysis method.
The Results page is refreshed to show results for the selected method.

Reports

Generating Reports	31
File Formats for Reports	35

Generating Reports

The software compiles mtDNA results into the following reports:

- **CODIS** — Variants from the selected samples for upload to CODIS.
- **EMPOP** — Specimen data, originating agency identifiers (ORIs), and other DNA profile information for upload to EMPOP.
- **Project** — Results for all samples in a project.
- **Sample** — Results for one sample.

The following sections describe how to generate each type of report. For report contents, see [File Formats for Reports \(page 35\)](#).

Create a CODIS Report

Create a CODIS report for one sample or multiple samples. Creating a report for one sample lets you specify extra specimen information and choose which positions to include in the report. When creating a report for multiple samples, the software automatically includes all positions for each sample.

Create a CODIS Report for One Sample

1. On the navigation rail, select  **Projects**.
2. Select a project, and then select a sample.
3. Select  **Reports**, and then select **mtDNA CODIS Report**.
The Create CODIS Report dialog box opens.
4. [Optional] Select **CODIS Report Defaults** to edit default report settings.
For more information, see [CODIS Report Defaults \(page 46\)](#).
5. Edit the following CMF header fields or accept the default values.
 - **Destination Laboratory ORI** — Identifier for the agency that processes the CODIS report.
 - **Source Laboratory ORI** — Identifier for the agency that processes the specimen.
 - **Submission User ID** — Identification of the person who generates the CODIS report.
 - **Batch ID** — Identifier that tracks where a DNA profile originated.
6. To include the name of the library prep kit in the report, select **Include Kit Name in Report**.
7. Select **Next**.
8. In the Specimen ID field, enter a unique identifier for the specimen or accept the default sample name.
9. In the Specimen Category list, select a preconfigured category or select **Custom** and enter the category name.
10. [Optional] Select **More Options**.
 1. **a.** Complete the following fields:

Universal Analysis Software – mtDNA Module Reference Guide

- **Source ID** — Specify whether the identity of the specimen contributor is known.
 - **NCIC Number** — Enter the unique number for storage in the NCIC system.
 - **ViCAP number** — Enter the unique number for storage in the ViCAP system.
 - **Case ID** — Enter the unique identifier for the case.
 - **Case ID** — Enter the law enforcement identifier associated with the specimen.
 - **Specimen Comment** — Enter any additional information about the specimen.
2. **b.** If the specimen has possible no call regions or other conditions requiring additional information to help with match resolution, select **Partial Profile**.
11. Select **Next**.
 12. Select the toggle key for each position you want to remove from the report.
By default, all positions are included.
 13. Select **Create** to generate the report.
 14. Select  **Reports** to open the Reports page.
The Reports page lists all reports for the project with the new report on top.
 15. In the Actions column, select **Download Report** to download the report as a .cmf file.

Create a CODIS Report for Multiple Samples

1. On the navigation rail, select  **Projects**.
2. Select a project, and then select  **Reports** to open the Reports page.
3. Select **Create Report**, and then select **mtDNA CODIS Report**.
The Create CODIS Report dialog box opens.
4. Select the checkbox of each sample you want to include in the report.
 - Scroll and use the arrows to move through the complete list of samples in the project.
 - Under Summary, select **Delete** to remove samples from the report.
5. Select **Next**.
6. [Optional] Select **CODIS Report Defaults** to edit default report settings.
For more information, see [CODIS Report Defaults \(page 46\)](#).
7. Edit the following CMF header fields or accept the default values:
 - **Destination Laboratory ORI** — Identifier for the agency that processes the CODIS report.
 - **Source Laboratory ORI** — Identifier for the agency that processes the specimen.
 - **Submission User ID** — Identification of the person who generates the CODIS report.
 - **Batch ID** — Identifier that tracks where a DNA profile originated.
8. To include the name of the library prep kit in the report, select **Include Kit Name in Report**.
9. Select **Next**.
10. In the Specimen ID field for each sample, enter a unique identifier for the specimen cited in the report.
11. Assign specimen categories per the following table. If you select Custom, enter a name for the custom specimen category.

Option	Instruction
Assign the same category to all samples.	In the Specimen Category list at the top of the dialog box, select a category, and then select Apply to All Specimens .
Assign a different category to each sample.	In the Specimen Category list for each sample, select the applicable category.

12. Select **Create** to generate the report.

The software returns to the Reports page, where the new report appears at the top of the list.

13. In the Actions column, select **Download Report** to download the report as a .cmf file.

Create an EMPOP Report

Create an EMPOP report for one or more samples in a project. When creating a report for one sample, you can start from the Results or Reports page.

Create an EMPOP Report for One Sample from Results

1. On the navigation rail, select  **Projects**.
2. Select a project, and then select a sample to open the Results page.
3. Select  **Reports**, and then select **mtDNA EMPOP Report** to generate the report.
4. Select  **Reports** to open the Reports page.
The Reports page lists the new report on top.
5. In the Actions column, select **Download** to download the report as a .txt file.

Create an EMPOP Report for One Sample from Reports

1. On the navigation rail, select  **Projects**.
2. Select a project, and then select  **Reports**.
The Reports page opens.
3. Select **Create Report**, and then select **mtDNA EMPOP Report**.
4. In the Create mtDNA EMPOP Report dialog box, select the checkbox of the applicable sample.
5. Select **Confirm** to generate the report.
The returns to the Reports page, where the new report appears at the top of the list.
6. In the Actions column, select **Download Report** to download the report as a .txt file.

Create an EMPOP Report for Multiple Samples

1. On the navigation rail, select  **Projects**.
2. Select a project, and then select  **Reports**.
The Reports page opens.
3. Select **Create Report**, and then select **mtDNA EMPOP Report**.
4. In the Create EMPOP Report dialog box, select the checkbox of each sample to include in the report.
 - Scroll and use the arrows to move through the complete list of samples in the project.
 - Under Summary, select **Delete** to remove samples from the report.

5. Select **Confirm** to generate the report.
The returns to the Reports page, where the new report appears at the top of the list.
6. In the Actions column, select **Download Report** to download the report as a .txt file.

Create a Project Report

1. On the navigation rail, select  **Projects**.
2. Select a project, and then select  **Reports** to open the Reports page.
3. Select **Create Report**, and then select **Project**.
If the software is analyzing samples, project report creation is disabled and a notification is displayed.
4. When prompted, select **Yes** to create the report.
The software generates the report and returns to the Reports page, where the new report appears at the top of the list.

If a project includes samples with different assay types, the software generates one project report per assay.

5. In the Actions column, select **Download** to download a report as an .xmlx file.

Create Sample Reports

Create one sample report at a time or simultaneously generate multiple sample reports, one for each selected sample in a project. When creating a report for one sample, you can start from the Results or Reports page.

Create One Sample Report from Results

1. On the navigation rail, select  **Projects**.
2. Select a project, and then select a sample to open the Results page.
3. Select  **Reports**, and then select **Sample Report** to generate the report.
4. Select  **Reports** to open the Reports page.
The Reports page lists the new report on top.
5. In the Actions column, select **Download** to download the report as an .xmlx file.

Create One Sample Report from Reports

1. On the navigation rail, select  **Projects**.
2. Select a project, and then select  **Reports** to open the Reports page.
3. Select **Create Report**, and then select **Sample**.
4. In the Create Sample Report dialog box, select the checkbox of the sample you want to create a report for. Scroll and use the arrows to see all samples.
5. Select **Confirm** to generate the report.
The software returns to the Reports page, where the new report appears at the top of the list.
6. In the Actions column, select **Download** to download a report as an .xmlx file.

Create Multiple Sample Reports

1. On the navigation rail, select  **Projects**.
2. Select a project, and then select  **Reports** to open the Reports page.
3. Select **Create Report**, and then select **Sample**.
4. In the Create Sample Report dialog box, select the checkbox of each sample you want to create a report for.
 - Scroll and use the arrows to move through the complete list of samples in the project.
 - Use the Delete icons under Summary to remove samples from the report.
5. Select **Confirm** to generate the reports.

The software returns to the Reports page, where the new reports appear at the top of the list.
6. In the Actions column, select **Download** to download a report as an .xml file.

File Formats for Reports

The following sections describe the format and contents of each type of report the software generates for mtDNA samples. For instructions on creating reports, see [Generating Reports \(page 31\)](#).

CODIS Report

A CODIS report is a .cmf file that includes a CMF header, reference sequence, and specimen profile. The specimen profile provides the following information for each sample in the report:

- Specimen ID
- Specimen category
- Source ID
- Fragments

EMPOP Report

An EMPOP report is a .txt file. For each sample, the report includes a sequence string in FASTA format, variant calls in rCRS format, and names for the following items:

- Project
- Sample
- Assay
- Analysis method
- Run

Universal Analysis Software – mtDNA Module Reference Guide

- Analysis method
- Run name
- Analytical and interpretation thresholds
- Minimum Q-score
- Minimum read count
- SNVs
- Insertion and deletion counts
- Call count for reagent blanks and negative controls
- No call count for samples and positive controls
- Manually changed calls
- Variant details:
 - Variant and reference calls
 - Total read count
 - Strand depth
 - Reads A, C, T, and G
 - Percent A, C, T, and G
 - Reads del
 - Percent del
 - Reads ref
 - Percent ref

Run, Project, and Sample Management

Create a Project	38
Edit a Project	38
Obtain the Analysis Folder Path	39
Printing	39
History	39

Create a Project

Create a project to group a selection of samples. For example: create a project of controls to compare controls from multiple projects, or group samples in a comparison to generate a project report of compared samples.

Although possible, combining samples with different assay types in a project is not recommended.

1. On the navigation rail, select  **Projects**.
2. Select **Create Project**.
3. In the Create Project dialog box, complete the following fields:
 - **Project Name** — Enter a preferred name to identify the project.
 - **Description** — Optionally enter a description of the project.
4. Select **Save** to create the project.
The new project opens to the Samples page.
5. Select **Add Samples**.
You can add any samples in the software to the project, analyzed or not.
6. In the Add Samples dialog box, in the **Search** field, enter a sample name.
7. In the search results, select the checkbox of each sample you want to add to the project.
Each selection appears in a summary list with the corresponding analysis method.
8. To remove a sample, clear the checkbox or select **X** in the summary.
9. Select **Add** to add the samples.
The Samples page lists the added samples. Runs the samples belong to appear in the Project sidebar.

Edit a Project

1. On the navigation rail, select  **Projects**.
2. Select a project, and then select  **Edit**.
3. In the Edit Project dialog box, do as follows.
 - Edit the project name.
 - Add or modify the optional project description.
4. Select **Save** to apply the changes.

Obtain the Analysis Folder Path

The analysis folder resides on the ForenSeq Universal Analysis Server and contains the Variant Call Format (VCF) and Binary Alignment Map (BAM) files for the analysis.

1. On the navigation rail, select  **Projects**.
2. Select a project, and then select a sample to open the Results page.
3. Select the Information icon.
4. The path is displayed under **Analysis Folder Path**.

Printing

The  **Print** icon converts the following pages into a printable format:

- **Run Details** — Prints the header of a Run Details page, the table listing samples in the run, run quality metrics, and sample representation.
- **Samples** — Prints the Samples page of a project.

Print a Run Details Page

1. On the navigation rail, select  **Runs**. Then select a run.
2. On the Run Details page, select  **Print**.
The software displays a printable view of the Run Details page.
3. Select **Print** to open the Print dialog box.
4. Choose printing options and print.
5. Select **Back to Run**.

Print a Samples Page

1. On the navigation rail, select  **Projects**. Then select a project.
2. On the Project sidebar, select  **Print**.
The software displays a printable view of the Samples page.
3. Select **Print** to open the Print dialog box.
4. Choose printing options and print.
5. Select **Back to Project**.

History

When the history setting is enabled, you can review system events and user actions recorded for a run, project, or sample and add comments. For example, you can comment on a sample to explain a call modification. For instructions on enabling the history setting, see [System Settings \(page 47\)](#).

Review Run Activity

1. On the navigation rail, select  **Runs**.

Universal Analysis Software – mtDNA Module Reference Guide

2. Select a run, and then select  **Run Activity**.
3. In the Run Activity dialog box, review system events and user actions for the run.
The dialog box lists system events and user actions from newest to oldest.
4. Navigate using the scroll bar and arrows. Narrow results using the following fields:
 - **Search** — Search the current history.
 - **Start Date and End Date** — View the history within a specified time frame.
 - **Show System Events** — See system events and user actions for the history or see user actions only.
5. When finished, select **Close**.

Review Project Activity

1. On the navigation rail, select  **Projects**.
2. Select a project, and then select  **Project Activity**.
3. In the Project Activity dialog box, review system events and user actions for the project.
The dialog box lists system events and user actions from newest to oldest.
4. Navigate using the scroll bar and arrows. Narrow results using the following fields:
 - **Search** — Search the current history.
 - **Start Date and End Date** — View the history within a specified time frame.
 - **Show System Events** — See system events and user actions for the history or see user actions only.
5. When finished, select **Close**.

Review Sample Activity

1. On the navigation rail, select  **Projects**.
2. Select a project, select a sample, and then select  **Sample Activity**.
3. In the Sample Activity dialog box, review system events and user actions for the sample.
The dialog box lists system events and user actions from newest to oldest.
4. Navigate using the scroll bar and arrows. Narrow results using the following fields:
 - **Search** — Search the current history.
 - **Start Date and End Date** — View the history within a specified time frame.
 - **Show System Events** — See system events and user actions for the history or see user actions only.
5. When finished, select **Close**.

Comment on a Run

1. On the navigation rail, select  **Runs**.
2. Select a run, and then select  **Run Activity**.
The Run Activity dialog box opens.
3. To add a comment to the history:
 1. **a.** In the Add New Comment field, type a comment.
 2. **b.** Select **Add** to save the comment.

Universal Analysis Software – mtDNA Module Reference Guide

The comment appears as an event at the top of the history.

4. To comment on an entry in the history:
 1. **a.** Select the system event or user action. Use the fields and scrolling to help locate the entry, if necessary.
 2. **b.** In the Edit Comment box, type a comment.
 3. **c.** Select **Save** to add the comment to the system event or user action.

The comment appears as the last item for the entry.

5. When finished, select **Close**.

Comment on a Project

1. On the navigation rail, select  **Projects**.
2. Select a project, and then select  **Project Activity**.

The Project Activity dialog box opens.

3. To add a comment to the history:
 1. **a.** In the Add New Comment field, type a comment.
 2. **b.** Select **Add** to save the comment.

The comment appears as an event at the top of the history.

4. To comment on an entry in the history:
 1. **a.** Select the system event or user action. Use the fields and scrolling to help locate the entry, if necessary.
 2. **b.** In the Edit Comment box, type a comment.
 3. **c.** Select **Save** to add the comment to the system event or user action.

The comment appears as the last item for the entry.

5. When finished, select **Close**.

Comment on a Sample

1. On the navigation rail, select  **Projects**.
2. Select a project, select a sample, and then select  **Sample Activity**.

The Sample Activity dialog box opens.

3. To add a comment to the history:
 1. **a.** In the Add New Comment field, type a comment.
 2. **b.** Select **Add** to save the comment.

The comment appears as an event at the top of the history.

4. To comment on an entry in the history:
 1. **a.** Select the system event or user action. Use the fields and scrolling to help locate the entry, if necessary.
 2. **b.** In the Edit Comment box, type a comment.
 3. **c.** Select **Save** to add the comment to the system event or user action.

The comment appears as the last item for the entry.

Universal Analysis Software – mtDNA Module Reference Guide

5. When finished, select **Close**.

Software Settings

Settings Overview	43
Analysis Methods	43
CODIS Report Defaults	46
User Management	47
System Settings	47
Enable or Disable Dark Mode	48
Display Software Version	48

Settings Overview

Use the following settings, which are available from the **Settings** icon on the navigation rail, to configure the software. All settings except Dark Mode are visible to administrators only.

- **Analysis Methods** — Create and manage analysis methods.
- **CODIS Report Defaults** — Set default values to prepopulate CODIS reports.
- **User Management** — Add users and assign permissions.
- **System Settings** — Set visibility into system events and user actions.
- **System Information** — View software version numbers for various bio-informatics analysis methods.
- **Dark Mode** — Change the color scheme of the user interface so it is color-blind compatible.

Analysis Methods

An analysis method is a collection of settings that inform how the software analyzes sequencing data and performs variant calling. The software ships with one or more default Verogen analysis methods.

- ForenSeq mtDNA Control Region — Verogen mtDNA Control Region Analysis Method
- ForenSeq mtDNA Whole Genome — Verogen mtDNA Whole Genome Analysis Method
- mtDNA Custom — Verogen mtDNA Custom Analysis Method

Create an Analysis Method

A default Verogen analysis method cannot be modified, but you can create an analysis method with the desired settings in one of two ways: copy an existing analysis method and modify the settings or start with a blank analysis method.

Performance is supported when using the default analysis methods. Perform internal validation studies as needed to develop settings for new analysis methods.

For information on how Verogen determined the default settings in the mtDNA analysis methods and help determining alternative settings, see *Analysis Settings in ForenSeq Universal Analysis Software v2.0 - Verogen analysis methods offer fast, flexible and accurate variant calling or mitochondrial DNA data (document # VD202004)*.

Copy an Analysis Method

1. On the navigation rail, select  **Settings**.
2. Select **Analysis Methods**.
3. Select **Create Analysis Method**.
The New Analysis Method dialog box opens.
4. In the New Analysis Method Options list, select **Copy an Existing Analysis Method**.
5. In the Source Analysis Method list, select an analysis method to use as a template.
6. In the New Analysis Method Name field, enter a unique name for the new analysis method.
7. Select **Save** to create the analysis method.
The new method appears with the default settings from the copied method.
8. Modify the default settings as desired.
For instructions, see [Modify Analysis Settings \(below\)](#).

Use a Blank Analysis Method

1. On the navigation rail, select  **Settings**.
2. Select **Analysis Methods**.
3. Select **Create Analysis Method**.
The New Analysis Method dialog box opens.
4. In the New Analysis Method Options list, select **New Blank Analysis Method**.
5. In the New Analysis Method Name field, enter a unique name for the new analysis method.
6. Select **Save** to create the analysis method.
The new method appears with blank settings.
7. Enter the desired settings for the analysis method.
For instructions, see [Modify Analysis Settings \(below\)](#).

Modify Analysis Settings

For a user-created mtDNA analysis method, you can modify the AT, IT, minimum Q-score, and minimum read count. If an analysis method is copied from the Verogen mtDNA Custom Analysis Method, you can also modify the library type and amplicons.

Edit an mtDNA Whole Genome or mtDNA Control Region Method

1. If necessary, navigate to the analysis method you want to edit:
 1. **a.** On the navigation rail, select  **Settings**.
 2. **b.** Select **Analysis Methods**.
 3. **c.** In the Select an Assay Type list, select the assay that corresponds to the applicable analysis method.
 4. **d.** Select the applicable analysis method.

Universal Analysis Software – mtDNA Module Reference Guide

2. To rename the analysis method, do as follows.
 1. **a.** Select  **Edit** to open the Edit Analysis Info dialog box.
 2. **b.** In the Analysis Method Name field, enter a new name.
 3. **c.** Select **Save** to apply the name change and close the dialog box.
3. Edit the following settings as desired:
 - In the Analytical Threshold and Interpretation Threshold boxes, adjust percentages.
 - In the Minimum Quality Score and Minimum Read Count boxes, adjust values.All percentages and values are greater than or equal to the value entered.
4. Select **Save** to apply the new settings.

Edit an mtDNA Custom Analysis Method

1. If necessary, navigate to the analysis method you want to edit:
 1. **a.** On the navigation rail, select  **Settings**.
 2. **b.** Select **Analysis Methods**.
 3. **c.** In the Select an Assay Type list, select the assay that corresponds to the applicable analysis method.
 4. **d.** Select the applicable analysis method.
2. To rename the analysis method, do as follows.
 1. **a.** Select  **Edit** to open the Edit Analysis Info dialog box.
 2. **b.** In the Analysis Method Name field, enter a new name.
 3. **c.** Select **Save** to apply the name change and close the dialog box.
3. Edit the following settings as desired:
 - In the Analytical Threshold and Interpretation Threshold boxes, adjust percentages.
 - In the Minimum Quality Score and Minimum Read Count boxes, adjust values.All percentages and values are greater than or equal to the value entered.
4. In the Library Type list, select one of the following options:
 - **Non-Directed Sequencing**—Either strand can be sequenced as Read 1 or Read 2. This option is the default.
 - **Primer-Directed Sequencing**—One strand is sequenced as Read 1 and the opposite strand as Read 2.The software lists amplicons available for the selected library type.
5. Modify the amplicons. (Required fields are indicated on-screen with an asterisk *.)
 - To remove an amplicon, in the Actions column, select **Delete**.
 - To edit an amplicon, in the Actions column, select **Edit**. In the Edit Amplicon dialog box, edit the desired fields, and then select **Save**.
 - To add an amplicon, select **Add Amplicon**. In the Add Amplicon dialog box, complete the necessary fields, and then select **Confirm**.
6. Select **Save** to apply the new settings.

Delete an Analysis Method

You can delete any user-created analysis method. The default Verogen analysis methods cannot be deleted.

Universal Analysis Software – mtDNA Module Reference Guide

1. On the navigation rail, select  **Settings**.
2. Select **Analysis Methods**.
3. In the Select an Assay Type list, select the assay that corresponds to the analysis method you want to delete.
4. Select the applicable analysis method.
5. In the upper-right corner of the analysis method, select **Delete Analysis Method**.
6. Select **OK** to confirm the deletion.

Change the Default Analysis Method

1. On the navigation rail, select  **Settings**.
2. Select **Analysis Methods**.
3. In the Select an Assay Type list, select the kit that corresponds to the analysis method you want to set as the default.
4. Select the applicable analysis method.
A green star identifies the current default.
5. Set the selected analysis method as the default for the assay type by selecting the empty star next to the analysis method name:

CODIS Report Defaults

To expedite report creation, define default values to appear in the CMF header of every CODIS report. You can override these values on a per-report basis.

Define CODIS Report Defaults

1. On the navigation rail, select  **Settings**.
2. Select **CODIS Report Defaults**.
3. Select **mtDNA CODIS Report** for mtDNA samples.
4. In each of the following fields, enter a default value:
 - **Destination Laboratory ORI** — Identifier for the agency that processes the CODIS report.
 - **Source Laboratory ORI** — Identifier for the agency that processes the specimen.
 - **Submission User ID** — Identification of the person who generates the CODIS report.
5. Select **Save** to apply the default values.

Set Specimen Categories

The Specimen Categories setting determines which specimen categories are available to assign to samples when creating a CODIS report. By default, all categories are available.

1. On the navigation rail, select  **Settings**.
2. Select **CODIS Report Defaults**.
3. Select **mtDNA CODIS Report** for mtDNA samples.
4. Select the Specimen Categories tab.
This tab lists all specimen categories, divided into CODIS and Custom columns.
5. Select a specimen category to make it available or unavailable for reports.

A blue toggle key indicates an available category.

Add Custom Specimen Categories

Add custom specimen categories to the software to supplement the predefined categories. Custom categories can be renamed, enabled or disabled, and deleted.

1. On the navigation rail, select  **Settings**.
2. Select **CODIS Report Defaults**.
3. Select **mtDNA CODIS Report** for mtDNA samples.
4. Select the **Specimen Categories** tab.
5. In the New Specimen Category Name field, enter a name of your choice.
6. Select **Add**.

The new category appears in the Custom column and is enabled by default.

User Management

The Admin Settings page lists each account by username. On this page, an administrator can manage permissions for each account, including granting access to new accounts and assigning access levels.

Enable and Disable Accounts

New accounts are disabled by default and must be enabled by an administrator.

1. On the navigation rail, select  **Settings**.
2. Select **User Management**.
3. Select **Enabled** to enable or disable an account.

Blue indicates an enabled account.

Assign Access Levels

1. On the navigation rail, select  **Settings**.
2. Select **User Management**.
3. Select **Administrator** to enable or disable administrator access for an account.

Blue indicates an administrator account.

System Settings

The System Settings page hosts the history setting. When enabled, this setting lets you view system events and user actions for runs, projects, and samples and includes these activities in project and sample reports.

By default, the history setting is enabled. However, the software **always records activity**, regardless of whether the setting is enabled. When enabling the setting, you can choose whether to store or delete activities that occurred since the last disabling.

Enable Run, Project, and Sample History

1. On the navigation rail, select  **Settings**.

2. Select **System Settings**.
3. Select **Enable or disable visibility of system events and user actions recorded by the software**.
4. In the dialog box, select whether to store or delete previous system events and user actions:
 - **Permanently Delete All Previous History** — Start recording and delete all previous activity. You can access only activities recorded from this point forward.
 - **Store Previous History** — Continue recording and grant access to all activity starting from the creation of a run, project, or sample.

An  Activity icon appears on the Run Details page, Project sidebar, and Results page. Selecting the icon opens the history.

Disable Run, Project, and Sample History

1. On the navigation rail, select  **Settings**.
2. Select **System Settings**.
3. Select **Enable or disable visibility of system events and user actions recorded by the software**.
The History icon disappears from the Run Details page, Project sidebar, and Results page.

Enable or Disable Dark Mode

1. On the navigation rail, select  **Settings**.
2. Select **Dark Mode** to darken or brighten the user interface.

Display Software Version

1. On the navigation rail, select  **Settings**.
2. Select **System Information**.

Troubleshooting

Troubleshoot a Run	49
Perform a Status Check	53

Troubleshoot a Run

The following sections provide recommendations for troubleshooting run problems. Most recommendations include library prep steps, which are documented in the library prep kit reference guides. Visit verogen.com/documentation to download the guide for your kit.

Low Intensity

Low intensity (number of reads) can indicate a problem with sample quantity or quality.

- Make sure that DNA concentrations are not below the targeted amount.
- Depending on the source, purify or dilute the DNA, prepare libraries, and resequence the sample.
- Resequence with fewer samples, maintaining a minimum of eight samples per run.

High Cluster Density

Review run metrics: If cluster density is high, but sample representation shows that reads are above the sample read count guideline, proceed with analysis. The following factors can result in insufficient data for analysis.

- During normalization, library from the Purified Library Plate or the Normalization Storage (LNB1) were overpipetted.
 - If the sample has a high number of alleles with quality control indicators, be sure to use correct volumes and reagents for normalization and resequence the sample.
 - If the sample volume is correct, further dilute libraries in the Pooled Normalized Libraries tube.
 - If normalization volumes or reagents were incorrect, repeat normalization and resequence.
- Too much adapter dimer in library prep.
 - If your laboratory has a Bioanalyzer or Fragment Analyzer, analyze 1–2 µl library from the Purified Library Plate to check the amount of primer dimer.
 - If the checked library has > 5% adapter dimer, repeat library purification using the 30 µl library remaining in each well of the Purified Library Plate. Adjust purification to use 30 µl Sample Purification Beads 2 (SPB2), add 32.5 µl Resuspension Buffer (RSB), and transfer 30 µl to a new 96-well PCR plate. Proceed with normalization.
 - For questions, contact Verogen Technical Support.
- An insufficient volume of Hybridization Buffer (HT1) was added, or too much normalized library was added to HT1.
 - If the sample has a high number of alleles with quality control indicators, resequence the sample. Make sure the normalized library is sufficiently diluted.

Low Cluster Density

Review run metrics: If cluster density for the run is low, but sample representation shows that reads are above the sample read count guideline, proceed with analysis. If the sample has an insufficient number of reads due to low density, the problem might be due to the following factors.

- DNA input was too low or too degraded, or the DNA dilution was incorrect. If possible, repeat library prep with more DNA.
- A critical reagent was not added during library prep. Repeat the prep.
- LNB1 was insufficiently pipetted during normalization. Consider diluting libraries in the Pooled Normalized Libraries tube less.
- A larger volume than expected of HT1 was pipetted during library prep. Repeat sequencing using the correct volume.
- Excessive sample was lost during purification. Repeat the library prep using best practices for bead handling.
- After heating, the bead-based normalized library pool was not immediately added to HT1. If possible, repeat the sequencing run following the protocol.
- The library was not otherwise denatured properly.
 - Make sure that the final library pool is denatured with freshly diluted 2 N NaOH (HP3) that has a pH > 12.5, if applicable.
 - Make sure that the final library pool is heat-denatured, if applicable.
- MiSeq FGx Maintenance Status:
 - Perform a maintenance wash and repeat sequencing.
 - Refill the wash tray and bottle for every wash.
 - Make all wash solutions daily.

Low Percentage of Clusters Passing Filter

If the percentage of clusters passing filter is too low, the problem might be due to the following factors.

- The cluster density is too high; see [High Cluster Density \(previous page\)](#). However, if sample representation shows that reads are above the sample count guideline, proceed with analysis.
- MiSeq FGx Maintenance Status:
 - Perform a maintenance wash and repeat sequencing.
 - Refill the wash tray and bottle for every wash.
 - Make all wash solutions daily.
- If the problem remains unresolved, reagents might not be performing as expected. Contact Verogen Technical Support.

High Phasing and Prephasing

The following factors can cause high phasing or prephasing.

- The cluster density is too high. See [High Cluster Density \(previous page\)](#).
- The ambient temperature is too high.

Universal Analysis Software – mtDNA Module Reference Guide

- Verify the appropriate temperature in the *MiSeq FGx Sequencing System Site Prep Guide (document # VD2018012)*.
- Check the MiSeq FGx System temperature logs to confirm whether the ambient temperature is too high. Contact Verogen Technical Support for help locating the logs.
- The reagent chiller temperature is too low or too high. The system displays an error.
- MiSeq FGx Maintenance Status:
 - Perform a maintenance wash and repeat sequencing.
 - Refill the wash tray and bottle for every wash.
 - Make all wash solutions daily.
- If the problem remains unresolved, reagents might not be performing as expected. Contact Verogen Technical Support.

Incomplete Run

The following factors can cause a run to stop prematurely.

- The HSC was not added during library prep. Add HSC per denature and dilute instructions and repeat the run.
- Cluster density is too low. See [Low Cluster Density \(previous page\)](#).
- If an error message is displayed, contact Verogen Technical Support.

Low Quality Scores for Read 1 or Read 2

The following factors can cause low quality scores for Read 1 and Read 2.

- Phasing and prephasing are too high. See [High Phasing and Prephasing \(previous page\)](#).
- The cluster density is too high. See [High Cluster Density \(page 49\)](#).
- If the problem remains unresolved, reagents might not be performing as expected. Contact Verogen Technical Support.

Low Quality Scores for Index 1 Read or Index 2 Read

The following factors can cause low quality scores for the Index 1 Read or the Index 2 Read.

- Phasing and prephasing are too high. See [Phasing and Prephasing \(page 55\)](#).
- The cluster density is too high. See [High Cluster Density \(page 49\)](#).
- Not enough samples were sequenced, so diversity in the index reads was low.
- If the problem remains unresolved, reagents might not be performing as expected. Contact Verogen Technical Support.

Low Reads per Sample

Review sample representation: If the samples exceed the sample read count (intensity) guideline, proceed with analysis. If the samples are below the guideline, the following factors are the likely cause.

- DNA input was too low or the DNA dilution was incorrect. If possible, repeat library prep with more DNA.
- A critical reagent was not added during library prep. Repeat the prep.
- LNB1 was insufficiently pipetted during normalization.
 - If the sample has a high number of alleles with quality control indicators, resequence the sample.

Universal Analysis Software – mtDNA Module Reference Guide

- If the sample has a low number of alleles with quality control indicators, consider diluting the libraries in the Pooled Normalized Libraries tube less and resequence the sample.
- If the run does not pass, repeat library normalization.
- HT1 was overpipetted during library prep. Repeat library normalization.
- The libraries were not heat-denatured during library prep. Perform the heat-denature step in the Diluted Normalized Libraries tube and repeat the run.

Low Reads per Sample for the HSC

Review sample representation: If the samples exceed the sample read count guideline, proceed with analysis. If the samples are below the guideline, the following factors are the likely cause:

- HSC was not added during library prep. Add HSC per denature and dilute instructions and repeat the run.
- HT1 was overpipetted during library prep. Repeat library normalization.
- HSC was not denatured with freshly diluted HP3. Make sure HSC is denatured with freshly diluted HP3 before adding it to the library, and then repeat the run.

Troubleshoot Analysis Errors

The following sections provide recommendations for troubleshooting potential analysis problems.

Analysis Error Message

Samples will error during analysis if there are no reads for the index combination.

1. Ensure that the correct index combination was used for the sample.
2. If the index combination is incorrect, follow steps in [Remove Samples from a Run \(page 25\)](#).
3. Add the sample back to the run with the correct index combination, following the steps in [Add Samples to a Run \(page 21\)](#).

The sample will automatically analyze.

If further errors occur, contact Verogen Technical Support

1. Navigate to the run folder, following the steps for obtaining the analysis folder path in [Review Run Activity \(page 39\)](#).
2. Navigate to the Log folder for that sample analysis.
3. Copy the log text file.
4. Navigate to D:\verogen\fuas\logs and copy the most recent log text file.

Quality Metrics Icons Are Unavailable

When all Quality Metrics icons are not available, the system might be unable to locate the data or the network where data are stored is disconnected. Use the following steps to troubleshoot.

1. Navigate to the run folder, following the steps for obtaining the analysis folder path in [Review Run Activity \(page 39\)](#).
2. Navigate up to the run folder (designated by a run number).
3. Ensure that the InterOp folder, the runInfo.xml file, and the runParameters.xml file are present.
 - If the folder is not present, copy this folder with the runInfo.xml and runParameters.xml files from the MiSeq, following the steps in [Review Run Metrics \(page 26\)](#).

Perform a Status Check

A status log in UAS aggregates information from a run, sample, or installation failure and a history of actions. Use the following instructions to export this information to send to Verogen when troubleshooting.

1. On the navigation rail, select  **Runs**. Then select a run.
2. On the Run Details page, select  **Run Details**, and then select **Run Status Check**.
3. In the Status Log Report dialog box, select **Copy**.
The copy command copies the content displayed under System Status and Most Recent Log to your clipboard.
4. Paste the copied content to an external file.
5. When finished, select **Close**.

Run Metrics

Accessing Run Metrics 54
 Quality Metrics 54
 Positive and Negative Controls 58
 Sample Representation 60

Accessing Run Metrics

Available in Projects after sequencing, run metrics provide run quality, positive and negative amplification control results, and sample representation for an overall view of run performance.


To open run metrics, select a QPN icon on the Project side sheet. A subset of these metrics, run quality and sample representation, are also available on Run Details pages.

Quality Metrics

Quality metrics monitor run quality during sequencing. After sequencing, the software preserves the quality metrics and makes them available in both Runs and Projects. In Runs, the quality metrics appear in a pop-up window. In Projects, they appear on a page.

Color-coded bars indicate results for each metric and show optimum ranges for cluster density, clusters passing filter, phasing, and prephasing. Short, vertical bars with corresponding icons indicate the quality of each read and the HSC.

Figure 8 Quality metrics overview



- A** Color-coded horizontal bars and values indicate run quality metrics.
- B** Thinner horizontal bars in dark gray indicate optimum ranges.
- C** Color-coded vertical bars and icons indicate the quality of each read and the HSC.

Run Quality Metrics

Run quality metrics include cluster density, clusters passing filter, phasing, and prephasing values. The following colors indicate overall quality.

Table 6 Colors for run quality metrics

Color	Indication
	Values are within the target range.
	Values require further evaluation of run quality.
	The target range for the metric.

Cluster Density

Cluster density is the number of clusters per square millimeter (K/mm²). For ForenSeq libraries, the target cluster density is 400–1650 K/mm². Values outside this range can still produce results sufficient for analysis. However, substantial deviations from the target range can impact other quality metrics and decrease the amount of usable data.

Clusters Passing Filter

Clusters passing filter is the percentage of clusters that passed the quality filter. This metric is based on the Illumina chastity filter, which measures quality and can detect low-quality base calls. Data appear after cycle 25.

The chastity of a base call is the ratio of the intensity of the greatest signal divided by the sum of the two greatest signals. If multiple base calls have a chastity value < 0.6 in the first 25 cycles, reads do not pass filter.

For ForenSeq libraries, the target value for a cluster passing filter is ≥ 80%. Values < 80% can still produce results that are sufficient for analysis. However, substantial deviations from the target can impact other quality metrics and decrease the amount of data.

Phasing and Prephasing

For both phasing and prephasing, lower percentages indicate higher quality run statistics. Phasing and prephasing values outside the target ranges can still produce results that are sufficient for analysis.

- Phasing shows the percentage of molecules in a cluster that fall behind the current cycle in Read 1. For ForenSeq libraries, the target phasing value is ≤ 0.25%.
- Prephasing shows the percentage of molecules in a cluster that jump ahead of the current cycle in Read 1. For ForenSeq libraries, the target prephasing value is ≤ 0.15%.

Read and Index Quality Metrics

The read and index quality metrics use the following colors and icons to indicate the status of each read and overall quality.

Table 7 Colors for read and index quality metrics

Color	Indication
	The average quality of assessed reads is within the target range.
	The average quality of assessed reads is outside the target range.
	The read has not yet occurred.

Table 8 Icons for read and index quality metrics

Icon	Indication
	The average quality of assessed reads is within the target range.
	The average quality of assessed reads is outside the target range.
	The read has not yet occurred.

Reads in a Run

A sequencing run completes up to four reads. Read 1 and Read 2 sequence the DNA template strands, and the Index 1 Read and Index 2 Read sequence the index adapters.

- **Read 1** — Read 1 sequencing primer is annealed to the template strand during cluster generation. RTA evaluates the first 50 cycles for quality.
- **Index 1 Read** — The Read 1 product is removed and the Index 1 sequencing primer is annealed to the same template strand as in Read 1. After index read preparation, the Index 1 Read is performed. RTA evaluates all eight cycles for quality.
- **Index 2 Read** — The Index 1 Read product is removed and the template anneals to the P5 primer grafted to the flow cell surface. The run proceeds through seven chemistry-only cycles without any imaging, followed by eight cycles of sequencing. RTA evaluates all eight cycles for quality.
- **Read 2** — The Index 2 Read product is extended to copy the original template strand. The original template strand is then removed and the Read 2 sequencing primer is annealed.

Cycles per Read

Quality metrics for each read appear after the read is complete. The number of cycles in each read depends on the assay.

Assay Type	Read	Cycles	Quality Metrics Appear
ForenSeq mtDNA Whole Genome	Read 1	1–201	Cycle 202
	Index 1 Read	202–209	Cycle 210
	Index 2 Read	210–217	Cycle 218
	Read 2	218–418	After cycle 418
ForenSeq mtDNA Control Region	Read 1	1–151	Cycle 152
	Index 1 Read	152–159	Cycle 160
	Index 2 Read	160–167	Cycle 168
	Read 2	168–318	After cycle 318

Human Sequencing Control

The software assesses the HSC and uses the following colors indicate the overall result.

Table 9 Colors for HSC results

Color	Indication
	The control meets or exceeds the criteria.
	The control does not meet at least one criterion.

When the HSC is orange, evaluate the run quality metrics to help determine whether to repeat the run. If the run quality metrics are green, indicating that values are within target ranges, you can interpret the sample with caution.

Loci

The HSC includes the loci in the following table. The software lists any loci that do not meet intensity or genotype concordance criteria.

Loci	151 or 201 Cycles	Loci	151 or 201 Cycles
D3S1358	✓	PentaE	✓
D5S818	✓	TPOX	✓
D7S820	✓	DYF387S1	
D8S1179	✓	DYS391	✓
D13S317	✓	DYS392	✓
D16S539	✓	DYS448	
D18S51	✓	DYS460	✓
D21S11		Y-GATA-H4	
FGA	✓	DXS7423	✓
PentaD	✓	DXS10074	
		DXS10103	✓

Positive and Negative Controls

The software assesses the same positive and negative control metrics for every run, regardless of assay type.

Table 10 Control assessments

Assay Type	Positive Amplification Control	Value
ForenSeq mtDNA Control Region	Control DNA HL60	Positions called
ForenSeq mtDNA Whole Genome	Control DNA HL60	Positions called

Positive Amplification Control

Positive amplification control metrics provide data for each sample identified as HL60. Every sample is analyzed and compared to this control. The software then lists each sample with values for how many reads passed filter and positions called. Colors indicate overall success. If applicable, a table displays discordant positions or loci.

Positions Called

Positions called shows the number of positions, including indels, with a base call out of the total number of positions for HL60. A short, vertical bar displays the following colors to indicate call success and concordance with HL60.

Table 11 Colors for positions for HL60

Color	Indication
	All positions have calls and are concordant.
	At least one position does not have a call or contains a discordant call.

Discordance Table

When a sample has calls that are discordant with the HL60 calls, the positive control metrics include a discordance table. For each discordant call, the table displays the position number or locus where the discordance occurred, the call observed for the sample, and the expected call.

The table displays up to 50 positions. When the number of discordant positions exceeds 50, use the scroll arrows to move through the table.

Pass Filter Reads

Pass filter reads displays the following colors to indicate whether the number of reads that passed filter meet the guideline of $\geq 50,000$ reads for ForenSeq mtDNA libraries. Reads below the guideline might still be usable. Review quality metrics and sample data to help make a determination.

Table 12 Colors for pass filter reads

Color	Indication
	The total number of pass filter reads exceeds the guideline.
	The total number of pass filter reads is below the guideline.
	The guideline for the number of reads that pass filter.

Negative Amplification Control

Negative control metrics provide data for each sample identified as a negative amplification control or reagent blank. The software lists each of these samples with the number of positions called and a color-coded indicator of overall success.

Positions Called

Positions called shows how many positions, including indels, a base call was made for. The short vertical bar uses the following colors to indicate call success.

Table 13 Colors for base calls of positions


Color	Indication
	No positions have a base call.
	At least one position has a base call.

Sample Representation

Sample representation is a metric that provides quantitative sample and run information. After sequencing, the software preserves sample representation data and makes it available in both Runs and Projects. In Runs, sample representation data appear in a pop-up window. In Projects, these data appear on a page.

Sample representation shows an index CV percentage for the run and the number of reads (intensity) and read distribution for each sample. These data are presented in a bar chart that also shows the sample read count guideline. Hovering over a bar displays the exact number of reads for the sample.

Figure 9 Sample Representation bar chart




- A** Sample name and exact number of reads for the sample
- B** Threshold of the read count guideline
- C** Bar indicating number of reads for the sample
- D** Exact number of reads for the sample (hovering)

Index CV

Above the sample representation bar chart, Index CV shows a percentage that represents the distribution of read counts for all samples in the run. CV is the coefficient of variation for the number of read counts across all indexes.

Figure 10 Index CV percentage



Sample Read Count Guideline

The sample read count guideline varies by assay type. As guidelines, these values are intended to help with quality reviews of the run and samples. Samples with fewer reads can still demonstrate complete coverage and provide enough data for interpretation. Review quality metrics and sample data to help make a determination.

Assay Type	Sample Read Count Guideline
ForenSeq mtDNA Control Region	50,000 reads per sample
ForenSeq mtDNA Whole Genome	400,000 reads per sample

The following colors indicate how sample read counts compare to the applicable sample read count guideline.

Table 14 Colors for sample read counts

Color	Indication
	The sample read count guideline.
	The sample meets or exceeds the guideline.
	The sample does not meet the guideline and requires careful interpretation.

Supporting Information

Special Handling of Sequencing Data 62

Special Handling of Sequencing Data

The following sections describe how the software handles special circumstances, such as ambiguous base call codes, when analyzing sequencing data from mtDNA samples.

Consensus Calls

The software supports variant calling for multiple bases at a DNA position. When multiple bases exceed the interpretation threshold, the software reports the International Union of Pure and Applied Chemistry (IUPAC) code as the consensus for the position.

IUPAC Code	Base
A	Adenine
C	Cytosine
G	Guanine
T	Thymine
R	A or G
Y	C or T
S	G or C
W	A or T
K	G or T
M	A or C
B	C, G, or T
D	A, G, or T
H	A, C, or T
V	A, C, or G

Ambiguous Calls

When data for the indel and a base call at a position exceed the IT, the software reports the ambiguous call code for an indel as the consensus for the position.

Base	Call
a	Adenine
c	Cytosine

Base	Call
g	Guanine
t	Thymine
r	a or g
y	c or t
s	g or c
w	a or t
k	g or t
m	a or c
b	c or g or t
d	a or g or t
h	a or c or t
v	a or c or g
n	Any base

Ambiguous Call Example

An example of an ambiguous call, the following figure shows an ambiguous deletion at position 501. Because an equal number of reads support a deletion call and a reference call of G, the software reports an ambiguous call of g as the consensus.

Figure 11 Ambiguous deletion at position 501

Coordinate	498	499	500	501	502	503
Call	A	C	T	g	A	C
Reads	A	C	T		A	C
	A	C	T		A	C
	A	C	T		A	C
	A	C	T		A	C
	A	C	T		A	C
	A	C	T		A	C
	A	C	T		A	C
	A	C	T		A	C
	A	C	T	G	A	C
	A	C	T	G	A	C
	A	C	T	G	A	C
	A	C	T	G	A	C
	A	C	T	G	A	C
	A	C	T	G	A	C

SWGDAM Nomenclature

The software follows SWGDAM nomenclature. For variant calling, differences between the reference and sample sequences are recorded as polymorphisms with the nucleotide position and DNA base difference indicated in the haplotype nomenclature (for example, 16182C).

Primer Trimming

The software trims primer binding sequences for ForenSeq mtDNA amplicons, which are based on a tiled primer design. This trimming allows unambiguous interrogation of mtDNA amplicon inserts without degenerate oligonucleotides. The mtDNA Variant Processor v1.0.0 application in Illumina BaseSpace Sequence Hub does not trim primers from amplicons generated from a tiled primer design.

Nuclear mtDNA Insertions

Nuclear DNA sequences homologous to mtDNA are often called nuclear mtDNA insertions (NUMTs). During library prep, these sections can be enriched alongside the respective mtDNA target. When performing alignment, UAS identifies reads that map closely to sequences in the Reference Human NumtS (RHNumtS) compilation and removes these reads before further analysis.^{5,6}

⁵ Mishmar, Dan, Eduardo Ruiz-Pesini, Martin Brandon, and Douglas C. Wallace, "Mitochondrial DNA-like Sequences in the Nucleus (NUMTs): Insights Into Our African Origins and the Mechanism of Foreign DNA Integration," *Human Mutation* 23 no. 2 (February 2004): 125–133, <https://doi.org/10.1002/humu.10304>.


⁶ Calabrese, Francesco Maria, Domenico Simone, Marcella Attimonelli "Primates and Mouse NumtS in the UCSC Genome Browser." Supplement, *BMC Bioinformatics* 13, no. S4 (March 2012): S15. <https://doi.org/10.1186/1471-2105-13-S4-S15>.

HVII C-Stretch Coverage

The C-stretch describes homopolymers of cytosine-stretches found in hypervariable region I (HVI) and hypervariable region II (HVII) of mtDNA. Primer pairs that amplify the region between position 262 and position 353 achieve coverage of the C-stretch (positions 303–315) in HVII. Sequences generated from the forward template strand have high accuracy and alignment. Sequences generated from the reverse template strand have high accuracy and alignment until reaching the C-stretch. Therefore, reads from the reverse strand that start sequencing at position 262 and do not meet alignment requirements are soft-clipped in UAS after position 303. The soft-clipped positions are not used for base calling. As a result, approximately half the coverage (read counts) are obtained for positions 304–353 relative to positions 262–303 (on the reverse strand).⁷⁻⁸

In the following figure, a Results page shows similar total read count and strand depth at positions 262–303 and 304–353. Positions 262–303, where reads from both strands are reported, have similar coverage. Total reads for positions 304–353 are about half that of positions 262–353.

Figure 12 Results for HL60 (100 pg gDNA) in UAS v2.1



⁷ Nakamura, Kensuke, Taku Oshima, Takuya Morimoto, Shun Ikeda, Hirofumi Yoshikawa, Yuh Shiwa, Shu Ishikawa, et al., "Sequence-specific error profile of Illumina sequencers," *Nucleic Acids Research* 39, no. 13 (July 2011): e90. <https://doi.org/10.1093/nar/gkr344>.

⁸ Allhoff, Manuel, Alexander Schönhuth, Marcel Martin, Ivan G. Costa, Sven Rahmann, and Tobias Marschall. Discovering motifs that induce sequencing errors," Supplement, *BMC Bioinformatics* 14, no. S5 (April 2013): S1. <https://doi.org/10.1186/1471-2105-14-S5-S1>.

Index

-
- A**
- activity setting 12, 47
 - adapter dimer 49
 - administrators
 - access 47
 - accounts 9
 - privileges 43, 47
 - settings 47
 - alignment 64
 - ambiguous calls 62
 - amplicons 64
 - modifying 44-45
 - analysis files 39
 - analysis methods 22, 29
 - default settings 43
 - editing 44
 - favorites 46
 - viewing results 29
 - analysis parameters 43, 61
 - analysis service 52
 - analysis status 16
 - analytical thresholds 14, 44
 - applications 8, 64
 - assay types 22, 43-44
 - AT 44-45
- B**
- BAM files 39
 - base call codes 62
 - base calls 14, 59-60, 62
 - modified 28
 - multiple 62
 - original 28
 - BaseSpace Sequence Hub 64
- C**
- C-stretch 64
 - Call and No Call filters 27
 - calls 14, 27
 - chastity filters 55
 - Chrome 9
 - cluster density 49-51, 54-55
 - clusters passing filter 50, 54-55
 - CMF headers 31-32, 35, 46
 - CODIS reports 31
- D**
- color coding 10
 - colors 25-27
 - commenting 39
 - comparison samples 27
 - concordance 59
 - configuration, software 43
 - contact information, Verogen 69
 - contributors 31-32
 - control software 8
 - controls 22, 26, 54
 - comparing 38
 - metrics 58
 - types 58
 - coordinates 14
 - coverage 14
 - coverage plot 14, 27
 - filters 27
 - custom assays 22
 - customer support 69
 - cycle numbers 57
- D**
- dark mode 43
 - dashboard 10
 - dashes, position viewer 14
 - data recording 47-48
 - database reports 31
 - default analysis settings 44
 - deletions 63
 - differences, samples 27
 - discordance 26, 59
 - diversity 51
 - DNA inputs 49-51
 - DNS 9
 - documentation 69
 - Domain Name System 9
 - downloading reports 13, 31-35
- E**
- EMPOP file format 33
 - EMPOP reports 31
 - enrichment 64
 - error messages 51
 - Excel workbooks 36
 - exceptions 62
- F**
- FASTA file format 35
 - filtering
 - reports 13
 - samples 27
 - filters 13, 23, 27
 - Call and No Call 27
 - flow cells 21
 - editing type 21
 - ForenSeq Analysis service 52
- G**
- genotype concordance 57-58
- H**
- help, technical 69
 - History icon 47-48
 - history setting 11-12, 39, 47
 - home-brew assays 22, 44
 - home page 10
 - hotspots 14
 - HSC 51-52
 - loci 58
 - quality metrics 25, 54
 - results 57
 - hyperlink run names 10
 - hyperlink sample names 10
 - hypervariable regions 64
- I**
- i5 read 56
 - i7 read 56
 - icons 10, 12, 16, 56
 - identical samples 27
 - importing samples 22
 - indels 62
 - Index 1 and Index 2 Reads, quality 56
 - index adapters 22, 56
 - index CV 60
 - Information icon 39
 - installing 9
 - intensity 51, 55, 57-58, 60
 - internet requirements 9
 - InterOp folder 26
 - interpretation thresholds 14, 44, 62

-
- IT 44-45
 - IUPAC codes 62
 - K
 - K/mm² 55
 - L
 - library prep kits 22
 - reference guides 49
 - library types 44-45
 - loci 58
 - logging in 10
 - logging out 10
 - lost passwords 9
 - M
 - matching samples 27
 - metrics, colors 25-26
 - minimum quality score 44-45
 - minimum read count 44-45
 - MiSeq FGx Control Software 8
 - MiSeq FGx System 8, 16
 - mtDNA navigator 14, 27
 - control region 14
 - filters 27
 - multiple bases called 62
 - N
 - navigation rails 10
 - negative amplification controls 13, 59
 - networks
 - requirements 9
 - troubleshooting 52
 - no calls 14, 27
 - non-directed sequencing 44
 - nucleotides 62
 - number of reads 15, 49
 - NUMTs 64
 - O
 - original calls 28
 - outputs 35
 - P
 - pass filter reads 26, 60-61
 - passwords, lost 9
 - paused runs 16
 - permissions 47
 - phasing and prephasing 26, 51, 54-55
 - plexities 51
 - position viewer 14, 28
 - dashes 14
 - filters 27
 - positions 14-15
 - selecting 14, 31
 - positions called 26
 - positive amplification controls 13, 58
 - primer-directed sequencing 44
 - primer pairs 64
 - primers 56
 - Print icon 39
 - printing 12
 - privileges, administrator 47
 - project summary, printing 39
 - projects 10
 - assigning samples 22, 24
 - creating 21-22, 24
 - descriptions 38
 - icons 12
 - naming 12, 38
 - summaries 10
 - Q
 - Q icon 10
 - QC indicators 13
 - QPN 12, 54
 - QPN icons 13
 - quality control 49, 61
 - quality filters 55
 - quality metrics 10-11, 13, 15, 25, 57
 - R
 - rails 10
 - rCRS 14, 35
 - read counts 14, 27, 61
 - read cycles 57
 - read definitions 56
 - read distributions 60
 - read numbers 27
 - Reads 1 and 2, quality 56
 - reagent blanks 59
 - reagent kits 21
 - reagent performance 50-52
 - Real-Time Analysis 8
 - reanalyzing samples 16, 29
 - recording data 47-48
 - reference sequences 35
 - reports
 - contents 35
 - downloading 13, 31-35
 - filtering 13
 - finding 13
 - types 31
 - worksheets 36
 - resolution 9
 - results 13
 - filtering 27
 - Results page 10
 - printing 39
 - RTA 56
 - run cycles 57
 - run details 10
 - run metrics 50, 57
 - run monitoring 11, 25, 54, 60
 - run status 16
 - runs 10
 - general information 11
 - metrics 49, 57
 - performance 11, 25
 - quality 26, 54
 - status 25
 - S
 - safety data sheets 69
 - sample comparisons, reporting 38
 - sample information 13, 21, 31, 36
 - Sample Overview section, printing 39
 - sample position, changing 14
 - sample read count guideline 26, 49-52, 60
 - sample representation 11, 13, 25-26, 49-52, 54
 - sample results 13
 - sample sheets 21-22
 - samples
 - adding information 21
 - analyzing 14
 - assigning to projects 21-22, 24, 29, 38
 - differences 27
 - filtering 27
 - importing 22
 - intensity 49
 - locating 23

Universal Analysis Software – mtDNA Module Reference Guide

- modifying 24
- no differences 27
- reanalyzing 29
- removing from runs 25
- searching 23
- status 16
- samples of interest 27
- Samples page
 - printing 39
- screen resolution 9
- SDS 69
- searching
 - reports 13
 - samples 23
- selected comparison samples 27
- sequences 14, 35
- sequencing processes 56
- servers 8, 39
- soft clipping 64
- software reset 9
- sorting 13, 23
- specimens 31-32, 35
- star icons 46
- stopped runs 16
- strand depth 14-15
- support, library prep 49
- system dashboard 10
- system events, definition 47
- T
- technical assistance 26, 69
- technical notes 43
- third-party assays 22, 44
- thresholds 43-45, 61
- tiling 64
- total read count 15
- total reads 14
- troubleshooting 52
- troubleshooting files 26
- U
- uninstalling 9
- user access 47
- user actions, definition 47
- V
- validations 44
- variant calls 43, 62

- variants 14, 35
- VCF files 39
- W
- washes 50
- web browser 9

Technical Support

For technical assistance, contact Verogen Technical Support.

Contact Information

Address 11111 Flintkote Avenue
San Diego CA 92121 USA

Website www.verogen.com

Email techsupport@verogen.com

Telephone +1.833.837.6436 toll-free (North America)
+1.858.285.4101 tel
+44 (0) 208 054 8706 (United Kingdom)

Safety data sheets (SDS)—Available for download from verogen.com/product-documentation.

Product documentation—Available for download from verogen.com/product-documentation.

Meet any challenge

Verogen is a dedicated developer of human identification products for sequencing and analysis of forensic genomic samples. Working closely with the forensics community, Verogen places exceptional value on flexible, scalable solutions that deliver results when you need them most.

Verogen

+1.833.837.6436 toll-free
(North America)

+1.858.285.4101
(outside North America)

techsupport@verogen.com

www.verogen.com