

Issue Date 01-Jun-2015

Revision Date 18-July-2018

Version 1.1

1. IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND OF THE COMPANY/UNDERTAKING**Product identifier****Product Name** LP#-LNS2, LIBRARY NORMALIZATION STORAGE BUFFER**Other means of identification****Product Code** 15041599**Synonyms** None**Recommended use of the chemical and restrictions on use****Recommended Use** reagent**Details of the supplier of the safety data sheet****Manufacturer Address**
Verogen, Inc.
11111 Flintkote Avenue
San Diego, CA 92121 USA
www.verogen.com**Company Phone Number**
1.833.837.6436
1.858.285.4101**E-mail address**
info@verogen.com**Emergency telephone number****24 Hour Emergency Phone Number** 1.858.285.4101**2. HAZARDS IDENTIFICATION****Classification****OSHA Regulatory Status**

This chemical is not considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.122).

Not a dangerous substance or mixture according to the Globally Harmonized System (GHS)

Label elements**Emergency Overview****Signal word**

None

Hazard statements

Not Hazardous

The product contains no substances which at their given concentration, are considered to be hazardous to health

Appearance colorless liquid**Physical state** Liquid**Odor** Odorless**Precautionary Statements - Prevention**

Handle in accordance with good industrial hygiene and safety practice

Precautionary Statements - Response

Specific treatment (see supplemental first aid instructions on this label)

Precautionary Statements - Storage

Store away from incompatible materials

Precautionary Statements - Disposal

Dispose of in accordance with federal, state and local regulations

Hazards not otherwise classified (HNOC)

Not Applicable

Other Information

Unknown Acute Toxicity 4.1 % of the mixture consists of ingredient(s) of unknown toxicity

3. COMPOSITION/INFORMATION ON INGREDIENTS

The product contains no substances which at their given concentration, are considered to be hazardous to health.

Chemical Name	CAS No	Weight-%	Trade Secret
Sodium hydroxide	1310-73-2	0.1 - 1	*

*The exact percentage (concentration) of composition has been withheld as a trade secret.

4. FIRST AID MEASURES**First aid measures****General advice**

Use first aid treatment according to the nature of the injury.

Eye contact

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Get medical attention if irritation develops and persists.

Skin Contact

Wash skin with soap and water. Get medical attention if irritation develops and persists.

Inhalation

If symptomatic, move to fresh air. Get medical attention if symptoms persist.

Ingestion

Get medical advice/attention.

Self-protection of the first aider

Ensure that medical personnel are aware of the material(s) involved and take precautions to protect themselves.

Most important symptoms and effects, both acute and delayed**Symptoms**

No information available.

Indication of any immediate medical attention and special treatment needed**Note to physicians**

Treat symptomatically.

5. FIRE-FIGHTING MEASURES**Suitable extinguishing media**

Water, Water spray (fog), Foam, Dry chemical, Carbon dioxide (CO2)

Unsuitable extinguishing media

None known.

Specific hazards arising from the chemical

None known.

Hazardous combustion products

Carbon oxides. Oxides of sulfur.

Explosion data

Sensitivity to Mechanical Impact None.
Sensitivity to Static Discharge None.

Protective equipment and precautions for firefighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency procedures

Personal precautions Use personal protective equipment as required.
For emergency responders Use personal protection recommended in Section 8. Keep unnecessary personnel away.

Environmental precautions

Environmental precautions See Section 12 for additional Ecological Information.

Methods and material for containment and cleaning up

Methods for containment Prevent further leakage or spillage if safe to do so. Dike far ahead of liquid spill for later disposal.
Methods for cleaning up Use appropriate personal protective equipment (PPE). Contain and/or absorb spill with inert material (specify sand, vermiculite or other appropriate material)*, then place in suitable container. Prevent runoff from entering drains, sewers or waterways.

7. HANDLING AND STORAGE

Precautions for safe handling

Advice on safe handling Use personal protective equipment as required. Handle in accordance with good industrial hygiene and safety practice.

Conditions for safe storage, including any incompatibilities

Storage Conditions Keep container tightly closed. Store away from incompatible materials. Store at -15° C to -25° C (5° F to -13° F).
Incompatible materials Strong oxidizing agents. Acids.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Control parameters

Exposure Guidelines

Chemical Name	ACGIH TLV	OSHA PEL	NIOSH IDLH
Sodium hydroxide 1310-73-2	Ceiling: 2 mg/m ³	TWA: 2 mg/m ³ (vacated) Ceiling: 2 mg/m ³	IDLH: 10 mg/m ³ Ceiling: 2 mg/m ³

NIOSH IDLH Immediately Dangerous to Life or Health

Other Information Vacated limits revoked by the Court of Appeals decision in AFL-CIO v. OSHA, 965 F.2d 962 (11th Cir., 1992).

Appropriate engineering controls

Engineering Controls Showers
Eyewash stations

Ventilation systems.

Individual protection measures, such as personal protective equipment

Eye/face protection	Wear safety glasses with side shields (or goggles).
Skin and body protection	Wear chemical-resistant gloves, footwear and protective clothing appropriate for risk of exposure. Contact glove manufacturer for specific information.
Respiratory protection	If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved respiratory protection should be worn. Positive-pressure supplied air respirators may be required for high airborne contaminant concentrations. Respiratory protection must be provided in accordance with current local regulations.
General Hygiene Considerations	Handle in accordance with good industrial hygiene and safety practice. Wash after handling this material and before eating, drinking and/or smoking. Regular cleaning of equipment, work area and clothing is recommended.

9. PHYSICAL AND CHEMICAL PROPERTIES

Information on basic physical and chemical properties

Physical state	Liquid
Appearance	colorless liquid
Odor	Odorless
Color	colorless

<u>Property</u>	<u>Values</u>	<u>Remarks • Method</u>
pH	Alkaline	
Melting point / freezing point	No data available	
Boiling point / boiling range		
Flash point		
Evaporation rate	No data available	
Flammability (solid, gas)	No data available	
Flammability Limit in Air		
Upper flammability limit:	No data available	
Lower flammability limit:	No data available	
Vapor pressure	No data available	
Vapor density	No data available	
Specific Gravity	No data available	
Water solubility	No data available	
Solubility in other solvents	No data available	
Partition coefficient	No data available	
Autoignition temperature	No data available	
Decomposition temperature	No data available	
Kinematic viscosity	No data available	
Dynamic viscosity	No data available	
Explosive properties	No unusual fire or explosion hazards noted.	
Oxidizing properties	No data available	

Other Information

Softening point	No data available
Molecular weight	No data available
VOC Content (%)	No data available
Density	No data available
Bulk density	No data available

10. STABILITY AND REACTIVITY

Reactivity

The product is stable and non-reactive under normal conditions of use, storage and transport.

Chemical stability

Stable under normal conditions. Stable under recommended storage conditions.

Possibility of Hazardous Reactions

None under normal processing.

Hazardous polymerization

Hazardous polymerization does not occur.

Conditions to avoid

None known.

Incompatible materials

Strong oxidizing agents. Acids.

Hazardous Decomposition Products

None known based on information supplied.

11. TOXICOLOGICAL INFORMATION

Information on likely routes of exposure**Product Information**

Low hazard for usual industrial or commercial handling by trained personnel

Inhalation

Under normal conditions of intended use, this material is not expected to be an inhalation hazard.

Eye contact

Direct contact with eyes may cause temporary irritation.

Skin Contact

Prolonged contact may cause redness and irritation.

Ingestion

No harmful effects expected in amounts likely to be ingested by accident.

Information on toxicological effects**Symptoms**

None known.

Delayed and immediate effects as well as chronic effects from short and long-term exposure**Skin corrosion/irritation**

Not classified.

Serious eye damage/eye irritation

Not classified.

Sensitization

Not classified.

Germ cell mutagenicity

Not classified.

Carcinogenicity

This product does not contain any carcinogens or potential carcinogens as listed by OSHA, IARC or NTP.

Reproductive toxicity

Not classified.

STOT - single exposure

Not classified.

STOT - repeated exposure

Not classified.

Aspiration hazard

Not classified.

Other adverse effects

No information available.

Numerical measures of toxicity - Product Information**Unknown Acute Toxicity**

4.1 % of the mixture consists of ingredient(s) of unknown toxicity

12. ECOLOGICAL INFORMATION

Ecotoxicity

4.2 % of the mixture consists of component(s) of unknown hazards to the aquatic environment

Chemical Name	Algae/aquatic plants	Fish	Crustacea
Sodium hydroxide 1310-73-2	-	45.4 mg/L: LC50 96 h Oncorhynchus mykiss static	-

Persistence and degradability

Not a PBT or vPvB substance or mixture.

Bioaccumulation

No data available

Mobility

No data available.

Other adverse effects

No data available

13. DISPOSAL CONSIDERATIONS

Waste treatment methods

Disposal of wastes Disposal should be in accordance with applicable regional, national and local laws and regulations.

Contaminated packaging Do not reuse container. Disposal should be in accordance with applicable regional, national and local laws and regulations.

This product contains one or more substances that are listed with the State of California as a hazardous waste.

Chemical Name	California Hazardous Waste Status
Sodium hydroxide 1310-73-2	Toxic Corrosive

14. TRANSPORT INFORMATION

DOT	Not regulated
TDG	Not regulated
MEX	Not regulated
ICAO (air)	Not regulated
IATA	Not regulated
IMDG	Not regulated
RID	Not regulated
ADR	Not regulated
ADN	Not regulated

15. REGULATORY INFORMATION

International Inventories

TSCA	Complies
DSL/NDSL	Complies
EINECS/ELINCS	Does not comply

ENCS Does not comply
IECSC Complies
KECL Does not comply
PICCS Complies
AICS Complies

Legend:

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory
DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List
EINECS/ELINCS - European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances
ENCS - Japan Existing and New Chemical Substances
IECSC - China Inventory of Existing Chemical Substances
KECL - Korean Existing and Evaluated Chemical Substances
PICCS - Philippines Inventory of Chemicals and Chemical Substances
AICS - Australian Inventory of Chemical Substances

US Federal Regulations

SARA 313

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372

SARA 311/312 Hazard Categories

Acute health hazard No
Chronic Health Hazard No
Fire hazard No
Sudden release of pressure hazard No
Reactive Hazard No

CWA (Clean Water Act)

This product contains the following substances which are regulated pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42)

Chemical Name	CWA - Reportable Quantities	CWA - Toxic Pollutants	CWA - Priority Pollutants	CWA - Hazardous Substances
Sodium hydroxide 1310-73-2	1000 lb	-	-	X

CERCLA

This material, as supplied, contains one or more substances regulated as a hazardous substance under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40 CFR 302)

Chemical Name	Hazardous Substances RQs	CERCLA/SARA RQ	Reportable Quantity (RQ)
Sodium hydroxide 1310-73-2	1000 lb	-	RQ 1000 lb final RQ RQ 454 kg final RQ

US State Regulations

California Proposition 65

This product does not contain any Proposition 65 chemicals

U.S. State Right-to-Know Regulations

Chemical Name	New Jersey	Massachusetts	Pennsylvania
Water 7732-18-5	-	-	X
Sodium hydroxide 1310-73-2	X	X	X

U.S. EPA Label Information

EPA Pesticide Registration Number Not Applicable

16. OTHER INFORMATION

<u>NFPA</u>	Health hazards 0	Flammability 0	Instability 0	Physical and Chemical Properties -
<u>HMIS</u>	Health hazards 0	Flammability 0	Physical hazards 0	Personal protection X

Prepared By Verogen Inc.
 Issue Date 01-Jun-2015
 Revision Date 18-July-2018

Revision Note
 Not Applicable.

Disclaimer

Although reasonable care has been taken in the preparation of this document, Verogen makes no warranties or representations of any kind as to the accuracy or completeness of the information contained herein, and Verogen assumes no responsibility regarding the suitability of this information or the product for the user's intended purposes or for the consequences of use of this information or the use of the product. The user of this (M)SDS and product should make a determination as to the suitability of the information and the product for their particular purpose(s). To the extent provided by applicable law, Verogen shall not be held liable for any damage or expenses of any kind resulting from handling, use, storage, or disposal of the product.

End of Safety Data Sheet