

VD2020019 LP#-SPB

Issuing Date 02-Oct-2018

Revision Date 18-May-2020

Revision A

1. IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND OF THE COMPANY/UNDERTAKING

Product identifier

Product Name LP#-SPB, SAMPLE PURIFICATION BEADS 10 ml

Other means of identification

Product Codes(s) 15041032

Synonyms None

Recommended use of the chemical and restrictions on use

Recommended Use reagent

Details of the supplier of the safety data sheet

Manufacturer Address

Illumina, Inc.
5200 Illumina Way
San Diego, CA 92122 USA

www.illumina.com

Company Phone Number

1-858-809-ILMN
1-858-202-4566

E-mail address

SDS@illumina.com

Emergency telephone number

24 Hour Emergency Phone Number 1-760-476-3962

2. HAZARDS IDENTIFICATION

GHS Classification

OSHA Regulatory Status

This chemical is not considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200).

Not a dangerous substance or mixture according to the Globally Harmonized System (GHS).

Label elements

Emergency Overview

Signal word

None

Hazard statements

Not Hazardous

The product contains no substances which at their given concentration are considered to be hazardous to health.

Appearance Brown liquid

Physical state Liquid

Odor Odorless

Precautionary Statements – Prevention

Handle in accordance with good industrial hygiene and safety practice.

Precautionary Statements – Storage

Store away from incompatible materials.

Precautionary Statements – Disposal

Dispose of in accordance with federal, state and local regulations.

Hazards not otherwise classified (HNOC)

Not applicable.

Other Information

Unknown acute toxicity 100 % of the mixture consists of ingredient(s) of unknown toxicity.

3. COMPOSITION/INFORMATION ON INGREDIENTS

The product contains no substances which at their given concentration are considered to be hazardous to health.

4. FIRST AID MEASURES**Description of first aid measures**

General Advice	Use first aid treatment according to the nature of the injury.
Eye Contact	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Get medical attention if irritation develops and persists.
Skin Contact	Wash skin with soap and water. Get medical attention if irritation develops and persists.
Inhalation	If symptomatic, move to fresh air. Get medical attention if symptoms persist.
Ingestion	Get medical advice/attention.
Self-protection of the first aider	Ensure that medical personnel are aware of the material(s) involved and take precautions to protect themselves.

Most important symptoms and effects, both acute and delayed**Symptoms** No information available.**Indication of any immediate medical attention and special treatment needed****Note to physicians** Treat symptomatically.**5. FIRE-FIGHTING MEASURES****Suitable extinguishing media**

Water, Water spray (fog), Foam, Dry chemical, Carbon dioxide (CO2).

Unsuitable extinguishing media

None known.

Specific hazards arising from the chemical

None known.

Hazardous combustion products

No information available.

Explosion data

Sensitivity to Mechanical Impact	None
Sensitivity to Static Discharge	None

Protective equipment and precautions for firefighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

6. ACCIDENTAL RELEASE MEASURES**Personal precautions, protective equipment and emergency procedures**

Personal precautions	Use personal protective equipment as required.
For emergency responders	Use personal protection recommended in Section 8. Keep unnecessary personnel away.

Environmental precautions

Environmental precautions Avoid discharge into drains, water courses, or onto the ground.

Methods and material for containment and cleaning up

Methods for containment Prevent further leakage or spillage if safe to do so. Dike far ahead of liquid spill for later disposal.

Methods for cleaning up Use personal protective equipment as required. Contain and collect spillage with non-combustible absorbent material, (e.g. sand, earth, diatomaceous earth, vermiculite) and place in container for disposal according to local / national regulations (see Section 13). Large spill. Prevent runoff from entering drains, sewers, or streams. Dike far ahead of liquid spill for later disposal.

Prevention of secondary hazards Clean contaminated objects and areas thoroughly observing environmental regulations.

7. HANDLING AND STORAGE

Precautions for safe handling

Advice on safe handling Use personal protective equipment as required. Handle in accordance with good industrial hygiene and safety practice.

Conditions for safe storage, including any incompatibilities

Storage Conditions Keep container tightly closed in a dry and well-ventilated place. Keep/store only in original container. Store away from incompatible materials.

Incompatible materials Strong oxidizing agents.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Control parameters**Exposure Guidelines**

Chemical Name	ACGIH TLV	OSHA PEL	NIOSH IDLH
Sodium azide 26628-22-8	Ceiling: 0.29 mg/m ³ NaN ₃ 0.11 ppm Hydrazoic acid vapor	(vacated) S* (vacated) Ceiling: 0.1 ppm HN ₃ (vacated) Ceiling: 0.3 mg/m ³ NaN ₃	Ceiling: 0.1 ppm HN ₃ Ceiling: 0.3 mg/m ³ NaN ₃

NIOSH IDLH Immediately Dangerous to Life or Health

Other Information Vacated limits revoked by the Court of Appeals decision in AFL-CIO v. OSHA, 965 F.2d 962 (11th Cir., 1992).

Appropriate engineering controls

Engineering Controls Showers
Eyewash stations
Ventilation systems

Individual protection measures, such as personal protective equipment

Eye/face protection Wear safety glasses with side shields (or goggles).

Skin and body protection Wear chemical-resistant gloves, footwear and protective clothing appropriate for risk of exposure. Contact glove manufacturer for specific information.

Respiratory protection If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved respiratory protection should be worn. Positive-pressure supplied air respirators may be required for high airborne contaminant concentrations. Respiratory protection must be provided in accordance with current local regulations.

General Hygiene Considerations Handle in accordance with good industrial hygiene and safety practice. Wash after handling this material and before eating, drinking and/or smoking. Regular cleaning of equipment, work area and clothing is recommended.

9. PHYSICAL AND CHEMICAL PROPERTIES

Information on basic physical and chemical properties

Physical state	Liquid
Appearance	Brown liquid
Odor	Odorless
Color	Brown
Odor threshold	No information available

<u>Property</u>	<u>Values</u>	<u>Remarks • Methods</u>
pH	8.4	
Melting point / freezing point	No data available	
Boiling point / boiling range		
Flash point		
Evaporation rate	No data available	
Flammability (solid, gas)	No data available	
Flammability Limit in Air		
Upper flammability or explosive limits:	No data available	
Lower flammability or explosive limits:	No data available	
Vapor pressure	No data available	
Vapor density	No data available	
Specific gravity	1.127	
Water solubility	Miscible	
Solubility in other solvents	No data available	
Partition coefficient	No data available	
Autoignition temperature	No data available	
Decomposition temperature	No data available	
Kinematic viscosity	No data available	
Dynamic viscosity	No data available	
Explosive properties	No unusual fire or explosion hazards noted.	
Oxidizing properties	No data available	

Other Information

Softening point	No data available
Molecular weight	No data available
VOC content (%)	No data available
Liquid density	No data available
Bulk density	No data available

10. STABILITY AND REACTIVITY

Reactivity

Chemical stability

Stable under normal conditions. Stable under recommended storage conditions.

Possibility of Hazardous Reactions

None under normal processing.

Hazardous polymerization

Hazardous polymerization does not occur.

Conditions to avoid

None known.

Incompatible materials

Strong oxidizing agents.

Hazardous decomposition products

None known based on information supplied.

11. TOXICOLOGICAL INFORMATION

Information on likely routes of exposure

Product Information	Low hazard for usual industrial or commercial handling by trained personnel.
Inhalation	Under normal conditions of intended use, this material is not expected to be an inhalation hazard.
Eye Contact	Direct contact with eyes may cause temporary irritation.
Skin Contact	Prolonged contact may cause redness and irritation.
Ingestion	No harmful effects expected in amounts likely to be ingested by accident.

Information on toxicological effects

Symptoms None known.

Delayed and immediate effects as well as chronic effects from short and long-term exposure

Skin corrosion/irritation	Not classified.
Serious eye damage/eye irritation	Not classified.
Skin sensitization	Not classified.
Respiratory sensitization	Not classified.
Germ cell mutagenicity	Not classified.
Carcinogenicity	This product does not contain any carcinogens or potential carcinogens as listed by OSHA, IARC or NTP.
Reproductive toxicity	Not classified.
STOT – single exposure	Not classified.
STOT – repeated exposure	Not classified.
Aspiration hazard	Not classified.
Other adverse effects	No information available.

Numerical measures of toxicity – Product Information

Unknown acute toxicity 100 % of the mixture consists of ingredient(s) of unknown toxicity.

12. ECOLOGICAL INFORMATION

Ecotoxicity

0 % of the mixture consists of components(s) of unknown hazards to the aquatic environment.

Chemical Name	Algae/aquatic plants	Fish	Crustacea
Sodium azide 26628-22-8	-	0.7 mg/L: LC50 96 h Lepomis macrochirus 0.8 mg/L: LC50 96 h Oncorhynchus mykiss 5.46 mg/L: LC50 96 h Pimephales promelas flow-through	-

Persistence and degradability

Not a PBT or vPvB substance or mixture.

Bioaccumulation

No information available.

Mobility

No data available.

Other adverse effects

No data available.

Ozone

Not classified.

Ozone depletion potential (ODP)

No information available.

13. DISPOSAL CONSIDERATIONS

Waste treatment methods

Disposal of wastes	Disposal should be in accordance with applicable regional, national and local laws and regulations.
Contaminated packaging	Do not reuse container. Disposal should be in accordance with applicable regional, national and local laws and regulations.
US EPA Waste Number	P105

This product contains one or more substances that are listed with the State of California as a hazardous waste.

14. TRANSPORT INFORMATION

DOT	Not regulated
TDG	Not regulated
MEX	Not regulated
ICAO (air)	Not regulated
IATA	Not regulated
IMDF	Not regulated
RID	Not regulated
ADR	Not regulated
ADN	Not regulated

15. REGULATORY INFORMATION

International Inventories

TSCA	Does not comply
DSL/NDSL	Does not comply
EINECS/ELINCS	Does not comply
ENCS	Does not comply
IECSC	Does not comply
KECL	Does not comply
PICCS	Does not comply
AICS	Does not comply
INSQ	This product complies with INSQ

Legend

TSCA—United States Toxic Substances Control Act Section 8(b) Inventory
DSL/NDSL—Canadian Domestic Substances List/Non-Domestic Substances List
EINECS/ELINCS—European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances
ENCS—Japan Existing and New Chemical Substances
IECSC—China Inventory of Existing Chemical Substances
KECL—Korean Existing and Evaluated Chemical Substances
PICCS—Philippines Inventory of Chemicals and Chemical Substances
AICS—Australian Inventory of Chemical Substances
INSQ—Mexican National Inventory of Chemical Substances

US Federal Regulations

SARA 313

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

SARA 311/312 Hazard Categories

Acute health hazard No
 Chronic Health Hazard No
 Fire hazard No
 Sudden release of pressure hazard No
 Reactive Hazard No

CWA (Clean Water Act)

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42).

CERCLA

This material, as supplied, contains one or more substances regulated as a hazardous substance under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40 CFR 302).

Chemical Name	Hazardous Substances RQs	CERCLA/SARA RQ	Reportable Quantity (RQ)
Sodium azide 26628-22-8	1000 lb	1000 lb	RQ 1000 lb final RQ RQ 454 kg final RQ

US State Regulations

California Proposition 65

This product does not contain any Proposition 65 chemicals.

U.S. State Right-to-Know Regulations

Chemical Name	New Jersey	Massachusetts	Pennsylvania
Sodium azide 26628-22-8	X	X	X

U.S. EPA Label Information

EPA Pesticide Registration Number Not applicable

16. OTHER INFORMATION

NFPA Health Hazards 1 Flammability 0 Instability 0 Physical and Chemical Properties -
 HMIS Health Hazards 1 Flammability 0 Physical Hazards 0 Personal Protection X

Prepared By Illumina Inc.
 Issuing Date 02-Oct-2018
 Revision Date 18-May-2020
 Revision Note Not applicable

Disclaimer

Although reasonable care has been taken in the preparation of this document, Illumina makes no warranties or representations of any kind as to the accuracy or completeness of the information contained herein, and Illumina assumes no responsibility regarding the suitability of this information or the product for the user's intended purposes or for the consequences of use of this information or the use of the product. The user of this (M)SDS and product should make a determination as to the suitability of the information and the product for their particular purpose(s). To the extent provided by applicable law, Illumina shall not be held liable for any damage or expenses of any kind resulting from handling, use, storage, or disposal of the product.

End of Safety Data Sheet